Resolute Onyx™

ZOTAROLIMUS-ELUTING CORONARY STENT SYSTEM

Resolute Onyx DES is designed to build on the market-leading clinical performance that Resolute DES demonstrated in the RESOLUTE All Comers clinical trial. Resolute Onyx DES offers the superior deliverability of the Integrity platform without compromising the drug elution characteristics of Resolute Integrity DES.

Expanded Size Matrix

Diameter (mm)	Stent Length (mm)									
	8	12	15	18	22	26	30	34	38	
2.00	8	12	15	18	22	26	30			
2.25	8	12	15	18	22	26	30	34	38	
2.6	2.50 ✓	8	12	15	18	22 ✓	26	30	34	38
2.75	8	12	15	18	22 ✓	26	30	34	38	
3.00	8	12	15	18	22	26	30	34	38	
3.50	8	12	15	18	22	26	30	34	38	
4.00	8	12	15	18	22	26	30	34	38	
4.50		12	15	18	22	26	30			
5.00		12	15	18	22	26	30			

Shaded cells highlight modifications to the size matrix: 2.00-mm stents and 34- and 38-mm lengths expand the size matrix; an 8-mm length replaces the 9-mm Resolute Integrity length; and a 15-mm length replaces the 14-mm Resolute Integrity length. 2.00-2.50 mm stents can be expanded to 3.25 mm; 2.75-3.00 mm stents can be expanded to 3.75 mm; 3.50-4.00 mm stents can be expanded to 4.75 mm; 4.50-5.00 mm stents can be expanded to 5.75 mm.

Technical Improvements


Technical Feature	Resolute Integrity DES	Resolute Onyx DES
Stent design	Integrity™ BMS Continuous Sinusoid Technology	Continuous Sinusoid Technology with Core Wire Technology
Drug and polymer	BioLinx and zotarolimus	No change
Delivery system	MicroTrac	Resolute Onyx
Catheter distal shaft O.D. (F)	2.7	2.7 (2.00-4.00 mm) 3.2 (4.50-5.00 mm)
Catheter distal shaft O.D. (in.)	0.036	0.036 (2.00-4.00 mm) 0.042 (4.50-5.00 mm)

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Changes in Selected Technical Specifications

Technical Feature	Resolute Integrity DES	Resolute Onyx DES
Crowns	Small vessel (2.25-2.75 mm): 7.5 Medium vessel (3.00-4.00 mm): 9.5	Small vessel (2.00-2.50 mm): 6.5 Medium vessel (2.75-3.00 mm): 8.5 Large vessel (3.50-4.00 mm): 9.5 Extra-large vessel (4.50-5.00 mm): 10.5
Stent material	✓ 2.2 Cobalt alloy conforming to ASTM F562-02	Shell: No change Core: Pt-Ir
Strut thickness	91 μ (0.0036 in.)	✓ 2.5 81 μ (0.0032 in.) (2.00-4.00 mm) 91 μ (0.0036 in.) (4.50-5.00 mm)
Polymer coating	BioLinx (C10:C19:PVP = 27:63:10)	No change
Base coating	Parylene C	No change
Drug density	2.4 ✓ Zotarolimus (~1.6 μg/mm ²)	No change
Delivery system	MicroTrac	Resolute Onyx
Balloon material	Fulcrum™	Improved
Platform	RX	No change
Catheter length (cm)	2.8 ✓ 140	No change
Nominal pressure (atm)	9	12 (2.00-5.00 mm)
Rated burst pressure (atm)	16	18 (2.00-4.00 mm) 16 (4.50-5.00 mm)
Crimping technology	Secure Technology	No change
Marker band material	Gold	Platinum iridium

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10 PATIENT SELECTION AND TREATMENT

See also **Section 5.5 - Use in Special Populations**. The risks and benefits described above should be carefully considered for each patient before use of the Resolute Onyx™ System. Factors to be utilized for patient selection should include an assessment of the risk of prolonged anticoagulation. Administration of P2Y₁₂ platelet inhibitor is recommended pre-procedure and for at least 6 months in stable ischemic heart disease patients and for at least 12 months in patients with acute coronary syndrome (ACS). In patients at higher risk of bleeding, DAPT discontinuation may be reasonable after 3 months in stable patients or 6 months in ACS patients (see **Section 5.1 - Pre- and Post-Procedure Antiplatelet Regimen**). Aspirin should be administered concomitantly with an approved antiplatelet medication and then continued indefinitely. The safety and effectiveness of the Resolute Onyx stent have not been evaluated in patients at high bleeding risk or those with contraindicated anticoagulation therapy.

11 PATIENT COUNSELING INFORMATION

Physicians should consider the following in counseling the patient about this product:

- Discuss the risks associated with stent placement
- Discuss the risks associated with a zotarolimus-eluting stent implant
- Discuss the risks/benefits issues for this particular patient
- Discuss alteration to current lifestyle immediately following the procedure and over the long term
- Discuss the risks of early discontinuation of the antiplatelet therapy

The following patient materials will be provided to physicians to educate their patients about the options available for treating coronary artery disease and provide contact information to the patient after their stent implant procedure:

- A Patient Guide which includes information on the Resolute Onyx™ Zotarolimus-Eluting Coronary Stent System, coronary artery disease, and the stent implantation procedure.
- A Stent Patient Implant Card that includes patient information, stent implant information and MRI guidelines. All patients should be instructed to keep this card in their possession at all times for procedure/stent identification.

12 HOW SUPPLIED

3.1 ✓

STERILE: This product is sterilized with ethylene oxide (EO) and is nonpyrogenic. Do not use if the package is opened or damaged. Do not resterilize. If the product or package is opened or damaged, return to Medtronic Returned Goods. Contact your local Medtronic Representative for return information.

CONTENTS: Package contains one (1) Resolute Onyx™ Zotarolimus-Eluting Coronary Stent mounted on either a Rapid Exchange (RX) or an Over-the-Wire (OTW) stent delivery system.

STORAGE: Store in the original container. Store at 25°C (77°F); excursions permitted to 15 - 30°C (59 - 86°F). Use by the "Use By" date noted on the package.

DISPOSAL INSTRUCTIONS: After use, dispose of product and packaging in accordance with hospital, administrative and/ or local government policy.

13 DIRECTIONS FOR USE

13.1 Access to Package Holding Sterile Stent Delivery System

Remove the stent delivery system from the package. Special care must be taken not to handle the stent or in any way disrupt its placement on the balloon. This is most important during catheter removal from packaging, placement over guidewire, and advancement

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
- Lubricious hydrophilic coating for reduced drag
- 0.02 mm (2.1 Fr) optimal stent
- PowerFlex™ technology enhances deliverability
- Resident hypertube for high stent column strength
- Enhanced tubular material malleability and bendability
- 0.01 mm (0.75 Fr) distal shaft
- 4.50-5.00 mm (3.2 Fr)
- Flattens in diameter for bundle or device visibility
- Reduced catheter profile under the stent enables lower crossing profile

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Nominal pressure
 *Does not include guidewire
 †Does not include guidewire