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


Designed for Durability

The Epic Stented Tissue Valve offers outstanding durability to help reduce the risk of reoperation. Identical in design to the Biocor Stented Tissue Valve, it features triple composite design, a unique pericardial shield and the FlexFit™ stent for ease of implant. The Epic Tissue Valve also includes Linx AC Technology, which is designed to improve long-term performance and valve durability.


Linx AC Technology

Tissue valve calcification is a common problem that impedes valve function. It is caused by a buildup of phospholipids, lipid deposits, and fibrous tissue. Calcification is a natural aging process. In animal studies, Linx AC Technology has been shown to reduce the risk of all these processes.


Triple Composite Design

Three separate porcine leaflets are matched to optimize leaflet coaptation and reduce stress.

Unique Pericardial Shield

The outflow edge of the Epic valve is covered with a pericardial shield, providing a tissue-to-tissue interface to help prevent the risk of abrasion.


FlexFit Stent

The FlexFit stent reduces leaflet stress.

Proven 20-year Durability

The Epic Valve is identical in design to the Biocor Valve, which delivers proven 20-year durability results in the aortic and mitral positions.

AVR: Study size = 1518 (female = 564, male = 964). Mean age = 70.8 years.
 MVR: Study size = 194 (female = 99, male = 95). Mean age = 64.9 years.

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Designed for Implantability

With the industry's lowest implant profile, the Epic Stented Tissue Valve with the FlexFit System offers optimal implantability in both the aortic and mitral positions.

Aortic Valves and Epic Supra

- Epic Supra provides a low profile, low ventricular protrusion, and low aortic protrusion.
- Epic Supra provides a low profile, low ventricular protrusion, and low aortic protrusion.
- Epic Supra provides a low profile, low ventricular protrusion, and low aortic protrusion.

Aortic Valve

9 mm ventricular protrusion

27 mm valve


Supra Valve

13 mm ventricular protrusion

23 mm valve

Protrusion Heights


24

Epic Mitral Valve

- Low profile reduces risk of LV outflow tract obstruction.
- Mitral retractor minimizes patient aortic annulus for suture looping.
- Lowest ventricular valve and holder height facilitates minimally-invasive surgical procedures.

Mitral Valve


9 mm ventricular protrusion

27 mm valve


Optimal Stent to Annulus Ratio: Epic Supra

The Epic Supra provides a larger stent-to-annulus ratio than the Epic Aortic valve.

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
Amas

The FlexFit System

The FlexFit System from St. Jude Medical works to help your patients by providing the most advanced mitral valve holder.

FlexFit Mitral Valve Holder

The FlexFit Mitral Valve Holder provides non-obstructive, low-profile, and atraumatic valve support for mitral valve surgery. It is designed to be easy to use.

FlexFit Aortic Valve Holder

The FlexFit Aortic Valve Holder improves visibility and cuff access for suturing.

FlexFit 20-second Rinse


Rinse time is reduced to 2 x 10 seconds.

FlexFit Click-in Handle

Both the aortic and mitral valve holders provide quick handle attachment and removal to stabilize the valve and minimize the risk of unseating.

FlexFit Valve Sizers and System Tray


FlexFit System sizers are optimized for valve selection, and the tray provides one convenient location for all valve accessories.

St. Jude Medical is committed to putting more control into the hands of physicians and thereby reducing risk. As part of a comprehensive portfolio that delivers on this promise, the Epic stented tissue valve with Linx AC Technology is designed to provide long-term durability and enhanced flexibility during implant to help reduce the risk of reoperation.

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Asst.
[Signature]

LONG-TERM DURABILITY BASED ON PROVEN DESIGN.

The Epic Stented Tissue Valve is identical in design to the Biocor™ Stented Tissue Valve, which has proven 20-year durability and more than 25 years of clinical experience. In addition, Epic Valve includes Linx AC Technology, a patented, proprietary valve treatment designed to improve long-term performance and valve durability.*

*There is no clinical data currently available that evaluates the long-term impact of anticalcification tissue treatment in humans.

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


Accessories

Model Number	Description
B1000	Bioprosthetic Heart Valve Sizer Set One (1) click-in holder handle Six (6) double-ended aortic/Supra sizes (19, 21, 23, 25, 27, 29) Five (5) mitral sizes (25, 27, 29, 31, 33) One (1) autoclavable tray for storage of components One (1) autoclavable tray cover
B1000-H	Replacement Handle One (1) click-in holder handle

3.5
4.3
3.9

Recommended Storage Temperature: 5° - 25° C (41° - 77° F)
Storage Solution: 0.5% formaldehyde solution.
Resterilization: Epic valves must not be resterilized.
Shelf Life: Four years.
Radiopaque Marker for X-ray: A stainless steel wire located under the sewing cuff allows the valve to be seen under X-ray.

1. Mykén PS, Bøch-Hansen O. A 20-year experience of 1712 patients with the Biocor porcine bioprosthesis. *J Thorac and Cardiovasc Surg.* 2009;137(1):76-81.
2. St. Jude Medical Instructions for Use: Stented Porcine Tissue Valves Epic/Epic Supra, 681936-001.
3. St. Jude Medical, Inc. Engineering Report 05-04-01.
4. Mosaic Marketing Materials, UC200103933b EN (2002). Hancock II Marketing Materials, UC200001201 EN (1999).
5. Carpentier-Edwards PERIMOUNT Magna Pericardial Bioprosthesis model 7000 TFX mitral Instructions for Use.
6. Reference Dimensions - St. Jude Medical Instructions for Use: Stented Porcine Tissue Valves Epic/Epic Supra, 681925-001. Internal diameter measurements from St. Jude Medical, Inc. Part number 300211.

* This size not available with all product enhancements. Please speak with a St. Jude Medical representative for further detail.

ATRIAL FIBRILLATION CARDIAC RHYTHM MANAGEMENT CARDIOVASCULAR NEUROMODULATION

Global Headquarters
One St. Jude Medical Drive
St. Paul, Minnesota 55117
USA
+1 651 756 2000
+1 651 756 3310 Fax

Cardiovascular Division
177 East County Road B
St. Paul, Minnesota 55117
USA
+1 651 490 4470
+1 651 490 4466 Fax

**St. Jude Medical
Coordination Center BVBA**
The Corporate Village
Avenue Da Vinci laan 11 Box F1
B-1935 Zaventem
Belgium
+32 2 774 68 11
+32 2 772 83 84 Fax

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St. Jude Medical (Hong Kong) Ltd.
Suite 1608, 16/F Exchange Tower
33 Wang Chiu Road
Kowloon Bay, Kowloon
Hong Kong SAR
+852 2956 7688
+852 2956 0622 Fax

St. Jude Medical Japan Co., Ltd.
3-1-30, Minami-Aoyama
Minato-ku
Tokyo 107 0052
Japan
+81 3 3423 6450
+81 3 3402 5586 Fax

ST. JUDE MEDICAL
MORE CONTROL. LESS RISK.

Biocor, Epic, FlexFit, Linx, ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2009 St. Jude Medical. All rights reserved.
Carpentier-Edwards PERIMOUNT is a registered trademark of Edwards Lifesciences Corporation. Mosaic and Hancock are registered trademarks of Medtronic, Inc. Mitrolflow Synergy is a registered trademark of St. Jude Medical Corp. **คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์**
Item 100015876

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
ST. JUDE MEDICAL
INSTRUCTIONS FOR USE
Stented Porcine Tissue Valves
Epic™ / Epic™ Supra

INDICATIONS FOR USE
 The Epic valve is indicated for patients requiring replacement of a diseased, damaged, or malfunctioning native aortic or mitral valve. It may also be used as a replacement for a previously implanted aortic and/or mitral prosthetic heart valve.

The Epic Supra valve is indicated for patients requiring replacement of a diseased, damaged, or malfunctioning native aortic heart valve. It may also be used as a replacement for a previously implanted aortic prosthetic heart valve.

CONTRAINDICATIONS
 None known.

WARNINGS AND PRECAUTIONS
Warnings

- Valve size selection is based on the size of the recipient annulus, and for supra-annular aortic placement, the size of the sinusoidal space. Implantation of an inappropriately large bioprosthesis may result in stent deformation, valvular incompetence, and/or damage to the surrounding tissue. The use of an inappropriately small bioprosthesis may result in suboptimal hemodynamics.
- Use only the St. Jude Medical Bioprosthetic Heart Valve Size Set Model B1600 with the Epic and Epic Supra valves.
- Accelerated deterioration due to calcific degeneration of the Epic and Epic Supra valves may occur in:
 - children, adolescents, or young adults;
 - patients with altered calcium metabolism (e.g., patients with hyperparathyroidism or chronic renal failure); or
 - individuals requiring hemodialysis.
- For single use only.
- Do not re-sterilize the valve by any method.
- Presence of a catheter or transvenous pacing lead through any bioprosthetic may damage the valve and is therefore not recommended.

Do not use if:

- the valve has been dropped, damaged, or contaminated in any way, or if there is any sign of deterioration;
- the expiration date has elapsed;
- the tamper-evident container seal is damaged, broken, or missing, or if fluid is leaking from the packaging; or
- the storage solution does not completely cover the valve.

Precautions

- The safety and effectiveness of the Epic and Epic Supra valves has not been established for the following specific populations:
 - patients who are pregnant
 - patients with chronic renal failure
 - patients with end-stage renal failure
 - patients with end-stage renal failure
 - patients with chronic kidney disease
 - patients with chronic renal insufficiency
 - patients requiring primary or secondary valve replacement
 - children, adolescents, or young adults
- Sizes are supplied non-sterile, and must be cleaned and sterilized prior to each use. Do not use cracked, deformed, or damaged size set components.
- Do not pass the flanged portion of the valve replica string into the annulus.
- Do not place the non-sterile exterior of the valve container in the sterile field.
- Do not expose the valve to solutions other than the formaldehyde valve storage solution in which it was shipped, the sterile isotonic saline solution used during the rinsing procedure, or the sterile isotonic saline solution used to irrigate the valve.
- Do not apply antibiotics to either the formaldehyde valve storage solution or the rinse solution.
- Do not apply antibiotics to the valve.
- Do not allow the valve tissue to dry. Place the valve in sterile isotonic saline solution immediately upon removal from the valve storage solution. Once removed from this solution, the valve should be periodically irrigated during implantation.
- Do not use the valve if shipping temperature indicators on the product carton have turned red, or if the valve has been improperly stored in temperature conditions outside of the 5 °C to 25 °C range.
- Do not implant the valve without thoroughly rinsing as directed.
- Do not increase the valve tissue. If a valve is damaged, the valve must be explanted and replaced.
- Do not attempt to repair a valve. Damaged valves must not be used.
- Do not use cutting edge needles, unprotected forceps, or sharp instruments as they may cause structural damage to the valve.
- Never handle the leader tissue.
- Position the mitral valve in a manner to avoid commissure avulsion of the left ventricular outflow tract, and minimize any potential of commissure contact with the ventricular wall.
- Position the aortic valve so that the stent posts do not obstruct the coronary ostia.
- Avoid prolonged contact with the formaldehyde storage solution. Immediately after contact, thoroughly flush any skin exposed to the solution with water. In case of contact with eyes, flush with water and seek appropriate medical care.

St. Jude Medical, a division of Medtronic, Inc., 515 North Dearborn Street, St. Paul, MN 55102, USA. Epic and Epic Supra are registered trademarks of St. Jude Medical, a division of Medtronic, Inc. © 2006 St. Jude Medical, a division of Medtronic, Inc. All rights reserved.

The valve is packaged in a formaldehyde storage solution. Store the valve in the upright position. The valve should be stored in temperatures from 5 °C to 25 °C (41 °F to 77 °F). Do not store the valve when significant temperature fluctuations may occur.

CAUTION: Do not implant the valve without thoroughly rinsing as directed.

CAUTION: Do not use the valve if shipping temperature indicators on the product carton have turned red, or if the valve has been improperly stored in temperature conditions outside the 5 °C to 25 °C (41 °F to 77 °F) range.

DIRECTIONS FOR USE
Pre-Implant Handling

The Epic and the Epic Supra valves are supplied in a storage container with a screw cap closure and tamper-evident seal. The contents of the container are sterile, and must be handled aseptically to prevent contamination.

Warnings

- Do not use the valve if the expiration date has elapsed.
- Do not use the valve if fluid is leaking from the packaging.
- Do not re-sterilize the valve by any method.

Removing the Valve from the Outer Packaging
Precautions

- Do not place the non-sterile exterior of the valve container in the sterile field.
- Do not expose the valve to solutions other than the formaldehyde valve storage solution in which it was shipped, the sterile isotonic saline solution used during the rinsing procedure, or the sterile isotonic saline solution used to irrigate the valve.
- Do not add antibiotics to either the formaldehyde valve storage solution or the rinse solution.
- Do not apply antibiotics to the valve.
- After rinsing, observe a valve of the appropriate size.
- Once the valve container has been removed from the outer packaging, examine the container for evidence of damage.

WARNING: The valve must not be implanted if the tamper-evident container seal is damaged, broken, or missing, or if fluid is leaking from the packaging.

WARNING: The valve must not be implanted if the storage solution does not completely cover the valve.

Verify the valve size and expiration date on the label. To remove the valve from the container, break the seal and remove the screw-top closure.

CAUTION: Avoid prolonged contact with the formaldehyde storage solution. Immediately after contact, thoroughly flush any skin exposed to the solution with water. In case of contact with eyes, flush with water and seek appropriate medical care.

Complete the medical device registration form and return it to St. Jude Medical. Place one of the pull-out labels with the designated model and serial number in the patient's chart.

Removing the Valve from the Storage Container

- Select the valve holder handle model B1600 H.
- With the circulating nurse holding the container, press the valve holder handle into the valve holder as shown in Figure 4, and remove the valve from the container.

CAUTION: Do not use cutting edge needles, unprotected forceps, or sharp instruments as they may cause structural damage to the valve. Never handle the leader tissue.

Inspect the valve for damage.

WARNING: Do not implant the valve if it has been dropped, damaged, or mishandled in any way, or if there is any sign of deterioration.

Rinse Procedure

CAUTION: Do not implant the valve without thoroughly rinsing as directed.

- Within the sterile field, prepare two sterile basins with a minimum of 500 ml of sterile isotonic saline solution in each basin.
- Holding the valve by the handle, fully immerse the valve support, the valve, the valve holder, and the portion of the holder handle that was submerged in the valve storage solution in the sterile isotonic saline solution in the first basin.
- Continually rinse the valve for 10 seconds, using a gentle back-and-forth motion.
- Repeat steps two and three in the second basin.
- After rinsing, leave the valve immersed in the basin until required by the surgeon for implantation.

CAUTION: Do not allow the valve tissue to dry. Place the valve in sterile isotonic saline rinse solution immediately upon removal from the valve storage solution.

Sizing Epic Standard Aortic Valves

Epic standard aortic valves are designed for supra-annular stent placement and supra-annular cuff placement. To determine the correct standard aortic valve size, use the Model B1600 aortic sizer. The Model B1600 aortic sizer is a double-ended tool, with a supra-valve replica end and an annular sizing end (Figure 5).

NOTE: Use only the aortic sizing end of the B1600 sizer to size standard aortic valves.

Identify the size that fits snugly in the annulus and select the corresponding valve size.

WARNING: Standard aortic valve size selection is based on the size of the recipient annulus. Implantation of an inappropriately large bioprosthesis may result in stent deformation, valvular incompetence, and/or damage to the surrounding tissues. The use of an inappropriately small bioprosthesis may result in suboptimal hemodynamics. Use only the St. Jude Medical Model B1600 sizers to size Epic standard aortic valves.

CAUTION: Sizers are supplied non-sterile, and must be cleaned and sterilized prior to each use. Do not use cracked, deformed, or damaged size set components.

Sizing Epic Supra Aortic Valves

Epic Supra valves are designed for supra-annular aortic placement. To determine the correct supra-aortic valve size, use the Model B1600 aortic sizer. The aortic sizer is a double-ended tool, with a supra-valve replica end and an annular sizing end. Use the annular sizing end to determine the size of the annulus. Insert the corresponding supra-valve replica end in the supra-annular space to confirm placement and fit of the valve (Figure 6).

WARNING: Supra-annular aortic valve size selection

Storage Solution: Sterile Isotonic Saline Solution	Expiration Date: 24 Months
Dimensions: 45mm	45mm
Weight: 150g	150g
Material: Porcine Tissue	Temperature: 5°C to 25°C (41°F to 77°F)
Manufacturer: St. Jude Medical	Do Not Re-sterilize
Product Name: Epic Supra	Pre-sterilized by Ethylene Oxide
Device Number: 271	Temperature Indicator: Not Applicable
Lot Number: 1234	Manufacturer: St. Jude Medical
REF: 1003394156	Device Number: 271
Do Not Re-sterilize	Manufacturer: St. Jude Medical
Autoclaved Disinfection Temperature: 121°C	Manufacturer: St. Jude Medical

DEVICE DESCRIPTION

The Epic™ valve (Figure 1) and the Epic™ Supra valve (Figure 2) are high-pressure latex valves manufactured from stents, porcine aortic valve cusps. The cusps are carefully matched for optimum leaflet coaptation and hemodynamics. Following tissue fixation, the tissue is mounted on the Flex™ polyester covered flexible acetate copolymer stent. The stent is a low-profile design with a scalloped shape. The stent is covered with knitted polyester fabric. The sewing cuff on Epic standard valves is formed from a knitted polyester fabric material that is covered by the knitted polyester fabric cover. The sewing cuff on Epic Supra valves is formed by enclosing a molded silicone elastomer within the same knitted polyester fabric as the Epic standard valve. The sewing cuff also contains three suture markers to facilitate valve placement. A suture steel wire is located under the sewing cuff for radiopacity.

A flexible pericardial tissue strip is attached to the outflow edge of the valve. The strip protects the leaflets as they open and close. The pericardial strip and the porcine valve cusps are preserved and cross-linked in a glutaraldehyde solution. Glutaraldehyde, formaldehyde, and ethanol are used in the valve sterilization process. The valve is attached to a holder and packaged in a sealed bag containing a formaldehyde solution. The Epic and Epic Supra valves are supplied sterile and non-proprietary.

Non-clinical testing has demonstrated that the Epic and Epic Supra valves are MR Conditional. Refer to the MRI Safety Information section of this manual for further information.

The Epic valves are available in the aortic and mitral size indicated in Table 1. Epic standard valves are designed to allow intra-annular placement of the inflow edge of the valve with supra-annular placement of the sewing cuff.

Epic Supra valves are available in the aortic sizes indicated in Table 1. Epic Supra valves are designed for supra-annular implantation of both the valve and the sewing cuff. The sewing cuff is the only difference between the Epic valve and the Epic Supra valve.

The Epic and Epic Supra valves are treated with the Flex™ antimicrobial protection.

Table 1: Catalog Number, Descriptions and Reference Dimensions

Catalog Number	Inlet Annulus Diameter (mm)	Aortic Ventricular Prosthesis (mm)	Total Height (mm)	Cuff Outer Diameter (mm)
Aortic Heart Valves, Epic				
E100-25A-00	21	9	14	25
E100-25A-01	21	9	15	27
E100-25A-02	25	16	16	26
E100-27A-03	27	11	17	31
E100-29A-06	29	12	19	33
Aortic Heart Valves, Epic Supra				
ES100-19-00	19	11	14	25
ES100-21-00	21	11	15	28
ES100-23-00	23	13	16	29
ES100-25-00	25	13	17	31
ES100-27-00	27	14	18	33
Mitral Heart Valves, Epic				
M100-25M-00	25	9	16	33
M100-27M-00	27	9	17	35
M100-29M-00	29	16	19	37
M100-31M-00	31	16	20	39
M100-33M-00	33	11	20	41

Table 1: Catalog Number, Descriptions and Reference Dimensions

Table 1 presents the number of patients implanted, cumulative follow-up, and late follow-up for each implant position.

Preoperative Patient Demographics
 Table 4 presents the preoperative patient demographics.

Efficacy Endpoints
 Quantitative data were collected throughout the study (i.e., NYHA functional classification, etc. parameters). Table 5 presents patient NYHA classification at one year follow-up. Table 6 and 7 present the hemodynamic follow-up results for the Epic aortic and mitral valve replacements.

Table 8 presents the hemodynamic follow-up results for the 25mm mitral Epic valve size augmented with data collected from Durofil for the Bicor 25mm valve.

INDIVIDUALIZATION OF TREATMENT
Anticoagulant and/or Antiplatelet Therapy
 Long-term, low-dose aspirin, unless contraindicated, is recommended for all patients with bioprosthetic valves. Long-term anticoagulant therapy, unless contraindicated, is recommended for all patients with bioprosthetic valves who have risk factors for thromboembolism.

PATIENT COUNSELING INFORMATION
 Long-term, low-dose aspirin, unless contraindicated, is recommended for all patients with bioprosthetic valves. Long-term anticoagulant therapy, unless contraindicated, is recommended for all patients with bioprosthetic valves who have risk factors for thromboembolism.

Patients with bioprostheses who undergo dental or other procedures that use potentially bactericidal should receive endocarditis prophylaxis as indicated.

St. Jude Medical publishes a patient brochure. Copies of this brochure are available through your SIM sales representative.

PACKAGING AND STORAGE
 Each Epic and Epic Supra valve is individually packaged in three layers of polyethylene. The Epic and Epic Supra valves are stored in the valve holder, and support facilitate handling and manipulation of the valve during removal from the container, rinsing, and implantation.

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is based on the size of the recipient annulus and the anatomy of the subannular space. Implantation of an inappropriately large bioprosthesis may result in stent deformation, valvular incompetence, and/or damage to the surrounding tissues. The use of an inappropriately small bioprosthesis may result in suboptimal hemodynamics. Use only St. Jude Medical Sizer Set Model B1000 to size Epic Supra valves.

CAUTION: Sizers are supplied non-sterile, and must be cleaned and sterilized prior to each use. Do not use cracked, deformed, or damaged sizer set components.

CAUTION: Do not pass the flanged portion of the supra valve replica sizing tool through the annulus.

Sizing Epic Mitral Valves

Epic mitral valves are designed for infra-annular stent placement and supra-annular cuff placement.

To determine the correct standard mitral valve size, use the Model B1030 mitral sizer (Figure 7). Identify the sizer that fits snugly in the annulus and select the corresponding valve size.

WARNING: Mitral valve size selection is based on the size of the recipient annulus. Implantation of an inappropriately large bioprosthesis may result in stent deformation, valvular incompetence, and/or damage to the surrounding tissues. The use of an inappropriately small bioprosthesis may result in suboptimal hemodynamics. Use only St. Jude Medical Sizer Set Model B1000 to size Epic mitral valves.

CAUTION: Sizers are supplied non-sterile, and must be cleaned and sterilized prior to each use. Do not use cracked, deformed, or damaged sizer set components.

Surgical Guidelines

The actual choice of surgical technique, modified in accordance with the instructions described herein, is left to the discretion of the individual surgeon.

When implanting supra-annular valves, non-everting mattress sutures are recommended.

Avoid any contact between the implantation raters and the leads.

Precautions:

- Do not allow the valve tissue to dry. Place the valve in isotonic sterile saline rinse solution immediately upon removal from the valve storage solution. Once removed from this solution, the valve should be periodically irrigated during implantation.
- Do not use cutting edge needles, unprotected forceps, or sharp instruments as they may cause structural damage to the valve.
- Never handle the leaflet tissue.
- Do not lacerate the valve tissue. If a valve is damaged, the valve must be explanted and replaced.
- Do not attempt to repair a valve. Damaged valves must not be used.

Aortic Valve Implantation

The Epic aortic valve is designed to allow infra-annular placement of the inflow edge of the valve with supra-annular placement of the sewing cuff. The Epic Supra valve is designed for supra-annular implantation of both the valve and the sewing cuff.

1. After sizing, choose a valve of the appropriate size.

CAUTION: Position the valve so that the stent posts do not obstruct the coronary ostia.

2. To facilitate implantation, the valve holder handle may be removed from the valve holder by depressing the release button on the valve holder (Figure 8).

3. To remove the valve holder from the valve, cut the three retaining sutures as shown in Figure 9, and pull the valve holder away from the valve. Examine the valve to ensure that there are no holder suture remnants.

Mitral Valve Implantation

1. After sizing, choose a valve of the appropriate size.

CAUTION: Position the mitral valve in a manner to avoid commissure obstruction of the left ventricular outflow tract, and minimize any potential of commissure contact with the semilunar wall.

2. To facilitate insertion of the mitral valve into the annulus, the mitral valve stent posts may be temporarily deflected inward during implantation. To deflect the valve stent posts inward, rotate the valve holder handle in the clockwise direction (Figure 10) while securing the valve holder. The valve stent posts will remain deflected until the valve retaining sutures are cut. To facilitate implantation, the valve holder handle may be removed from the valve holder by depressing the release button on the valve holder.

NOTE: Take care to avoid looping or entangling sutures around the commissural posts, as this may result in a compromise of leaflet function.

3. To remove the valve holder from the valve, cut the three retaining sutures as shown in Figure 11, and pull the valve holder away from the valve. Examine the valve to ensure that there are no holder suture remnants.

INTRA-OPERATIVE ASSESSMENT

The suggested method for assessing competence of the Epic and the Epic Supra valve is with intra-operative Doppler echocardiography.

PATIENT REGISTRATION

A medical device registration form is included with each device. After implantation, please complete all requested information, and return the original form to the address indicated on the medical device registration form. Tracking by manufacturers is mandatory within the United States and Canada. For countries outside of the United States and Canada, please disregard any request for patient information if this contradicts your local legal or regulatory requirements regarding patient privacy.

LIMITED WARRANTY

St. Jude Medical (SJM) warrants that reasonable care has been used in the manufacturing of this device. THIS WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES NOT EXPRESSLY SET FORTH HEREIN, WHETHER EXPRESS OR IMPLIED BY OPERATION OF LAW OR OTHERWISE INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, since handling, storage, cleaning, and sterilization of this device as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond SJM's control directly affect this device and the results obtained from its use. SJM SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE directly or indirectly arising from the use of this device other than the replacement of all or part of it. SJM neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.

Some states in the United States do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you. This limited warranty gives you specific legal rights, and you may have other rights which vary from jurisdiction to jurisdiction.

Descriptions of specifications, appearing in SJM literature, are meant solely to generally describe the device at the time of manufacture and do not constitute any express warranties.

Edwards LH, Clark EE, Cohn LH, Goehriester CL, Miller CM, White RD. Guidelines for Separating Mediastinal and Marital after Cardiac Valvular Operations. Ann Thorac Surg 1994;62:202-5


Figure 1: Epic Aortic and Epic Mitral Valves


Figure 2: Epic Supra Aortic Valve


Figure 3: Place the valve in the valve holder.

To release the valve support ring, the valve cuff depresses the stent posts at the base of the valve support, follow the support ring.


Figure 4: Release the valve holder.

Place the supra valve replica end in the supra-annular space to confirm placement and fit of the Epic supra valve.


Figure 5: Mitral Sizer


Figure 8: Release the valve holder handle by depressing the button.


Figure 9: Cut three sutures to remove aortic holder.

Figure 10: Rotate the valve holder handle to deflect the mitral stent posts inward.


Figure 11: Cut three sutures to remove mitral holder.

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TABLES

Table 1: Observed Adverse Event Rates

All subjects entered into study; N=762, cumulative follow-up 7172.4 Late Patient-years

Adverse Event	Early Events n (%)	Late Events n (%) (p<0.05)	Dependent Proportion Ratio (a)	Freedom From Event 1 Year % (95% CI)
Myocardial Infarction	2 (0.3)	1 (0.1)	0.164	99.6% (98.6%, 99.6%)
Structural Deterioration	0 (0.0)	2 (0.3)	0.324	100.0% (100.0%, 100.0%)
Paravalvular Leak	2 (0.3)	11 (1.5)	1.263	98.2% (96.7%, 99.6%)
Embolism	20 (2.6)	16 (2.2)	1.175	94.6% (92.8%, 96.3%)
Valve Thrombosis	1 (0.1)	0 (0.0)	0.014	99.5% (97.8%, 100.0%)
Major Bleeding, Pericardial Effusion and/or Anticoagulant-related (including Stroke)	76 (10.0)	13 (1.8)	1.377	93.2% (91.0%, 94.5%)
Myocardial Infarction, Myocardial Ischemia*	27 (3.5)	7 (1.0)	0.85	95.2% (93.2%, 96.4%)
Endocarditis	1 (0.1)	0 (0.0)	0.002	99.5% (97.1%, 99.2%)
Reoperation	1 (0.1)	11 (1.5)	1.457	98.3% (96.8%, 99.4%)
Mortality - Valve Related	2 (0.3)	5 (0.7)	0.834	99.2% (98.1%, 99.7%)

1. Late events are those occurring on or after 30 days post-implant. The early adverse event rate (a) is calculated as the ratio of early adverse events divided by the total number of subjects enrolled in the trial.
 2. Late events are those occurring 31 days post-implant and thereafter.
 3. Dependent proportion ratio is the ratio of the dependent event rate to the independent event rate.
 4. Freedom from adverse event is calculated using only completed therapy.

Table 3: Patient Numbers, and Cumulative and Late Patient Follow-Up

All subjects entered into study; N=762

Mean, SD, Min, and Max are represented in "Patient-years"

Implant Position	Number of Subjects	Total Patient-years	Mean	SD	Minimum	Maximum
Cumulative Patient-years						
Isolated Aortic	757	587.15	1.05	0.70	0.00	3.09
Isolated Mitral	176	165.76	0.94	0.73	0.00	3.10
Double	29	27.43	0.95	0.62	0.00	2.01
All Implants	762	772.34	1.02	0.72	0.00	3.10
Late Patient-years*						
Isolated Aortic	374	210.22	1.14	0.61	0.00	3.10
Isolated Mitral	157	156.12	1.00	0.65	0.00	3.01
Double	22	20.34	0.85	0.60	0.00	1.93
All Implants	644	712.48	1.11	0.62	0.00	3.01

* Late patient-years are determined from 31 days post-implant to the remaining event.

Table 4: Preoperative Patient Demographics

All subjects entered into study; N=762

Variable	Isolated Aortic (N=552)	Isolated Mitral (N=176)	Double (N=29)	All (N=762)
Age	74.4 ± 9.1 (24, 92)	72.1 ± 8.9 (44, 91)	75.9 ± 8.3 (55, 92)	73.9 ± 9.2 (24, 92)
Gender (Male)	61.6%	44.3%	34.5%	58.2%
Preoperative NYHA				
I	9.2%	8.0%	0.0%	8.5%
II	34.6%	30.1%	34.5%	33.6%
III	43.1%	41.2%	44.5%	42.8%
IV	12.4%	16.2%	11.6%	14.6%
Unknown	0.5%	6.0%	0.0%	0.5%

Table 5: Effectiveness Outcomes, NYHA Functional Classification: 1 year follow-up*

Subjects with both preoperative and 1 year NYHA measurements, N=606; n = number per subgroup

NYHA Class	Isolated Aortic n=358				Isolated Mitral n=173				Double n=14			
	Preoperative		1 Year		Preoperative		1 Year		Preoperative		1 Year	
	n	% (n/n)	n	% (n/n)	n	% (n/n)	n	% (n/n)	n	% (n/n)	n	% (n/n)
I	55	9.9	244	67.1	6	6.5	66	73.0	0	0	12	85.7
II	121	27.4	56	27.8	30	32.3	23	24.2	6	43.0	1	7.1
III	146	31.1	11	3.1	40	43.0	5	5.1	4	28.6	1	7.1
IV	28	6.1	9	2.5	17	18.1	0	0	4	28.6	0	0
All	151	18.6	153	18.6	83	10.0	93	11.0	14	10.0	14	10.0

* NYHA with both preoperative and one year NYHA measurements available are included in this table.

Table 6: Effectiveness Outcomes at One Year Follow-Up Visit, Hemodynamic Results - All Epic Aortic Valves

All aortic subjects entered into study; N=568

Hemodynamic Parameter	21mm	23mm	25mm	27mm	29mm*
Mean Gradient (mmHg)	n=49	n=120	n=121	n=36	n=10
- Mean ± SD	19.1 ± 8.2	13.9 ± 6.0	12.1 ± 5.1	11.4 ± 4.1	2.5 ± 3.3
- Min, Max	1.1, 43.5	1.2, 33.0	1.2, 34.3	0.5, 26.5	2.7, 12.7
EOA (cm²)	n=46	n=118	n=121	n=25	n=10
- Mean ± SD	10.4 ± 0.3	14.3 ± 0.5	15.2 ± 0.5	16.4 ± 0.4	24.1 ± 1.1
- Min, Max	0.5, 2.3	0.5, 1.0	0.1, 2.9	0.8, 2.7	1.2, 4.6
Regurgitation (n, %)	n=56	n=120	n=128	n=26	n=10
- None	47 (84.0%)	101 (79.9%)	92 (72.5%)	25 (100%)	9 (90.0%)
- Trivial	7 (12.5%)	17 (13.3%)	29 (22.5%)	0 (0%)	1 (10%)
- Mild	2 (3.6%)	6 (4.7%)	6 (4.7%)	1 (4%)	0 (0%)
- Moderate	0 (0%)	0 (0%)	1 (0.8%)	0 (0%)	0 (0%)
- Severe	0 (0%)	0 (0%)	0 (0%)	0 (0%)	0 (0%)
- Unknown	0 (0%)	0 (0%)	0 (0%)	0 (0%)	0 (0%)

* Number 2 subjects at 29mm were not included in this table.

Table 7: Effectiveness Outcomes at One Year Follow-Up Visit, Hemodynamic Results - All Epic Mitral Valves

All mitral subjects entered into study; N=297

Hemodynamic Parameter	21mm*	23mm	25mm	27mm*
Mean Gradient (mmHg)	n=30	n=41	n=25	n=24
- Mean ± SD	4.1 ± 2.0	5.5 ± 1.2	4.8 ± 1.3	4.1 ± 1.6
- Min, Max	2.7, 14.0	2.8, 10.0	2.6, 6.3	1.5, 7.9
EOA (cm²)	n=16	n=26	n=15	n=22
- Mean ± SD	14.1 ± 0.2	19.4 ± 0.5	16.9 ± 0.4	15.4 ± 0.3
- Min, Max	0.6, 2.1	0.6, 2.8	1.1, 2.4	1.1, 2.2
Regurgitation (n, %)	n=30	n=41	n=25	n=25
- None	27 (90.0%)	41 (100%)	21 (84.0%)	23 (92.0%)
- Trivial	0 (0%)	1 (2.4%)	0 (0%)	1 (4.0%)
- Mild	0 (0%)	1 (2.4%)	1 (4.0%)	0 (0%)
- Moderate	0 (0%)	0 (0%)	0 (0%)	1 (4.0%)
- Severe	0 (0%)	0 (0%)	0 (0%)	0 (0%)
- Unknown	0 (0%)	0 (0%)	0 (0%)	0 (0%)

* Number 2 subjects at 27mm were not included in this table.
 * Number of subjects included.

Table 8: Bicuspid/Epic Effectiveness Outcomes, 25mm Mitral Hemodynamic results (n=13 months); All mitral 25mm subjects; N=15


Bicuspid/Epic Study (n=13) Epic Valve (n=2)	
Hemodynamic Parameters	
Mean Gradient (mmHg)	n=15
- Mean ± SD	4.8 ± 2.5
- Min, Max	2.7, 12.9
EOA (cm²)	
- Mean ± SD	17.1 ± 0.3
- Min, Max	6.3, 1.7
Regurgitation (n, %)	
- None	11 (73.3%)
- Trivial	0 (0.0%)
- Mild	4 (26.7%)
- Moderate	0 (0.0%)
- Severe	0 (0.0%)
- Unknown	0 (0.0%)

n = number of subjects included.

Manufacturer:
St. Jude Medical
197 County Road E, East
St. Paul, Minnesota 55117 USA
+1 800 544 1464
+1 651 463 4410
Minnesota

Manufacturing Facility:
St. Jude Medical Curitiba Ltda.
Edifício 464
Calle O, Ave. 2
Zona Franca Croyal
E. Croyal, Alajurá
Curitiba, Brazil

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Bicuspid, Epic, Hefix, Hefix ST, MDR MEDICAL, the nine squares symbol and MDRP CONTROL LESS EPIC are trademarks and service marks of St. Jude Medical, Inc. and its related companies.

Brazilian patent 8402133-0; US patent 5,246,275

Copyright © 2010 St. Jude Medical, Inc. All rights reserved.

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ASMA

By Royal Charter

EC Design-Examination Certificate

Directive 93/42/EEC on Medical Devices, Annex II Section 4

No. CE 617865
Issued To: **St. Jude Medical**
177 County Road B East
St Paul
Minnesota
55117
USA

A.2

In respect of:
Biocor™, Epic™ and Epic™ Supra Heart Valves

BSI has performed a design examination on the above devices in accordance with the Council Directive 93/42/EEC Annex II Section 4 and Regulation 722/2012. The design conforms to the requirements of this directive and regulation. For marketing of these products an additional Annex II excluding Section 4 certificate is required.

For and on behalf of BSI, a Notified Body for the above Directive (Notified Body Number 2797):

Gary E Slack, Senior Vice President Medical Devices

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Expiry Date: 2024-05-26

...making excellence a habit.™

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.
 This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Bay Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780
 BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.
 A member of BSI Group of Companies.

EC Design-Examination Certificate

Supplementary Information to CE 617865

Issued To: **St. Jude Medical**
177 County Road B East
St Paul
Minnesota
55117
USA

Biocor Valve

Catalogue Number	Device Name	Model Type	Intended Purpose	Classification
B100-21A	Biocor Valve	Aortic Model Number	Biocor™ and Epic™ valves are indicated for patients requiring replacement of a diseased, damaged, or malfunctioning native aortic and/or mitral heart valve. Biocor and Epic valves may also be used as replacements for previously implanted aortic and/or mitral prosthetic heart valves.	Class III
B100-23A	Biocor Valve	Aortic Model Number		Class III
B100-25A	Biocor Valve	Aortic Model Number		Class III
B100-27A	Biocor Valve	Aortic Model Number		Class III
B100-29A	Biocor Valve	Aortic Model Number		Class III
B100-25M	Biocor Valve	Mitral Model Number		Class III
B100-27M	Biocor Valve	Mitral Model Number		Class III
B100-29M	Biocor Valve	Mitral Model Number		Class III
B100-31M	Biocor Valve	Mitral Model Number		Class III
B100-33M	Biocor Valve	Mitral Model Number		Class III

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...making excellence a habit.™

Page 2 of 6

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Bay Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780

BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.

A member of BSI Group of Companies.

Handwritten signature

Handwritten signature

EC Design-Examination Certificate

Supplementary Information to CE 617865

Issued To: **St. Jude Medical**
177 County Road B East
St Paul
Minnesota
55117
USA

Epic Valve

Catalogue Number	Device Name	Model Type	Intended Purpose	Classification
E100-21A	Epic Valve	Aortic Model Number	Biocor™ and Epic™ valves are indicated for patients requiring replacement of a diseased, damaged, or malfunctioning native aortic and/or mitral heart valve. Biocor and Epic valves may also be used as replacements for previously implanted aortic and/or mitral prosthetic heart valves.	Class III
E100-23A	Epic Valve	Aortic Model Number		Class III
E100-25A	Epic Valve	Aortic Model Number		Class III
E100-27A	Epic Valve	Aortic Model Number		Class III
E100-29A	Epic Valve	Aortic Model Number		Class III
E100-25M	Epic Valve	Mitral Model Number		Class III
E100-27M	Epic Valve	Mitral Model Number		Class III
E100-29M	Epic Valve	Mitral Model Number		Class III
E100-31M	Epic Valve	Mitral Model Number		Class III
E100-33M	Epic Valve	Mitral Model Number		Class III

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Expiry Date: **2024-05-26**

...making excellence a habit.™

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

Handwritten signature

Handwritten signature

Date 2nd August 2021

**Letter of Authorization of Distributor
Thailand**

To Whom It May Concern:

This is to certify that St. Jude Medical, LLC ("St. Jude Medical"), having a location at One St. Jude Medical Drive, St. Paul, Minnesota, 55117, USA, authorizes its affiliate, Abbott Medical (Thailand) Co., Ltd (previously known as "St. Jude Medical (Thailand) Co. Ltd."), located at Rasa Tower II, 17th Floor, Unit 1704, 555 Phaholyothin Road, Kwang Chatuchak, Chatuchak, Bangkok 10900, in Thailand, to distribute products of St. Jude Medical in Thailand until further notice.

As such, Abbott Medical (Thailand) Co., Ltd, is authorized to undertake the necessary actions to participate in all auctions and tenders regarding the supply of medical equipment and products from any of the divisions of St. Jude Medical identified above, and, in connection with such activity, to conclude contracts for the supply of such products, to carry out deliveries according to the conditions stipulated in the contracts and tender offers, to appoint sub-distributors to distribute such products in Thailand, and also to carry out any maintenance service, training of staff and to fulfill guarantee obligation, as needed and appropriate.

St. Jude Medical's authorization is valid for all the product lines of the St. Jude Medical Divisions until St. Jude Medical provides further written notification confirming that the services of Abbott Medical (Thailand) Co., Ltd. for the above matters are no longer needed.


Yours sincerely,

Name: John McCoy Jr.
Position: Vice President and Treasurer

STATE OF ILLINOIS
COUNTY OF LAKE

Signed and attested before me on 2nd August, 2021 by John McCoy Jr.

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Beneficiary Name : Abbott Medical (Thailand) Co., Ltd.
Incorporated in U.S.A with Limited Liability