Elecsys total P1NP

REF	REF	Σ	SYSTEM
03141071190	03141071500	100	cobas e 411 cobas e 601 cobas e 602

English

System information

For cobas e 411 analyzer: test number 190
For cobas e 601 and cobas e 602 analyzers: Application Code Number 044

Intended use

Immunoassay for the in vitro quantitative determination of total P1NP in human serum and plasma.
This assay is intended for use in monitoring therapy following the diagnosis of osteoporosis^{1,2,3} in post-menopausal women and in patients diagnosed with Paget's disease of the bone.^{4,5,6}
The electrochemiluminescence immunoassay "ECLIA" is intended for use on Elecsys and cobas e immunoassay analyzers.

Summary

Type I collagen is an important component of bone matrix and osteoblasts secrete its precursor procollagen molecule during bone formation.⁷
Type 1 procollagen contains both N-(amino) and C-(carboxy) terminal extensions. These extensions (propeptides) are cleaved by enzymes during bone matrix formation and released into the circulation. The propeptide measured by Elecsys total P1NP assay is from the amino terminal, hence P1NP procollagen-type 1 N-terminal-propeptide. P1NP is released during type 1 collagen formation and its subsequent incorporation into the bone matrix, and thus may be defined as a true bone formation marker.^{8,9} P1NP appears to be released as a trimeric structure (derived from the trimeric collagen structure) but is rapidly broken down to a monomeric form by thermal degradation effects.^{10,11} This Elecsys P1NP assay detects both fractions present in blood and is therefore called total P1NP.

Circulating levels of total P1NP demonstrate a significant change upon anti-resorptive^{12,13} as well as anabolic¹³ therapy within few months from treatment start. Suboptimal response to treatment may indicate non-compliance or the presence of secondary causes of osteoporosis which may need addressing.

Serum P1NP has been selected by the IOF-IFCC Bone Marker Standards Working Group as marker for bone formation, mainly based on the following criteria:

- it has been evaluated both for fracture prediction and monitoring osteoporosis therapies
- the assay is widely available, suitable for serum or plasma samples, with well documented sample handling and stability.⁷

Test principle

Sandwich principle. Total duration of assay: 18 minutes.

- 1st incubation: 20 µL of sample and a biotinylated monoclonal P1NP-specific antibody are incubated together.
- 2nd incubation: After addition of streptavidin labeled microparticles and a monoclonal P1NP-specific antibody labeled with a ruthenium complex¹, a sandwich complex is formed which becomes bound to the solid phase via interaction of biotin and streptavidin.
- The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell/ProCell M. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier.
- Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the reagent barcode or e-barcode.

a) Tris(2,2'-bipyridyl)ruthenium(II)-complex (Ru(bpy)₃²⁺)

Reagents - working solutions

The reagent rackpack is labeled as TP1NP.

- M Streptavidin-coated microparticles (transparent cap), 1 bottle, 6.5 mL: Streptavidin-coated microparticles 0.72 mg/mL; preservative.
- R1 Anti-P1NP-Ab-biotin (gray cap), 1 bottle, 10 mL: Biotinylated monoclonal anti-P1NP antibody (mouse) 2.5 mg/L; phosphate buffer 100 mmol/L, pH 7.2; preservative.
- R2 Anti-P1NP-Ab-Ru(bpy)₃²⁺ (black cap), 1 bottle, 8 mL: Monoclonal anti-P1NP antibody (mouse) labeled with ruthenium complex 2.5 mg/L; phosphate buffer 100 mmol/L, pH 7.2; preservative.

Precautions and warnings

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious or microbial waste:
Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:
Apply all relevant local disposal regulations to determine the safe disposal. Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

Warning

H317 May cause an allergic skin reaction.

Prevention:

- P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
- P272 Contaminated work clothing should not be allowed out of the workplace.
- P280 Wear protective gloves.

Response:

- P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.
- P362 + P364 Take off contaminated clothing and wash it before reuse.

Disposal:

- P501 Dispose of contents/container to an approved waste disposal plant.

Product safety labeling follows EU GHS guidance.

Contact phone: all countries: +49-621-7590

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

Reagent handling

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

All information required for correct operation is read in from the respective reagent barcodes.

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Elecsys total P1NP

Do not freeze.

Store the Elecsys reagent kit upright in order to ensure complete availability of the microparticles during automatic mixing prior to use.

Stability:	
unopened at 2-8 °C	up to the stated expiration date
after opening at 2-8 °C	8 weeks
on the analyzers	8 weeks

Specimen collection and preparation

Only the specimens listed below were tested and found acceptable. Serum collected using standard sampling tubes or tubes containing separating gel.

Li-heparin, K₂EDTA and K₃EDTA plasma.

Li-heparin and K₂EDTA plasma tubes containing separating gel can be used.

Criterion: Method comparison serum versus plasma. Slope 0.9-1.1 + intercept within ± 5 ng/mL + coefficient of correlation ≥ 0.95 .

Stable for 24 hours at 15-25 °C, 5 days at 2-8 °C, 6 months at -20 °C (± 5 °C). Samples may be frozen and thawed up to 5 times without adverse effects.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay. Do not use samples and controls stabilized with azide.

Ensure the samples, calibrators and controls are at 20-25 °C prior to measurement.

Due to possible evaporation effects, samples, calibrators and controls on the analyzers should be analyzed/measured within 2 hours.

Note: Avoid hemolysis! Samples showing visible signs of hemolysis may cause interference.

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

- [REF] 03141090190, total P1NP CalSet, for 4 x 1.0 mL
- [REF] 05618860190, PreciControl Varia, for 4 x 3.0 mL
- [REF] 11732277122, Diluent Universal, 2 x 16 mL sample diluent or [REF] 03183971122, Diluent Universal, 2 x 36 mL sample diluent
- General laboratory equipment

cobas e analyzer

Additional materials for the cobas e 411 analyzer:

- [REF] 11662988122, ProCell, 6 x 380 mL system buffer
- [REF] 11662970122, CleanCell, 6 x 380 mL measuring cell cleaning solution
- [REF] 11930346122, Elecsys SysWash, 1 x 500 mL washwater additive
- [REF] 11933159001, Adapter for SysClean
- [REF] 11706802001, AssayCup, 60 x 60 reaction cups
- [REF] 11706799001, AssayTip, 30 x 120 pipette tips
- [REF] 11800507001, Clean-Liner

Additional materials for cobas e 601 and cobas e 602 analyzers:

- [REF] 04880340190, ProCell M, 2 x 2 L system buffer
- [REF] 04880293190, CleanCell M, 2 x 2 L measuring cell cleaning solution
- [REF] 03023141001, PC/CC-Cups, 12 cups to prewarm ProCell M and CleanCell M before use

- [REF] 03005712190, ProbeWash M, 12 x 70 mL cleaning solution for run finalization and rinsing during reagent change

- [REF] 12102137001, AssayTip/AssayCup, 48 magazines x 84 reaction cups or pipette tips, waste bags

- [REF] 03023150001, WasteLiner, waste bags

- [REF] 03027651001, SysClean Adapter M

Additional materials for all analyzers:

- [REF] 11298500316, ISE Cleaning Solution/Elecsys SysClean, 5 x 100 mL system cleaning solution

Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

Resuspension of the microparticles takes place automatically prior to use. Read in the test-specific parameters via the reagent barcode. If in exceptional cases the barcode cannot be read, enter the 15-digit sequence of numbers.

Bring the cooled reagents to approximately 20 °C and place on the reagent disk (20 °C) of the analyzer. Avoid foam formation. The system automatically regulates the temperature of the reagents and the opening/closing of the bottles.

Calibration

Traceability: This method has been standardized against reference standards precisely defined by weighing native P1NP into an analyte-free human serum matrix.

Every Elecsys reagent set has a barcoded label containing specific information for calibration of the particular reagent lot. The predefined master curve is adapted to the analyzer using the relevant CalSet.

Calibration interval may be extended based on acceptable verification of calibration by the laboratory.

Renewed calibration is recommended as follows:

- after 1 month (28 days) when using the same reagent lot
- after 7 days (when using the same reagent kit on the analyzer)
- as required: e.g. quality control findings outside the defined limits

Quality control

For quality control, use PreciControl Varia.

In addition, other suitable control material can be used.

Controls for the various concentration ranges should be run individually at least once every 24 hours when the test is in use, once per reagent kit, and following each calibration.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

If necessary, repeat the measurement of the samples concerned.

Follow the applicable government regulations and local guidelines for quality control.

Calculation

The analyzer automatically calculates the analyte concentration of each sample (either in $\mu\text{g/L}$ or ng/mL).

Limitations - Interference

The effect of the following endogenous substances and pharmaceutical compounds on assay performance was tested. Interferences were tested up to the listed concentrations and no impact on results was observed:

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Elecsys total P1NP

Compound	Concentration tested
Rheumatoid factors	≤ 1000 IU/mL

Criterion: Recovery of ± 1.5 ng/mL of initial value for samples ≤ 15 ng/mL and within ± 10 % of initial value for samples > 15 ng/mL.

Samples should not be taken from patients receiving therapy with high biotin doses (i.e. > 5 mg/day) until at least 8 hours following the last biotin administration.

There is no high-dose hook effect at P1NP concentrations up to 3900 µg/L (ng/mL).

In vitro tests were performed on 28 commonly used pharmaceuticals. No interference with the assay was found.

The pharmaceuticals tested included drugs common to the management of osteoporosis, i.e. ibandronate, Risedronate and Alendronate as well as oestrogen-based therapies ("HRT") and also calcium and vitamin D supplementation therapies; none showed any evidence of assay interference.

In rare cases, interference due to extremely high titers of antibodies to analyte-specific antibodies, streptavidin or ruthenium can occur. These effects are minimized by suitable test design.

Bone metabolism may be affected by the use of cytotoxic agents. Results obtained from patients treated with such therapies should be interpreted with caution.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Limits and ranges

Measuring range

5-1200 µg/L or ng/mL (defined by the lower detection limit and the maximum of the master curve). Values below the lower detection limit are reported as < 5 µg/L (ng/mL). Values above the measuring range are reported as > 1200 µg/L (ng/mL).

Lower limits of measurement

Lower detection limit of the test

Lower detection limit: < 5 ng/mL (µg/L)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying two standard deviations above that of the lowest standard (master calibrator, standard 1 + 2 SD, repeatability study, n = 21).

Dilution

Samples with P1NP concentrations above the measuring range can be diluted with Diluent Universal. The recommended dilution is 1:2 (either automatically by the analyzers or manually). The concentration of the diluted sample must be > 100 µg/L (ng/mL).

After manual dilution, multiply the result by the dilution factor.

After dilution by the analyzers, the software automatically takes the dilution into account when calculating the sample concentration.

Non-linear dilution behaviour may be seen when using sera from patients diagnosed with renal insufficiency.

Expected values


Sera taken from 573 healthy female volunteers who had been enrolled in a study of determinants of bone loss (OFELY^{14,15}) were measured for total P1NP levels. The following results were obtained (µg/L or ng/mL):¹⁶

	Post-menopausal			Pre-menopausal
	All	HRT ^{b)} yes	HRT no	All
N	444	154	290	129
5 th percentile	16.27	14.28	20.25	15.13
Median	37.09	28.48	42.94	27.80
Mean	40.43	31.74	45.05	30.10
95 th percentile	73.87	58.92	76.31	58.59

b) HRT = patients receiving hormone replacement therapy


Below are frequency plots showing total P1NP concentration ranges with normal, untreated pre- versus untreated post-menopausal women (figure 1) and below (figure 2) P1NP ranges in post-menopausal women receiving HRT therapy versus those not receiving therapy.

Figure 1: Frequency of total P1NP concentrations (µg/L or ng/mL) observed in normal, untreated, pre- (n = 129) and post- (n = 290) menopausal women

x: total P1NP (µg/L or ng/mL) - - : pre-menopausal
y: Frequency (%) - - : post-menopausal

Figure 2: The effect of hormone replacement therapy on total P1NP (µg/L or ng/mL) concentration (frequency distribution in treated ("HRT yes"; n = 154) and untreated ("HRT no"; n = 290) post-menopausal women

x: total P1NP (µg/L or ng/mL) - - : HRT yes
y: Frequency (%) - - : HRT no

The measurement of total P1NP shows minimal circadian or seasonal variation (approx. 6 %) ^{17,18} and food intake or diet show no detectable influence upon serum levels. ^{19,20}

Significantly elevated serum total P1NP levels are associated with the presence of metastatic bone disease and may also be seen in patients with renal insufficiency. ^{21,22,23} Diseases associated with secondary bone disease may have an effect upon levels of total P1NP. ^{24,25}

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using Elecsys reagents, samples and controls in a protocol (EP5-A2) of the CLSI (Clinical and Laboratory Standards Institute): 2 runs per day in duplicate each for 21 days (n = 84). The following results were obtained:

Elecsys total P1NP

cobas e 411 analyzer					
Sample	Mean µg/L (ng/mL)	Repeatability		Intermediate precision	
		SD µg/L (ng/mL)	CV %	SD µg/L (ng/mL)	CV %
Human serum 1	12.8	0.340	2.7	0.531	4.1
Human serum 2	57.2	1.04	1.8	1.34	2.4
Human serum 3	527	6.93	1.3	11.7	2.2
Human serum 4	33.4	0.891	2.7	0.960	2.9
Human serum 5	1140	34.1	3.0	37.6	3.3
PreciControl Varia 1	30.5	0.488	1.6	0.629	2.1
PreciControl Varia 2	166	2.06	1.2	2.83	1.7

cobas e 601 and cobas e 602 analyzers					
Sample	Mean µg/L (ng/mL)	Repeatability		Intermediate precision	
		SD µg/L (ng/mL)	CV %	SD µg/L (ng/mL)	CV %
Human serum 1	14.4	0.275	1.9	0.527	3.7
Human serum 2	57.9	1.15	2.0	1.47	2.5
Human serum 3	496	8.42	1.7	11.3	2.3
Human serum 4	34.2	0.816	2.4	1.10	3.2
Human serum 5	1090	34.8	3.2	37.2	3.4
PreciControl Varia 1	30.7	0.524	1.7	0.797	2.6
PreciControl Varia 2	158	2.29	1.5	3.04	1.9

Method comparison

A comparison of the Elecsys total P1NP assay (y) with a commercially available P1NP test (x) using clinical samples gave the following correlations (µg/L or ng/mL):

Number of samples measured: 76

Passing/Bablok²⁶ Linear regression
 $y = 1.01x + 1.73$ $y = 1.02x - 5.55$
 $r = 0.879$ $r = 0.981$

The sample concentrations were between 19.7 and 319.5 µg/L (ng/mL).

Analytical specificity

Cross-reactivities of < 1% were seen with the following analytes: β-CrossLaps, N-MID Osteocalcin, parathyroid hormone (PTH), and 25-OH vitamin D.

References

- Fink E, Cormier C, Steinmetz P, et al. Differences in the capacity of several biochemical bone markers to assess high bone turnover in early menopause and response to alendronate therapy. *Osteoporos Int* 2000;11(4):295-303.
- Garnero P, Stevens RE, Ayres SA, et al. Short-term effects of new synthetic conjugated estrogens on biochemical markers of bone turnover. *J Clin Pharmacol* 2002;42(3):290-296.
- Sharp CA, Evans SF, Risteli L, et al. Effects of low- and conventional-dose transcutaneous HRT over 2 years on bone metabolism in younger and older post-menopausal women. *Eur J Clin Invest* 1996;26:763-771.
- Alvarez L, Guanabens N, Peris P, et al. Usefulness of biochemical markers of bone turnover in assessing response to treatment of Paget's disease. *Bone* 2001;29(5):447-452.
- Alvarez L, RicOs C, Peris P, et al. Components of biological variation of biochemical parameters in Paget's bone disease. *Bone* 2000;26(6):571-576.

- Reid IR, Davidson JS, Wattie D, et al. Comparative responses of bone turnover markers to bisphosphonate therapy in Paget's disease of bone. *Bone* 2004;35:224-230.
- Vasikaran S, Eastell R, Bruyère O, et al. Markers of bone turnover for the prediction of fracture risk and monitoring of osteoporosis treatment: a need for international reference standards. *Osteoporos Int* 2011;22:391-420
- Orum O, Hansen M, Jensen CH, et al. Procollagen type 1 N-terminal Propeptide (P1NP) as an indicator of type 1 collagen metabolism: ELISA development, reference interval, and hypovitaminosis D induced hyperparathyroidism. *Bone* 1996;19(2):157-163.
- Brandt J, Frederiksen JK, Jensen CH, et al. The N- and C-terminal propeptides of human procollagen type 1 (P1NP and P1CP): molecular heterogeneity and assay technology. Pgs. 73-81 In *Bone Markers Biochemical and Clinical perspectives*. Eds. Eastell R, Baumann M, Hoyle NR and Wiecezorek L. Dunitz, London 2001.
- Jensen CH, Hansen M, Brandt J, et al. Quantification of the N-terminal propeptide of human procollagen type 1 (P1NP): Comparison of ELISA and RIA with respect to different molecular forms. *Clin Chim Acta* 1998;269(1):31-41.
- Brandt J, Krogh TH, Jensen CH, et al. Thermal instability of the trimeric structure of the N-terminal propeptide of human procollagen type 1 in relation to assay technology. *Clin Chem* 1999;45(1):47-53.
- Grey A, Bolland M, Wattie D, et al. Prolonged antiresorptive activity of zoledronate: a randomized, controlled trial. *J Bone Miner Res* 2010;25(10):2251-5.
- Garnero P, Vergnaud P, Hoyle N. Evaluation of a Fully Automated Serum Assay for Total N-Terminal Propeptide of Type I Collagen in Postmenopausal Osteoporosis. *Clin Chem* 2008;54(1):188-96.
- Garnero P, Sornay-Rendu E, Chapuy MC, et al. Markers of bone turnover predict postmenopausal forearm bone loss over four years: The OFELY study. *J Bone Miner Res* 1999;14:1614-1621.
- Garnero P, Sornay-Rendu E, Clauserat B, et al. Biochemical markers of bone turnover, endogenous hormones and the risk of fractures in postmenopausal women: The OFELY study. *J Bone Miner Res* 2000;15:1526-1536
- Garnero P. Personal communication. Data on file at Roche Diagnostics GmbH, Jan 2004.
- Ahmad AM, Hopkins MT, Fraser WD, et al. Parathyroid hormone secretory pattern, circulating activity, and effect upon bone turnover in adult growth hormone deficiency. *Bone* 2003;32(2):170-179.
- Blumsohn A, Naylor KE, Timm W, et al. Absence of marked seasonal change in bone turnover: A longitudinal and multicentre cross-sectional study. *J Bone Miner Res* 2003;18(7):1274-1281.
- Clowes JA, Hannon RA, Yap TS, et al. Effect of feeding on bone turnover markers and its impact on biological variability of measurements. *Bone* 2002;30(6):886-890.
- Clowes JA, Robinson RT, Heller SR, et al. Acute changes of bone turnover and PTH induced by insulin and glucose: Euglycemic and hyperglycaemic hyperinsulinemic clamp studies. *J Clin Endocrinol Metab* 2002;87:3324-3329.
- Marin L, Koivula MK, Jukkola-Vuorinen A, et al. Comparison of total and intact aminoterminal propeptide of type I procollagen assays in patients with breast cancer with or without bone metastases. *Ann Clin Biochem* 2011;48(5):447-451.
- Oremek GH, Sauer-Eppel H and Klepzig M. Total procollagen type 1 amino-terminal propeptide (total P1NP) as a bone metastasis marker in gynecological carcinomas. *Anticancer Res* 2007;27(4a):1961-2.
- Kolzumi M, Takahashi S, and Ogata E. Comparison of serum bone resorption markers in the diagnosis of skeletal metastasis. *Anticancer Res* 2003;23(5b):4095-9.
- Eastell R. and Szulc P. Use of bone turnover markers in postmenopausal osteoporosis. *Lancet Diabetes Endocrinol* 2017;5(11):908-923.
- Cavalier E, Lukas P, Callisi A, et al. Aminoterminal propeptide of type I procollagen (P1NP) in chronic kidney disease patients: the assay matters. *Clin Chim Acta* 2013;425:117-8.

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Elecsys total P1NP

26 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets, the product information and the Method Sheets of all necessary components (if available in your country).

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard (for USA: see dialog.roche.com for definition of symbols used):

- CONTENT** Contents of kit
- SYSTEM** Analyzers/instruments on which reagents can be used
- REAGENT** Reagent
- CALIBRATOR** Calibrator
- Volume after reconstitution or mixing
- GTIN** Global Trade Item Number

COBAS, COBAS E, ELECSYS and PRECICONTROL are trademarks of Roche. INTRALIPID is a trademark of Fresenius Kabi AB.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

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
Convenience

cobas e packs for convenient handling

- ▶ cobas e reagent packs are liquid, ready-to-use and simplify handling
- ▶ All-in-one cassette packaging for each parameter combined with calibrators streamlines logistics
- ▶ Automatic pack opening/closing ensures longer onboard stability
- ▶ Programming-by-loading concept ensures seamless, precise data handling
- ▶ Data review screen provides quick result traceability

Enhanced data review screen


Performance

9 min applications for speed and quality

- ▶ Over 75 assays offer comprehensive coverage of 7 indication areas
- ▶ High onboard stability and shelf-life allow continuous availability of both routine and esoteric parameters
- ▶ Disposable tips and assay cups for reaction eliminate risk of carryover
- ▶ Integrated safety features ensure result integrity

PTH Precision Monitoring NEQAS 2006

Elecsys assay performance in external ring-trials: PTH; Source: <http://www.ukneqas.org.uk/>; July 2008


Confidence

ECL technology: The gold standard

- ▶ High analytical sensitivity enables large measuring ranges and low sample volume
Measuring ranges: Estradiol 4-4,300 pg/ml; Progesterone 0-6 ng/ml; Total PSA 0-100 ng/ml;
Sample Volume: β -hCG requires 10 μ l; 3rd Gen TSH requires 50 μ l
- ▶ Activation with voltage for a controlled reaction allows high precision to avoid unnecessary repeats and reruns
- ▶ Short incubation time enables the timely availability of innovative applications

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
CD-ROM, USB & printer ports enable efficient data management

cobas[®] touchscreen user interface is standardized to other cobas systems

Disposable tips and cups eliminate carryover

Measuring cell provides unmatched detection using ECL technology

Also available rack version with 75 sample positions

30 continuous random access sample positions with STAT priority

18 channels enable use of up to 18 parameters simultaneously

Easy removal of disposables via solid waste container

Consumables and liquid waste are easily accessible

Full workstation consolidation

- ▶ The cobas e 411 analyzer is designed to deliver accurate results on account of proven ECL technology
- ▶ It consolidates routine endocrinology with cardiac, cancer, bone and infectious disease markers

Serum work area integration

- ▶ Combined with the cobas c 311 analyzer, the system processes up to 95% of your serum work area workload
- ▶ Work area management solutions integrate sample and data management across systems

Seamless laboratory networking

- ▶ As part of the cobas modular platform, the cobas e 411 analyzer offers consistent results with other cobas e systems
- ▶ Universal reagent carriers and analytical components enable flexible network solutions

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Technical specifications

System	Fully automated, continuous loading and random access system for immunoassay analysis based on ECL technology. Available as a disk and rack system. System components: Analytical module including touchscreen PC	
Test throughput	Up to 86 tests/hr	
Reagent concept	cobas e reagent pack: Onboard reagent capacity: Evaporation protection:	Ready-to-use reagents 18 channels/reagent slots for up to 18 different assays Automatic closure of reagent lids
Programmable parameters	Max 60 assays definable via 2D-barcode including all necessary information (programming by loading)	
Sample types	Serum, Plasma, Urine and others	
Onboard sample capacity	Load/unload capacity: Rack types: STAT handling:	Disk: 30 samples Rack: 75 samples on 15 RD standard 5 position rack Routine, STAT, Control, Calibrator Disk: Any unoccupied position on the sample disk Rack: Dedicated STAT port on rack feeder
Sample container types	Primary tubes: Sample cup: Cups on tube	5-10 ml; 16x100, 16x75, 13x100, 13x75 mm 2.5 ml
Sample volume	10 to 50 µl per test, depending on assay protocol	
Sample probe	Disposable tips for individual sample to avoid carry over; Sample clot detection; Liquid level detection	
Sample barcode types	PDF417, Code 128, Codabar (NW 7), Interleaved 2 of 5, Code 39	
System interfaces	Touchscreen PC, RS 232 serial interface, bi-directional Standard PC ports (USB, ethernet, serial etc) for other communication devices	
Sample database	2,000 results for routine, STAT and controls	
Test methods	Pre-defined assay protocols (sandwich, competitive, titration)	
Calibration methods	2-point calibration per lot or per cobas e reagent pack	
QC methods	Individual QC; Cumulative QC; Preventive QC after calibration	
Rerun/repeat/reflex	Rerun, repeat and reflex supported	
Electrical requirements	Power requirements: Power consumption:	100-120 VAC 50/60 Hz single phase or 200-240 VAC 50/60 Hz single phase 1000 VA
Water/waste requirements	Water supply: Water quality: Liquid waste: Water consumption:	3l container ≤10µS/cm or ≥0.1 mega Ohm, bacteria-free 4l container with overflow detection or direct drain Approx. 3l for 250 tests/approx. 12 ml/cycle
Regulatory compliance	GS, CE, UL, C-UL, CB-report and certificate	
Operating conditions	Ambient temperature: Ambient humidity: Noise output:	8 to 32°C (64.4°F to 89.6°F) 20% to 80% 60 dbA (stand-by mode), 63 dbA (avg. during operation)
Physical dimensions	Width: Depth: Height: Weight:	120 to 170 cm (disk/rack) 73 to 95 cm (disk/rack) 56 cm (w/o PC unit) Approx. 170 kg (disk) and 210 kg (rack)

COBAS, COBAS E, and LIFE NEEDS ANSWERS are trademarks of Roche.

Other brand or product names are trademarks of their respective holders.

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© 2020, Roche Diagnostics