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
Dräger Technology for LIFE


Illumination is a critical factor during surgical procedures.

The new Dräger Polaris® 600 surgical light combines advanced technologies to support precision and efficiency. Intuitive controls adapt the light to each procedure and individual surgeon's need, so you can focus on what matters most – your patient.


UNSURPASSED SURGICAL LIGHTING PERFORMANCE

With the special optical system, the homogeneous light field can be adjusted incrementally to quickly provide the right focus and optimal illumination at the surgical site, helping to eliminate shadows and providing maximum performance.


The Polaris 600 is extensively tested with our unique computerized testing methods. The LED lighting, specially designed for the Polaris 600, ensures 50,000 hours of uninterrupted, high-quality performance.

EASY CONTROL PUTS EVERY PROCEDURE IN THE RIGHT LIGHT

Unique handle with Sterile Touch Control (STC) allows the surgeon to adjust the light intensity and field diameter easily and independently from non-sterile personnel right at the fingertips.

INDUSTRY LEADING WIRELESS TECHNOLOGY

Optional HD camera with wireless video transfer, automatically adjusts to the current light conditions and provides an unrestricted view of the surgical area.

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

INDUSTRY LEADING WIRELESS TECHNOLOGY

Optional HD camera with wireless video transfer, integrated in the light head, automatically adjusts to the current light conditions and provides an unrestricted view of the surgical area.


INTUITIVE LIGHT TOUCH CONTROL PANEL

Control panel on the light head allows for complete lighting control for the entire OR team.

POLARIS 600 WALL CONTROL

Adjusts all aspects of the color coded Polaris lights as well as the HD Camera

OR INTEGRATION



Polaris 600 is equipped with an integrated OR interface, which allows you to control multiple devices via a central system.

FLEXIBLE TO MEET YOUR NEEDS

The Polaris 600 is available in several different configurations depending on your surgical needs and requirements.

A FAMILY OF PRODUCTS

More than a single product, Dräger also offers a full range of Polaris lighting.

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

ENHANCED OPTIONS 2.1.6

An intuitive light touch control panel for complete lighting control that benefits the entire OR Team.

STERILE CONTROL 2.1.8.3

The unique handle with Sterile Touch Control (STC) allows the surgeon to adjust the light easily and independently from non-sterile personnel. Touch-sensitive technology adjusts the light field size and light intensity. All can be easily operated with surgical gloves.

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ADJUSTABLE COLOR TEMPERATURE

From warm white to cold white (3800, 4400, 5000 and 5600 Kelvin) you can switch between four modes to reduce fatigue or glare during longer procedures and increase individual perception for each surgical situation.

2.1.8.13

3800 K

4400 K

5000 K

5600 K

EASE OF MOVEMENT

Movement of the Polaris 600 is simple and fluid. It easily moves to where it's needed allowing you to focus on what's important.

STREAMLINED DESIGN

Flat and easy to clean surfaces and laminar flow compatible design.

2.1.5

WIRELESS HD CAMERA

The optional HD camera with wireless video transfer can be added at any time. It automatically adjusts to the current light conditions.

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
91 00 371 16 02 11 Marketing Communications | AX | Printed in the USA | Client: (line) | Subject to Modifications | © 2016 Draeger USA

CORPORATE HEADQUARTERS
 Drägerwerk AG & Co. KGaA
 Moislinger Allee 53-55
 28558 Lübeck, Germany

Locate your Regional Sales Representative at:
www.draeger.com/contact

USA
 Draeger Medical, Inc.
 3135 Quarry Road
 Telford, PA 18969-1042, USA
 Tel +1 215 721 5400
 Toll-free +1 800 437 2437
 Fax +1 215 723 5935
info.usa@draeger.com

CANADA
 Draeger Medical Canada Inc.
 2425 Skymark Avenue, Unit 1
 Mississauga, Ontario, L4W 4Y6
 Tel +1 905 212 6600
 Toll-free +1 866 343 2273
 Fax +1 905 212 6601
Canada.support@draeger.com

Manufacturer:
 Drägerwerk AG & Co. KGaA
 Moislinger Allee 53-55
 23558 Lübeck, Germany


The quality management system at Draeger Medical Systems, Inc. is certified according to ISO 13485, ISO 9001 and Annex II.3 of Directive 93 / 42 / EEC (Medical devices).

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
WARNING

Surgical light


To properly use this medical device,
read and comply with these
instructions for use.

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Intended use of Polaris 600 light

Polaris product family

The Polaris 600 light system is classified as a surgical light system in accordance with IEC 60601-2-41 and is intended to be used for the local illumination of the operating field and examination field on the patient in operating and treatment rooms.

The intended uses remain unchanged when the Polaris 600 light is combined with a Polaris 100 or 200 light.

Ceiling-mounted versions

Single light

As a single light, the Polaris 600 light is intended to be used as a surgical or examination light in operating and treatment rooms for diagnostic and treatment purposes which can be interrupted without posing a hazard to the patient if the light fails.

As a single light, the Polaris 600 light is classified as a "small surgical light" in accordance with IEC60601-2-41.

Double or triple light combination

The combination of 2 or more Polaris 600 lights as a surgical light system enables use in operating and treatment rooms and meets the requirements for a fail-safe surgical light system.

Intended use of the Dräger disposable handle

The Dräger disposable handle provides a sterile cover for the associated inner handle of Polaris systems. The Dräger disposable handle is intended solely for single use.

Intended use of the MedView camera 2.1.8.16

The MedView camera is intended for recording the operating field for the purpose of patient documentation, teaching, and training.

In the Polaris 600 light system, the MedView camera is used together with a medical display as a video system. The MedView camera must only be used in conjunction with the Polaris 600 light.

The recordings must only be used for non-diagnostic purposes.

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Application

Environments of use

Polaris 600 lights are intended for use in rooms used for medical purposes, especially operating rooms and treatment rooms.

Description

Light

The Polaris 600 light is certified and approved as a surgical light. The Polaris 600 light uses LED bulbs. This ensures low energy consumption and a long service life of the bulb. The small size of the LED bulbs makes it possible to use a modern design which facilitates ease of handling and cleaning while also maintaining good flow behavior under laminar flow ceilings. 2.1.5

Color temperature

2.1.8.13

The color temperature can be adjusted in 4 steps, which the surgeon can select depending on tissue structure, operating field, and individual perception.

Cold white light increases concentration, while warm white light reduces dazzle effects and increases the contrast.

OR light mode 2.1.8.1

The Polaris 600 light combines the illumination from 92 white LED bulbs to create a homogeneous light column with deep cavity illumination and minimal shadowing. The adjustable color temperature, good color reproduction, and adjustable illuminance ensure optimum working conditions for demanding procedures within the operating field and the examination field. The illuminance can be adjusted in the OR light mode.

Ambient light mode (Endo light) 2.1.4

In the ambient light mode, the surroundings can be illuminated for minimally invasive procedures. In this light mode, the illuminance is significantly lower than in the OR light mode. The illuminance cannot be adjusted in the ambient light mode.

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

Multimedia system with 2 Polaris 600 lights (configuration example)

A Ceiling cover

B IR receiver for Polaris 600 remote control for MedView camera

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

Polaris 600 light with sterilizable handle E 2.1.1

- A Gimbal mounting
- B Integrated handle
- C Light control panel
- D Sterilizable handle E
- E Lens with LED
- F Bottom glass

Polaris 600 light with sterilizable handle E and MedView camera 2.1.8.16


- A Gimbal mounting
- B Integrated handle
- C Light control panel
- D Sterilizable handle E
- E Lens with LED
- F Bottom glass
- G MedView camera

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
Light control panel 2.1.6

- A key: Light on/standby
 - B key: Switch synchronization on/off 2.1.7
 - C key: Reduce illuminance
 - D key: Increase illuminance
 - E key: Set color temperature
 - F key: Set light field diameter
 - G key: Switch ambient light (Endo light) on/off
- For further information on controlling the lighting, see "Operation" on page 35.

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
2.1.7

Wall-mounted control panel for 2 Polaris 600 lights and a MedView camera (configuration example)

For a Polaris 600 light with MedView camera, the camera control area is on the right-hand side of the shared wall-mounted control panel.

Light control area

- A key: Light on/standby
- B key: Reduce illuminance
- C key: Increase illuminance
- D key: Set color temperature
- E key: Set light field diameter
- F key: Switch ambient light (Endo light) on/off
- G key: Switch synchronization on/off

Camera control area

- H key: MedView camera on/standby
- I key: Zoom out video image (Wide)
- J key: Zoom in video image (Tele)
- K key: Switch automatic exposure on/off
- L key: Reduce exposure (manually)
- M key: Increase exposure (manually)
- N AF key: Switch autofocus on/off
- O key: Focus on distant object (manually)
- P key: Focus on near object (manually)
- Q key: Switch freeze frame on/off

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
Operating concept

	Polaris 600 remote control for MedView camera
Light field diameter	-
Illuminance	-
Color temperature	-
Synchronization	-
Ambient light (Endo light)	-
MedView camera	X

Sterilizable handle E 2.1.8.3

illuminate and the light field diameter.

NOTE

The sterilizable handle E will lose its functionality if disposable handles and the associated handle adapter for disposable handles from manufacturers other than Dräger are used. It will no longer be possible to set the illuminance or the light field diameter in the operating section.

The operating section of the handle is touch-sensitive, i.e., it is sufficient to lay a finger on the surface. It is not necessary to exert any pressure. The better the contact between finger and surface, the easier it is to adjust the light.

The operating section of the handle can be operated when surgical gloves are worn.

- A wiping motion with one finger across the surface of the operating section at any location will adjust the illuminance, see "Setting illuminance on the sterilizable handle E" on page 43.
- Double-tapping with one finger on the surface of the operating section at any location will adjust the light field diameter, see "Setting illuminance on the sterilizable handle E" on page 43.

The handle is used both for positioning and for adjusting the light.

It consists of the operating section (A) and the handle section (B).

The handle section is used for sterile positioning.


The operating section is used for setting the

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
Correct operation of the sterilizable handle E


The touch-sensitive part of the handle is designed so that, for correct operation, only one finger must touch the operating section at any one time. If several fingers touch the operating section simultaneously, the operating signal will not be implemented.

Correct operation

- One finger at a time is touching the operating section

Operating signal will not be implemented

- Several fingers are touching the operating section at the same time

Wall-mounted control panel (option) 2.1.7

The wall-mounted control panel can incorporate up to 3 light control panels. Optionally it can include a camera control panel.

The keys on the wall-mounted control panel are pressure-sensitive, i.e., a certain pressure must be applied by the finger.

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