

cobas[®] 4800 HPV Test



FOR *IN VITRO* DIAGNOSTIC USE.

cobas [®] 4800 System Sample Preparation Kit	c4800 SMPL PREP	960 Tests 240 Tests	P/N: 05235804190 P/N: 05235782190
cobas [®] 4800 HPV Amplification/Detection Kit	c4800 HPV AMP/DET	960 Tests 240 Tests	P/N: 05235910190 P/N: 05235901190
cobas [®] 4800 HPV Controls Kit	c4800 HPV CTLS	10 Sets	P/N: 05235855190
cobas [®] 4800 System Liquid Cytology Preparation Kit	c4800 LIQ CYT	960 Tests 240 Tests	P/N: 05235839190 P/N: 05235812190
cobas [®] 4800 System Wash Buffer Kit	c4800 WB	960 Tests 240 Tests	P/N: 05235871190 P/N: 05235863190

NOTICE: The purchase of this product allows the purchaser to use it for amplification and detection of nucleic acid sequences by polymerase chain reaction (PCR) and related processes for human *in vitro* diagnostics. No general patent or other license of any kind other than this specific right of use from purchase is granted hereby.

INTENDED USE

คุณสมบัติทั่วไป และคุณลักษณะเฉพาะข้อ 2.1

The cobas[®] 4800 Human Papillomavirus (HPV) Test is a qualitative *in vitro* test for the detection of Human Papillomavirus in patient specimens. The test utilizes amplification of target DNA by the Polymerase Chain Reaction (PCR) and nucleic acid hybridization for the detection of 14 high-risk (HR) HPV types in a single analysis. The test specifically identifies HPV16 and HPV18 while concurrently detecting the other high risk types (31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68) at clinically relevant infection levels. Specimens are limited to cervical cells collected in Roche Cell Collection Medium (Roche Molecular Systems, Inc.), cobas[®] PCR Cell Collection Media (Roche Molecular Systems, Inc.), PreservCyt[®] Solution (Hologic Corp.) and SurePath[™] Preservative Fluid (BD Diagnostics[™] InPath[™].)

คุณลักษณะเฉพาะข้อ 2.2

Indications for use of the cobas[®] 4800 HPV Test are:

- The cobas[®] 4800 HPV Test is indicated for use in screening patients with ASC-US (atypical squamous cells of undetermined significance) cervical cytology results to determine the need for referral to colposcopy.
- The cobas[®] 4800 HPV Test is indicated for use in screening patients with ASC-US cervical cytology results to assess the presence or absence of high-risk HPV genotypes 16 and 18.
- The cobas[®] 4800 HPV Test is indicated for use adjunctively with cervical cytology to assess the presence or absence of high risk HPV types.
- The cobas[®] 4800 HPV Test is indicated for use adjunctively with cervical cytology to assess the presence or absence of HPV genotypes 16 and 18.
- The cobas[®] 4800 HPV Test is indicated for use as a first-line primary screening test to identify women at increased risk for the development of cervical cancer or presence of high-grade disease.
- The cobas[®] 4800 HPV Test is indicated for use as a first-line primary screening test to assess the presence or absence of HPV genotypes 16 and 18.

The results from the cobas[®] HPV Test, together with the physician's assessment of cytology history, other risk factors, and professional guidelines, may be used to guide patient management. The results of the cobas[®] HPV Test are not intended to prevent women from proceeding to colposcopy.

SUMMARY AND EXPLANATION OF THE TEST

Persistent infection with human papillomavirus (HPV) is the cause of cervical cancer and its precursor cervical intraepithelial neoplasia (CIN)¹⁻³. The presence of HPV has been implicated in greater than 99% of cervical cancers, worldwide³. HPV is a small, non-enveloped, double-stranded DNA virus, with a genome of approximately 8000 nucleotides. There are more than 118 different types of HPV^{4,5}, and approximately 40 different HPVs that can infect the human anogenital mucosa^{6,7}. However, only a subset of 13 to 18 of these types is considered high-risk for the development of cervical cancer and its precursor lesions^{3,9-13}. In an analysis of data from the International Agency of Research on Cancer (IARC) multi-center case-control study, the pooled OR (Odds Ratio) for squamous-cell cervical cancer with HPV infection was 158.2 when the analysis was restricted to studies using well validated HPV detection techniques¹². In this study, the odds ratios for cervical cancer ranged from 109 to 278 in studies from different parts of the world¹².

Although persistent infection with high-risk (HR) HPV is a necessary cause of cervical cancer and its precursor lesions, a very small percentage of infections progress to these disease states. Sexually transmitted infection with HPV is extremely common, with estimates of up to 75% of all women experiencing exposure to HPV at some point¹⁴. However, > 90% of infected women will mount an effective immune response and clear the infection in 6 to 24 months without any long term health consequences¹⁵⁻²⁰. An infection with any HPV type can produce cervical intraepithelial neoplasia (CIN) although this also usually resolves once the HPV infection has been cleared²¹.

In developed countries with cervical cancer screening programs, the Pap smear has been used since the mid-1950s as the primary tool to detect early precursors to cervical cancer. Although it has decreased the death rates due to cervical cancer dramatically in those countries, the Pap smear requires interpretation by highly trained cytopathologists and is a relatively inaccurate test with a high rate of false negatives. Cytological abnormalities observed in the Pap smear are primarily due to infection with HPV; however, various inflammatory or sampling variations can result in false positive Pap results. Triage of an abnormal Pap smear involves repeat testing, colposcopy and biopsy. A histologically confirmed high-grade lesion must be surgically removed in order to prevent the development of invasive cervical cancer.

Papillomavirus is extremely difficult to culture *in vitro*, and not all patients infected with HPV have a demonstrable antibody response. Nucleic acid (DNA) testing by PCR is a non-invasive method for determining the presence of a cervical HPV infection. The implementation of HPV DNA testing has increased the efficiency of cervical cancer screening programs by detecting high-risk lesions earlier in women 30 years and older with NILM cytology and by reducing the need for unnecessary colposcopy and treatment in patients 21 and older with ASC-US cytology. Furthermore, the superior

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

The Document Revision Information section is located at the end of this document.
05641225001-21EN

๑. ลงชื่อ.....ประธานกรรมการ Doc Rev. 19.0
๒. ลงชื่อ.....กรรมการ
๓. ลงชื่อ.....กรรมการ

sensitivity of HPV testing over Pap smears for the detection of high grade disease in a screening population has been well documented^{22,23}. With superior sensitivity established, HPV DNA testing as a first-line primary screening test has been proposed and adopted in some screening programs.

PRINCIPLES OF THE PROCEDURE ข้อ 3.2 และข้อ 3.3

The **cobas**[®] 4800 HPV Test is based on two major processes: (1) automated specimen preparation to simultaneously extract HPV and cellular DNA; (2) PCR amplification²⁴ of target DNA sequences using both HPV and β -globin specific complementary primer pairs and real-time detection of cleaved fluorescent-labeled HPV and β -globin specific oligonucleotide detection probes. The concurrent extraction, amplification and detection of β -globin in the **cobas**[®] 4800 HPV Test monitors the entire test process.

The Master Mix reagent for the **cobas**[®] 4800 HPV Test contains primer pairs and probes specific for the 14 high-risk HPV types and β -globin DNA. The detection of amplified DNA (amplicon) is performed during thermal cycling using oligonucleotide probes labeled with four different fluorescent dyes. The amplified signal from twelve high-risk HPV types (31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68), is detected using the same fluorescent dye, while HPV16, HPV18 and β -globin signals are each detected with their own dedicated fluorescent dye.

Specimen Preparation ข้อ 3.1 และข้อ 3.8

Specimen preparation for the **cobas**[®] 4800 HPV Test is automated with the use of the **cobas** x 480 instrument. Cervical specimens collected in Roche Cell Collection Medium, **cobas**[®] PCR Cell Collection Media, PreservCyt[®] Solution or SurePath[™] Preservative Fluid are digested under denaturing conditions at elevated temperatures and then lysed in the presence of chaotropic reagent. Released HPV nucleic acids, along with the β -globin DNA serving as process control, are purified through absorption to magnetic glass particles, washed and finally separated from these particles, making them ready for PCR amplification and detection.

PCR Amplification

Target Selection

The **cobas**[®] 4800 HPV Test uses primers to define a sequence of approximately 200 nucleotides within the polymorphic L1 region of the HPV genome. A pool of HPV primers present in the Master Mix is designed to amplify HPV DNA from 14 high-risk types (16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68)^{3,9,13,25}. Fluorescent oligonucleotide probes bind to polymorphic regions within the sequence defined by these primers.

An additional primer pair and probe target the human β -globin gene (330 bp amplicon) to provide a process control. ข้อ 3.5.1 Internal cellular control

Target Amplification

EagleZ05 DNA Polymerase²⁶, a chemically modified version of *Thermus species* Z05 DNA polymerase²⁷, is utilized for "hot start" amplification of the HPV targets and the β -globin control. First, the PCR reaction mixture is heated to activate EagleZ05 DNA Polymerase, to denature the viral DNA and genomic DNA and to expose the primer target sequences. As the mixture cools, the upstream and downstream primers anneal to the target DNA sequences. The EagleZ05 DNA Polymerase, in the presence of divalent metal ion and excess dNTPs, extends the primer(s), and a second DNA strand is synthesized. This completes the first cycle of PCR, yielding a double-stranded DNA copy of the target region of the HPV genome and β -globin gene. The DNA Polymerase extends the annealed primers along the target templates to produce an approximately 200-base pair double-stranded HPV target DNA molecule or a 330 base pair β -globin DNA molecule termed an amplicon. This process is repeated for a number of cycles, each cycle effectively doubling the amount of amplicon DNA. Amplification occurs only in the region of the HPV genome and/or β -globin gene between the appropriate primer pair. The entire genome is not amplified.

Automated Real-time Detection

The **cobas**[®] 4800 HPV Test utilizes real-time^{28,29} PCR technology. Each oligonucleotide probe in the reaction is labeled with a fluorescent dye that serves as a reporter, and with a quencher that quenches fluorescent emissions from the dye in an intact probe. As amplification progresses, probes that are complementary to the amplicon bind to specific single-stranded DNA sequences and are cleaved by the 5' to 3' nuclease activity of the EagleZ05 DNA Polymerase. Once the reporter dye is separated from the quencher by this nuclease activity, it emits fluorescence of a characteristic wavelength when excited by the proper spectrum of light. This characteristic wavelength for each dye allows HPV16 amplicon, HPV18 amplicon, other HR amplicons (31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68) and the β -globin control to be measured independently because the probes specific for these sequences are labeled with different dyes.

Selective Amplification

ข้อ 3.6

Selective amplification of target nucleic acid from the clinical specimen is achieved in the **cobas**[®] 4800 HPV Test by the use of AmpErase (uracil-N-glycosylase) enzyme and deoxyuridine triphosphate (dUTP). AmpErase enzyme recognizes and catalyzes the destruction of DNA strands containing deoxyuridine²⁸, but not DNA containing deoxythymidine. Deoxyuridine is not present in naturally occurring DNA, but is always present in amplicon due to the use of deoxyuridine triphosphate in place of thymidine triphosphate as one of the dNTPs in the Master Mix reagent; therefore, only amplicon contain deoxyuridine. Deoxyuridine renders contaminating amplicon susceptible to destruction by AmpErase enzyme prior to amplification of the target DNA. AmpErase enzyme, which is included in the Master Mix reagent, catalyzes the cleavage of deoxyuridine-containing DNA at the deoxyuridine residues by opening the deoxyribose chain at the C1-position. When heated in the first thermal cycling step, the amplicon DNA chain breaks at the position of the deoxyuridine, thereby rendering the DNA non-amplifiable. AmpErase enzyme is inactive at temperatures above 55°C, i.e., throughout the thermal cycling steps, and therefore does not destroy target amplicon. AmpErase enzyme in the **cobas**[®] 4800 HPV Test has been demonstrated to inactivate at least 10³ copies of deoxyuridine-containing HPV amplicon per PCR.


คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑. ลงชื่อ..... ประธานกรรมการ

๒. ลงชื่อ..... กรรมการ


๓. ลงชื่อ..... กรรมการ

REAGENTS

cobas® 4800 System Sample Preparation Kit (c4800 SMPL PREP) 240 Tests (P/N: 05235782190)			
Kit components	Reagent ingredients	Quantity per kit	Safety symbol and warning*
MGP (cobas® 4800 System Magnetic Glass Particles)	Magnetic glass particles 93% Isopropanol ^b	10 x 4.6 mL	 <p>DANGER H225: Highly flammable liquid and vapour. H319: Causes serious eye irritation. H336: May cause drowsiness or dizziness. P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P233: Keep container tightly closed. P261: Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray. P280: Wear protective gloves/ eye protection/ face protection. P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. P370 + P378: In case of fire: Use dry sand, dry chemical or alcohol-resistant foam to extinguish. 67-63-0 Propan-2-ol</p>
EB (cobas® 4800 System Elution Buffer)	Tris buffer 0.09% Sodium azide	10 x 18 mL	N/A

* Product safety labeling primarily follows EU GHS guidance

^b Hazardous substance

cobas® 4800 System Sample Preparation Kit (c4800 SMPL PREP) 960 Tests (P/N: 05235804190)			
Kit components	Reagent ingredients	Quantity per kit	Safety symbol and warning*
MGP (cobas® 4800 System Magnetic Glass Particles)	Magnetic glass particles 93% Isopropanol ^b	10 x 13.5 mL	 <p>DANGER H225: Highly flammable liquid and vapour. H319: Causes serious eye irritation. H336: May cause drowsiness or dizziness. P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P233: Keep container tightly closed. P261: Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray. P280: Wear protective gloves/ protective clothing / eye protection/ face protection. P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water. P370 + P378: In case of fire: Use dry sand, dry chemical or alcohol-resistant foam to extinguish. 67-63-0 Propan-2-ol</p>
EB (cobas® 4800 System Elution Buffer)	Tris buffer 0.09% Sodium azide	10 x 18 mL	N/A

* Product safety labeling primarily follows EU GHS guidance

^b Hazardous substance

คณะกรรมการพิจารณาผลการประกาศราคาอิเล็กทรอนิกส์

๑. ลงชื่อ..... ประธานกรรมการ



๒. ลงชื่อ..... กรรมการ

๓. ลงชื่อ..... กรรมการ



cobas® 4800 System Wash Buffer Kit (c4800 WB) 240 Tests (P/N: 05235863190)			
Kit components	Reagent Ingredients	Quantity per kit	Safety symbol and warning
WB (cobas® 4800 System Wash Buffer)	Sodium citrate dihydrate 0.05% N-Methylisothiazolone HCl	10 x 55 mL	N/A

cobas® 4800 System Wash Buffer Kit (c4800 WB) 960 Tests (P/N: 05235871190)			
Kit components	Reagent Ingredients	Quantity per kit	Safety symbol and warning
WB (cobas® 4800 System Wash Buffer)	Sodium citrate dihydrate 0.05% N-Methyl isothiazolone HCl	10 x 200 mL	N/A

cobas® 4800 System Liquid Cytology Preparation Kit (c4800 LIQ CYT) 240 Tests (P/N: 05235812190)			
Kit components	Reagent Ingredients	Quantity per kit	Safety symbol and warning*
PK (cobas® 4800 Proteinase K)	Tris buffer < 0.05% EDTA Calcium chloride Calcium acetate Glycerol < 2% Proteinase K ^b	10 x 0.9 mL	 DANGER H317: May cause an allergic skin reaction. H334: May cause allergy or asthma symptoms or breathing difficulties if inhaled. P261: Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray. P280: Wear protective gloves. P284: Wear respiratory protection. P304 + P340: IF INHALED: Remove person to fresh air and keep comfortable for breathing. P333 + P313: If skin irritation or rash occurs: Get medical advice/ attention. P342 + P311: If experiencing respiratory symptoms: Call a POISON CENTER/doctor. 39450-01-5 Proteinase, Tritirachium album serine
SDS (cobas® 4800 System SDS Reagent)	Tris buffer 0.2% SDS 0.09% Sodium azide	10 x 3 mL	N/A
LVS (cobas® 4800 System Lysis Buffer)	Tris buffer 37% (w/w) Guanidine HCl ^b < 6% Polydocanol ^b	10 x 10 mL	 DANGER H302: Harmful if swallowed. H315: Causes skin irritation. H318: Causes serious eye damage. P264: Wash skin thoroughly after handling. P270: Do not eat, drink or smoke when using this product. P280: Wear protective gloves/ eye protection/ face protection. P301 + P312 + P330: IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell. Rinse mouth. P305 + P351 + P338 + P310: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor. P501: Dispose of contents/ container to an approved waste disposal plant. 50-01-1 Guanidinium chloride 9002-92-0 Polydocanol

* Product safety labeling primarily follows EU GHS guidance

^b Hazardous substance

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์



๑. ลงชื่อ..... ประธานกรรมการ

๒. ลงชื่อ..... กรรมการ

๓. ลงชื่อ..... กรรมการ

Doc Rev. 19.0

cobas® 4800 System Liquid Cytology Preparation Kit (c4800 LIQ CYT)
860 Tests (P/N: 05235839190)

Kit components	Reagent Ingredients	Quantity per kit	Safety symbol and warning*
PK (cobas® 4800 Proteinase K)	Tris buffer < 0.05% EDTA Calcium chloride Calcium acetate Glycerol < 2% Proteinase K ^b	20 x 1.2 mL	 DANGER H317: May cause an allergic skin reaction. H334: May cause allergy or asthma symptoms or breathing difficulties if inhaled. P261: Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray. P280: Wear protective gloves. P284: Wear respiratory protection. P304 + P340: IF INHALED: Remove person to fresh air and keep comfortable for breathing. P333 + P313: If skin irritation or rash occurs: Get medical advice/ attention. P342 + P311: If experiencing respiratory symptoms: Call a POISON CENTER/doctor. 39450-01-6 Proteinase, Tritirachium album serine
SDS (cobas® 4800 System SDS Reagent)	Tris buffer 0.2% SDS 0.09% Sodium azide	10 x 9 mL	N/A
LYS (cobas® 4800 System Lysis Buffer)	Tris buffer 37% (w/v) Guanidine HCl ^b < 5% Polydocanol ^b	10 x 36 mL	 DANGER H302: Harmful if swallowed. H315: Causes skin irritation. H318: Causes serious eye damage. P264: Wash skin thoroughly after handling. P270: Do not eat, drink or smoke when using this product. P280: Wear protective gloves/ eye protection/ face protection. P301 + P312 + P330: IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell. Rinse mouth. P305 + P351 + P338 + P310: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor. P501: Dispose of contents/ container to an approved waste disposal plant. 50-01-1 Guanidinium chloride 9002-92-0 Polydocanol

* Product safety labeling primarily follows EU GHS guidance


^b Hazardous substance

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑.ลงชื่อ.....ประธานกรรมการ

๒.ลงชื่อ.....กรรมการ

๓.ลงชื่อ.....กรรมการ



cobas® 4800 HPV Amplification/Detection Kit (c4800 HPV AMP/DET) 240 Tests (P/N: 05235901190)			
Kit components	Reagent ingredients	Quantity per kit	Safety symbol and warning
HPV MMX (cobas® 4800 HPV Master Mix)	Tricine buffer Potassium acetate Potassium hydroxide Glycerol < 0.13% dATP, dCTP, dGTP, dUTP < 0.01% Upstream and downstream HPV primers < 0.01% Upstream and downstream β-globin primers < 0.01% Fluorescent-labeled HPV probes < 0.01% Fluorescent-labeled β-globin probes < 0.10% EagleZ05 DNA polymerase (microbial) < 0.10% AmpErase (uracil-N-glycosylase) enzyme (microbial) 0.09% Sodium azide	10 x 0.5 mL	N/A
HPV Mg/Mn (cobas® 4800 HPV Mg/Mn Solution)	Magnesium acetate Manganese acetate < 0.02% Glacial acetic acid 0.09% Sodium azide	10 x 1.0 mL	N/A

cobas® 4800 HPV Amplification/Detection Kit (c4800 HPV AMP/DET) 960 Tests (P/N: 05235910190)			
Kit components	Reagent ingredients	Quantity per kit	Safety symbol and warning
HPV MMX (cobas® 4800 HPV Master Mix)	Tricine buffer Potassium acetate Potassium hydroxide Glycerol < 0.13% dATP, dCTP, dGTP, dUTP < 0.01% Upstream and downstream HPV primers < 0.01% Upstream and downstream β-globin primers < 0.01% Fluorescent-labeled HPV probes < 0.01% Fluorescent-labeled β-globin probes < 0.10% EagleZ05 DNA polymerase (microbial) < 0.10% AmpErase (uracil-N-glycosylase) enzyme (microbial) 0.09% Sodium azide	20 x 1.0 mL	N/A
HPV Mg/Mn (cobas® 4800 HPV Mg/Mn Solution)	Magnesium acetate Manganese acetate < 0.02% Glacial acetic acid 0.09% Sodium azide	10 x 1.0 mL	N/A

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑. ลงชื่อ..... ประธานกรรมการ

๒. ลงชื่อ..... กรรมการ

๓. ลงชื่อ..... กรรมการ

cobas® 4800 HPV Controls Kit (c4800 HPV CTLs) 10 Sets (P/N: 05235855190)			
Kit components	Reagent ingredients	Quantity per kit	Safety symbol and warning
HPV (+) C (cobas® 4800 HPV Positive Control)	Tris buffer EDTA 0.05% Sodium azide < 0.00001% Poly rA RNA (synthetic) < 0.00001% Non-infectious plasmid DNA (microbial) containing HPV 18, 18, 39 sequences < 0.00001% Non-infectious plasmid DNA (microbial) containing human β-globin sequences	10 x 0.5 mL	N/A
(-) C (cobas® 4800 System Negative Control)	Tris buffer EDTA 0.05% Sodium azide < 0.00001% Poly rA RNA (synthetic)	10 x 0.5 mL	N/A

NOTE: Product safety labeling primarily follows EU GHS guidance.

WARNINGS AND PRECAUTIONS

A. FOR *IN VITRO* DIAGNOSTIC USE.

- B. This test is for use with cervical specimens collected using Roche Cell Collection Medium, cobas® PCR Cell Collection Media, PreservCyt® Solution and SurePath™ Preservative Fluid.
- C. Do not pipette by mouth.
- D. Do not eat, drink or smoke in laboratory work areas. Wear protective disposable gloves, laboratory coats and eye protection when handling specimens and kit reagents. Wash hands thoroughly after handling specimens and test reagents.
- E. Avoid microbial and DNA contamination of reagents.
- F. Dispose of unused reagents and waste in accordance with country, federal, state and local regulations.
- G. Do not use reagents after their expiration dates.
- H. Do not pool reagents.
- I. Material Safety Data Sheets (MSDS) are available on request from your local Roche office.
- J. Gloves must be worn and must be changed between handling specimens and cobas® 4800 reagents to prevent contamination.
- K. Specimens should be handled as infectious using safe laboratory procedures such as those outlined in *Biosafety in Microbiological and Biomedical Laboratories*³¹ and in the CLSI Document M29-A3³².
- L. **LYS** contains guanidine hydrochloride. **Do not allow direct contact between guanidine hydrochloride and sodium hypochlorite (bleach) or other highly reactive reagents such as acids or bases. These mixtures can release a noxious gas.** If liquid containing guanidine hydrochloride is spilled, clean with suitable laboratory detergent and water. If a spill occurs with potentially infectious agents, **FIRST** clean the affected area first with laboratory detergent and water, and then with 0.5% sodium hypochlorite.
- M. **MGP** contains isopropanol and is highly flammable. Keep away from open flames and potential spark producing environments.
- N. **EB, SDS, HPV MMX, HPV Mg/Mn, (-) C, and HPV (+) C** contain sodium azide. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. While disposing of sodium azide containing solutions down laboratory sinks, flush the drains with a large volume of cold water to prevent azide buildup.
- O. Wear eye protection, laboratory coats and disposable gloves when handling any reagents. Avoid contact of these materials with the skin, eyes or mucous membranes. If contact does occur, immediately wash with large amounts of water. Burns can occur if left untreated. If spills occur, dilute with water before wiping dry.
- P. All disposable items are for one time use. Do not reuse.
- Q. Do not use sodium hypochlorite solution (bleach) for cleaning the cobas x 480 instrument or cobas z 480 analyzer. Clean the cobas x 480 instrument or cobas z 480 analyzer according to procedures described in the cobas® 4800 System - User Assistance.
- R. For additional warnings, precautions and procedures to reduce the risk of contamination for the cobas x 480 instrument or cobas z 480 analyzer, consult the cobas® 4800 System - User Assistance.

STORAGE AND HANDLING REQUIREMENTS

- A. **Do not freeze reagents.**
- B. Store MGP, EB, PK, SDS, LYS, HPV MMX, HPV Mg/Mn, HPV (+) C and (-) C at 2-8°C. These reagents are stable until the expiration date indicated.
- C. Store WB at 15-25°C. This reagent is stable until the expiration date indicated.

05641225001-21EN

๑.ลงชื่อ.....ประธานกรรมการ
 ๒.ลงชื่อ.....กรรมการ
 ๓.ลงชื่อ.....กรรมการ

Doc Rev. 19.0

MATERIALS PROVIDED

A. cobas® 4800 System Sample Preparation Kit (P/N: 05235782190)	c4800 SMPL PREP	240 Tests
MGP (cobas® 4800 System Magnetic Glass Particles) EB (cobas® 4800 System Elution Buffer)		
B. cobas® 4800 System Sample Preparation Kit (P/N: 05235804190)	c4800 SMPL PREP	960 Tests
MGP (cobas® 4800 System Magnetic Glass Particles) EB (cobas® 4800 System Elution Buffer)		
C. cobas® 4800 System Wash Buffer Kit (P/N: 05235863190)	c4800 WB	240 Tests
WB (cobas® 4800 System Wash Buffer)		
D. cobas® 4800 System Wash Buffer Kit (P/N: 05235871190)	c4800 WB	960 Tests
WB (cobas® 4800 System Wash Buffer)		
E. cobas® 4800 HPV Amplification/Detection Kit (P/N: 05235901190)	c4800 HPV AMP/DET	240 Tests
HPV MMX (cobas® 4800 HPV Master Mix) HPV Mg/Mn (cobas® 4800 HPV Mg/Mn Solution)		
F. cobas® 4800 HPV Amplification/Detection Kit (P/N: 05235910190)	c4800 HPV AMP/DET	960 Tests
HPV MMX (cobas® 4800 HPV Master Mix) HPV Mg/Mn (cobas® 4800 HPV Mg/Mn Solution)		
G. cobas® 4800 System Liquid Cytology Preparation Kit (P/N: 05235912190)	c4800 LIQ CYT	240 Tests
PK (cobas® 4800 Proteinase K) SDS (cobas® 4800 System SDS Reagent) LYS (cobas® 4800 System Lysis Buffer)		
H. cobas® 4800 System Liquid Cytology Preparation Kit (P/N: 05235839190)	c4800 LIQ CYT	960 Tests
PK (cobas® 4800 Proteinase K) SDS (cobas® 4800 System SDS Reagent) LYS (cobas® 4800 System Lysis Buffer)		
I. cobas® 4800 HPV Controls Kit (P/N: 05235855190)	c4800 HPV CTL5	10 Sets
HPV (+) C (cobas® 4800 HPV Positive Control) (-) C (cobas® 4800 System Negative Control)		

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑. ลงชื่อ..... ประธานกรรมการ

๒. ลงชื่อ..... กรรมการ

๓. ลงชื่อ..... กรรมการ

MATERIALS REQUIRED BUT NOT PROVIDED

Specimen and Reagent Handling

- **cobas**® PCR Cell Collection Media (Roche P/N 05619637190, optional)
- Roche Cell Collection Medium (Roche P/N 07994745190, optional)
- Roche Cell Collection Medium Replacement Caps (Roche P/N 08037230190, optional)
- Cervical Collection Brush (Roche P/N 08399832190, optional)
- CO-RE Tips, 1000 µL, rack of 96 (Roche P/N 04639642001 or Hamilton P/N 235905)
- 50 mL Reagent Reservoir (Roche P/N 05232732001)
- 200 mL Reagent Reservoir (Roche P/N 05232759001)
- For HPV ASAP v2.0.1 use **cobas**® 4800 System Extraction (deep well) plate 1.6 mL (Roche P/N 05232716001)
- For HPV ASAP v2.1 use **cobas**® 4800 System Extraction (deep well) plate 2.0 mL (Roche P/N 06884008001)
- **cobas**® 4800 System AD (microwell) plate 0.3 mL and Sealing Film (Roche P/N 05232724001)
- Solid waste bag [Roche P/N 05530873001 (small) or 04691989001 (large)]
- Hamilton STAR Plastic Chute (Roche P/N 04639669001)
- Tubes 13 mL Round Base (Roche: P/N 07958048190) for use as secondary sample tubes
- Caps, neutral color (Roche P/N 07958056190; for recapping post-run specimens in 13 mL Round Base tubes)
- Disposable gloves, powderless

Instrumentation and Software

- **cobas** x 480 instrument
- **cobas** z 480 analyzer
- **cobas**® 4800 System Control Unit with System Software version 2.2 or higher
- **cobas**® 4800 System **cobas**® HPV AP software version 2.0. or higher

Optional Equipment and Materials

- **cobas**® Sample Prep Buffer (Roche P/N.06526985190; Tris buffered detergent)*
- Pipettes: capable of delivering 1000 µL
- Aerosol barrier DNase-free tips: capable of delivering 1000 µL
- Centrifuge equipped with a swinging bucket rotor with minimum RCF of 1500
- Stand-alone magnetic plate (Roche P/N 05440777001)
- Vortex Mixer (single tube)
- Multi-tube vortexer [e.g. VWR P/N 58816-116]
- Heat-resistant barcode labels (RACO Industries; Cat # RAC-225075-9501)
- Thermometer -20/150°C (VWR Cat# 89095-600) or equivalent
- Digital Heater Block 120V (VWR Cat# 75838-294) or equivalent
- 12-Hole Heat Block Module 16mm (VWR Cat# 13259-162) or equivalent

* An open bottle of **cobas**® Sample Prep Buffer (CSPB) may be stored at ambient temperature (15-30°C) for up to 21 days and up to 4 separate uses for the pre-analytic treatment of SurePath™ samples.

SPECIMEN COLLECTION, TRANSPORT AND STORAGE

NOTE: Handle all specimens as if they are capable of transmitting infectious agents.

ข้อ 3.7

A. Specimen Collection

Cervical specimens collected in Roche Cell Collection Medium, **cobas**® PCR Cell Collection Media, PreservCyt® Solution and SurePath™ Preservative Fluid have been validated for use with the **cobas**® 4800 HPV Test. Follow the manufacturer's instructions for collecting cervical specimens.

B. Specimen Transport

Cervical specimens collected in Roche Cell Collection Medium, **cobas**® PCR Cell Collection Media, PreservCyt® Solution and SurePath™ Preservative Fluid can be transported at 2-30°C. Transportation of HPV specimens must comply with country, federal, state and local regulations for the transport of etiologic agents³³.

C. Specimen Storage

Cervical specimens collected in Roche Cell Collection Medium, **cobas**® PCR Cell Collection Media and PreservCyt® Solution may be stored at 2-30°C for up to 6 months after the date of collection. Cervical specimens collected in SurePath™ Preservative Fluid may be stored at 2-8°C for up to 6 months or may be stored at 15-30°C for up to 6 weeks after the date of collection provided that SurePath™ Preservative Fluid matrix-induced crosslinks are reversed through treatment with **cobas**® Sample Prep Buffer prior to HPV testing.

INSTRUCTIONS FOR USE

NOTE: All reagents except HPV MMX and HPV Mg/Mn must be at ambient temperature prior to loading on the **cobas** x 480 instrument. The HPV MMX and HPV Mg/Mn may be taken directly from 2-8°C storage as they will equilibrate to ambient temperature on board the **cobas** x 480 instrument by the time they are used in the process.

๑.ลงชื่อ.....ประธานกรรมการ

๒.ลงชื่อ.....กรรมการ

๓.ลงชื่อ.....กรรมการ

NOTE: Specimens in Roche Cell Collection Medium, cobas® PCR Cell Collection Media, PreservCyt® Solution and SurePath™ Preservative Fluid must be at ambient temperature before loading on the cobas x 480 instrument.

NOTE: Refer to the cobas® 4800 System - User Assistance for detailed operating instructions.

Run Size: ข้อ 3.5 ตารางวิเคราะห์ได้สูงสุด 96 เทสต์ต่อรอบการตรวจรวม Control
The cobas® 4800 System is designed to support the cobas® 4800 HPV Test with run sizes from 1 to 94 specimens plus controls (up to 96 assays per run). Each cobas® 4800 System Sample Preparation Kit, cobas® 4800 System Liquid Cytology Preparation Kit and cobas® 4800 System Wash Buffer Kit contains reagents sufficient for 10 runs of either 24 tests (240 tests per kit) or 96 tests (960 tests per kit). Each cobas® 4800 HPV Amplification/Detection Kit contains reagents sufficient for 10 runs of either 24 tests (240 tests per kit) or 96 tests (960 tests per kit); multiple 240 Test Kits can be used to optimize reagent usage for 48 or 72 tests. The cobas® 4800 HPV Controls Kit contains reagents sufficient for a total of 10 runs (10 sets per kit). The minimum run size on the cobas® 4800 System is 1 specimen plus controls. One replicate of the cobas® 4800 System Negative Control [(-) C] and one replicate of the cobas® 4800 HPV Positive Control [HPV (+) C] are required to perform each test run (see "Quality Control" section).

Workflow:

NOTE: Although not an optimal use of reagents, a System Sample Preparation 960 Test Kit can be used for a 24 sample run and an HPV Amplification/Detection 960 Test Kit can be used for a 24, 48, or 72 sample run.

The cobas® 4800 HPV Test can be run using either of two workflows, referred to as "Full Workflow" or "Recovery Workflow" within the cobas® 4800 Software.

HPV Full Workflow:

The "HPV Full Workflow" consists of sample preparation on the cobas x 480 instrument followed by amplification/detection on the cobas z 480 analyzer. Run size can be a 24-test format (from 1 to 22 specimens plus 2 controls) or a 96-test format (from 1 to 94 specimens plus 2 controls). Refer to the "Performing a Full Workflow" section below and the cobas® 4800 System - User Assistance for details.

HPV Recovery Workflow:

The "HPV Recovery Workflow" consists of manual PCR plate setup using eluate from the processed deep well plate followed by amplification/detection on the cobas z 480 analyzer. Refer to the "Performing a Recovery Workflow" section below and the cobas® 4800 System - User Assistance for details.

Specimens:

There are four cervical specimen types that can be assayed using the cobas® 4800 HPV Test: a) cervical specimens in Roche Cell Collection Medium, b) cervical specimens in PreservCyt® Solution, c) cervical specimens in cobas® PCR Cell Collection Media and d) cervical specimens in SurePath™ Preservative Fluid. Roche Cell Collection Medium, PreservCyt® Solution, and cobas® PCR Cell Collection Media specimens may be processed directly out of their primary containers with a proper barcode or out of a properly barcoded 13 mL round-based tube on the cobas x 480 instrument. SurePath™ specimens must be transferred into a properly barcoded 13 mL round-based tube for specimen treatment (See Treatment of SurePath™ primary specimens section) and processing on the cobas x 480 instrument. Consult the cobas® 4800 System - User Assistance for proper barcoding procedures and the list of acceptable barcodes for the cobas® 4800 System.

NOTE: SurePath™ specimens must be treated with cobas® Sample Prep Buffer to reverse matrix-induced cross-links prior to HPV testing on the cobas® 4800 System.

Treatment of SurePath™ primary specimens with cobas® Sample Prep Buffer to reverse matrix-induced crosslinks

NOTE: Heat-resistant barcodes are required for tubes used to reverse matrix-induced cross-links (see the Optional Equipment and Materials section).

NOTE: It is recommended that steps B, C, G and H below are done in a biological hood to minimize possible cross-contamination.

- A. Prepare a barcoded 13 mL round-based tube with 0.5 mL of cobas® Sample Prep Buffer for each SurePath™ specimen to be tested. An open bottle of cobas® Sample Prep Buffer (CSPB) may be stored at ambient temperature (15-30°C) for up to 21 days and up to 4 separate uses for the pre-analytic treatment of SurePath™ samples.
- B. Vortex SurePath™ specimens for 10 seconds prior to transfer. Transfer 0.5 mL of each SurePath™ specimen into a 13 mL round-based tube prepared in step A. Re-cap each tube before moving to the next. Always change pipet tips for each specimen.
- C. Vortex each tube for 1 second.
- D. Transfer tubes to the heating unit set at 120°C (see Optional Equipment and Materials section). Up to 48 tubes can be processed per batch.
- E. Heat for 20 minutes.
- F. After heating, remove tubes to a collection rack and cool at ambient temperature for 10 minutes.
- G. Vortex each tube for 5 seconds.
- H. Transfer tubes to 24 position cobas® 4800 specimen racks, discard caps and process on the cobas® 4800 System for HPV testing.

SurePath™ specimens treated with cobas® Sample Prep Buffer can be stored for future HPV testing if, for example, cytology evaluation is required first. The following procedure should be followed:

- A. Follow the treatment procedure above to step G.
- B. Store tubes with SurePath™ specimens treated with cobas® Sample Prep Buffer at 2-30°C for up to 4 weeks prior to HPV testing on the cobas® 4800 System.

NOTE: The minimum volume required in the Roche Cell Collection Medium, cobas® PCR Cell Collection Media, PreservCyt® Solution primary containers is 3.0 mL. When using 13 mL round-based secondary tubes, fill to a minimum volume of 1.0 mL and a maximum volume of 10 mL.

NOTE: Use only Roche Cell Collection Medium, cobas® PCR Cell Collection Media, PreservCyt® Solution and SurePath™ Preservative Fluid to collect cervical specimens for the cobas® 4800 HPV Test. The cobas® 4800 HPV Test has not been validated with other media types. Using the cobas® 4800 HPV Test with other media types could lead to false negative, false positive and/or invalid results.

๑.ลงชื่อ.....ประธานกรรมการ
๒.ลงชื่อ.....กรรมการ
๓.ลงชื่อ.....กรรมการ

NOTE: To avoid cross-contamination of processed specimens, additional caps for vials (see Materials Required But Not Provided section) should be used to recap specimens after processing. Re-cap tightly. Store and ship the vials in upright orientation.

NOTE: It may be necessary to aliquot specimens into barcoded 13 mL round-based tubes for processing on the cobas x 480 instrument. Use pipettors with aerosol-barrier or positive-displacement tips to handle specimens. To avoid cross-contamination of processed specimens, additional caps for these tubes in an alternate color (neutral; see Materials Required But Not Provided section) should be used to recap these specimens after processing.

NOTE: Use caution when transferring specimens from primary containers to 13 mL round-based secondary tubes. Vortex primary specimens prior to transfer. Change pipet tips after each specimen.

NOTE: Do not process specimens which appear bloody or have a dark brown color.

A single run can have any combination of specimens (Roche Cell Collection Medium, cobas[®] PCR Cell Collection Media, PreservCyt[®] Solution and/or SurePath[™] Preservative Fluid) and each specimen can be tested with either the HPV High Risk or HPV High Risk Plus Genotyping sub-tests.

Workflows

Performing a Full Workflow:

- A. The cobas[®] 4800 HPV Test may be used for 1 to 94 specimens plus one cobas[®] 4800 System negative control and one cobas[®] 4800 HPV positive control.
- B. Perform the system startup and maintenance procedures by following the instructions in the cobas[®] 4800 System - User Assistance.
- C. Start a new run by clicking the "New run" button.
- D. In the Selection test window, select Workflow type "Full" then select the Test "HPV".
- E. Enter a run name or leave as the default run name, then click "OK" to proceed.
- F. Follow the software wizard guide to load specimens.

NOTE: Specimens can be loaded in barcoded primary or secondary tubes in any order.

NOTE: If primary containers for Roche Cell Collection Medium, cobas[®] PCR Cell Collection Media, and PreservCyt[®] Solution specimens are used for processing, vortex prior to loading.

- G. Select a Specimen type for each specimen.
 - Choose "PC" for ordering Roche Cell Collection Medium, PreservCyt[®] Solution, or cobas[®] PCR Cell Collection Media specimens.
 - Choose "SP" for ordering SurePath[™] Preservative Fluid specimens.
- H. Select the Request result for each specimen.
 - Choose Requested result "HPV High Risk Panel" to report any one of, or combination of high risk HPV types: 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68 test results.
 - Choose Requested result "HPV High Risk Panel Plus Genotyping" to report any one of, or combination of high risk HPV types: 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68 and to separately report high risk HPV Type 16 and high risk HPV Type 18 test results.
- I. Follow the software wizard guide to load all consumables.
- J. Follow the software wizard guide to load all reagents.

NOTE: Controls [HPV (+) C and (-) C] are not loaded together with specimens. They are loaded onto the reagent carrier during reagent loading. Two positions (A1 and B1) on each of the Extraction plate and microwell plate are reserved for the HPV (+) and (-) controls, respectively.

NOTE: The cobas[®] 4800 System has an internal clock to monitor the length of time the reagents are on-board. Once the WB is scanned, 1 hour is allowed to complete the loading process and click on the Start button. A countdown timer is displayed on the Workplace Tab. The system will not allow the run to start if the on-board time has expired.

NOTE: To assure the accurate transfer of MGP, vortex or vigorously shake the MGP vial prior to dispensing into the reagent reservoir.

- K. Load the sample preparation reagents (WB, MGP, EB, SDS and LYS) into the barcoded reagent reservoirs using the "scan-scan-pour-place" method:
- Scan the reagent bottle barcode ข้อ 3.9
 - Scan the reagent reservoir barcode
 - Pour the reagent into the reservoir
 - Place the filled reagent reservoir into the designated position on the reagent carrier

L. The reagent reservoirs are available in two sizes: 200 mL and 50 mL. Follow the software wizard guide to select the appropriate reagent reservoir sizes. The reagent reservoir barcodes must face to the right of the carrier.

NOTE: Amplification/detection reagents (HPV MMX and HPV Mg/Mn), Controls [HPV (+) C and (-) C] and PK are loaded directly onto the reagent carrier and scanned by the cobas x 480 instrument automatically.

NOTE: All reagents and reagent reservoirs are barcoded and designed for one time use. The cobas[®] 4800 Software tracks the use of the reagents and reagent reservoirs and rejects previously used reagents or reagent reservoirs. The software also verifies that sufficient reagents are loaded on the instrument.

NOTE: The cobas[®] 4800 Software tracks the expiration date of all reagents. Reagents that are beyond their expiration date will not be accepted for use on the cobas[®] 4800 System.

M. Start sample preparation by clicking on "Start Run".

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

05641225001-21EN

๑.ลงชื่อ.....ประธานกรรมการ

๒.ลงชื่อ.....กรรมการ

๓.ลงชื่อ.....กรรมการ

Doc Rev. 19.0

- N. After successful completion of sample preparation, click "Unload" to unload the plate carrier.
- ** The status of sample preparation can be reviewed at this point, prior to clicking "Unload". See the **cobas**[®] 4800 System - User Assistance.
- O. Follow the instructions in the **cobas**[®] 4800 System - User Assistance to seal the microwell plate, transport the plate to the **cobas z 480** analyzer and start the amplification and detection run.

NOTE: The **cobas[®] 4800 System has an internal clock to monitor the length of time after addition of the prepared samples to working master mix. Amplification and detection should be started as soon as possible but no later than 90 minutes after the end of the **cobas x 480** instrument run. A countdown timer is displayed on the Workplace Tab.**

- P. When the amplification and detection run is completed, unload the microwell plate from the **cobas z 480** analyzer.
- Q. Follow the instructions in the **cobas**[®] 4800 System - User Assistance to review and accept results.

Performing a Recovery Workflow

NOTE: The Recovery Workflow is available as a recovery option in the event that the full workflow cannot be completed due to circumstances beyond the user's control (e.g. power failure during amplification/detection run).

NOTE: Only samples successfully processed on the **cobas x 480 instrument can be amplified/detected using the Recovery run. System surveillance for reagents and consumables is limited during the Recovery run. No sample position tracking is provided when using the Recovery workflow – the end user must ensure that the actual position of a sample on the microwell plate corresponds to the one designated in the Recovery Plate Layout Report Work Order file. Extreme care must be exercised while preparing the microwell plate to ensure proper PCR set-up and to avoid contamination.**

NOTE: Samples processed on the **cobas x 480 instrument have limited stability. They must be amplified/detected using the Recovery workflow within 24 hours if stored at 2°C to 30°C.**

- A. Start a Recovery run by clicking the New run button.
- B. In the Test Selection window, select "Recovery" then select test type "HPV".
- C. Enter a run name or leave as the default run name, then click OK to proceed.
- D. Select a run to recover.
- E. If using HPV ASAP v2.1, scan the original DWP ID from the full workflow.
- F. Enter the new MWP ID.
- G. Enter the Master Mix and Metal Ions IDs for all Amplification/Detection reagent vials in the kit.
- H. Prepare the **cobas**[®] 4800 HPV working master mix:
1. For a 240 Test Kit, add 240 µL of HPV Mg/Mn to one vial of HPV MMX (0.5 mL vial from 240 Test Kit).
 2. For a 960 Test Kit, add 450 µL of HPV Mg/Mn to each of the two vials of HPV MMX (1.0 mL vials from 960 Test Kit).

NOTE: The Recovery run must be started within 90 minutes of addition of HPV Mg/Mn to the HPV MMX. The system does not monitor the length of time after addition of the prepared samples to working master mix in the Recovery workflow. The end user must ensure that amplification and detection is started within the allotted time.

- I. Thoroughly mix working master mix by carefully inverting the vial(s). Do not vortex the working master mix.
- J. Transfer 25 µL of working Master Mix to the required wells in the microwell plate.
- K. Place the Extraction plate from the run to be repeated onto the stand-alone magnetic plate.
- L. Manually transfer 25 µL of eluate from the Extraction plate wells to the corresponding wells in the microwell plate. Ensure that well positions are maintained (e.g. eluate in A1 well in Extraction plate is transferred to A1 on the microwell plate). Ensure that no MGP is carried over to the microwell plate.
- M. Follow the instructions in the **cobas**[®] 4800 System - User Assistance to seal the microwell plate.
- N. Centrifuge the microwell plate using a swinging bucket rotor for at least 5 seconds at 1500 RCF.
- O. Transport the plate to the **cobas z 480** analyzer and start the amplification and detection run.
- P. When the amplification and detection run is completed, unload the microwell plate from the **cobas z 480** analyzer.
- Q. Follow the instructions in the **cobas**[®] 4800 System - User Assistance to review and accept results.

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑. ลงชื่อ..... ประธานกรรมการ

๒. ลงชื่อ..... กรรมการ

๓. ลงชื่อ..... กรรมการ

Interpretation of Results

NOTE: All assay and run validation is performed by the cobas® 4800 Software.

NOTE: A valid run may include both valid and invalid specimen results.

For a valid run, specimen results are interpreted as shown in Tables 1 and 2:

Table 1
Result Interpretation of the cobas® 4800 HPV Test for Presence of HPV DNA

cobas® 4800 HPV Test	Result Report and Interpretation
Requested Result "HPV High Risk Panel":	
POS HR HPV	High Risk HPV Positive Specimen is positive for the DNA of any one of, or combination of, the following high risk HPV types: 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68.
NEG HR HPV	High Risk HPV Negative* HPV types 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68 DNA were undetectable or below the pre-set threshold.
Invalid HR HPV	High Risk HPV Invalid The results for HR HPV are invalid. For PreservCyt® specimens, the original specimen should be re-tested no more than two times to obtain valid results. If the results are still invalid a new specimen should be obtained. For SurePath™ specimens the original specimen should be retested if sufficient volume remains. If the results are still invalid a new specimen should be obtained.
Failed	No Result for Specimen Consult the cobas® 4800 System - User Assistance for instructions to review run flags and recommended actions. Original specimen should be re-tested to obtain valid result.
Requested Result "HPV High Risk Panel Plus Genotyping":	
POS Other HR HPV	Other High Risk HPV Positive Specimen is positive for the DNA of any one of, or combination of, the following high risk HPV types: 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68.
NEG Other HR HPV	Other High Risk HPV Negative* HPV types 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68 DNA were undetectable or below the pre-set threshold.
Invalid Other HR HPV	Invalid Other High Risk HPV The result for Other HR HPV is invalid. For PreservCyt® specimens, the original specimen should be re-tested no more than two times to obtain valid results. If the results are still invalid a new specimen should be obtained. For SurePath™ specimens the original specimen should be retested if sufficient volume remains. If the results are still invalid a new specimen should be obtained.
POS HPV16	HPV16 Positive Specimen is positive for HPV type 16 DNA.
NEG HPV16	HPV16 Negative* HPV type 16 DNA was undetectable or below the pre-set threshold.
Invalid HPV16	Invalid HPV16 The result for HPV16 is invalid. For PreservCyt® specimens, the original specimen should be re-tested no more than two times to obtain valid results. If the results are still invalid a new specimen should be obtained. For SurePath™ specimens the original specimen should be retested if sufficient volume remains. If the results are still invalid a new specimen should be obtained.
POS HPV18	HPV18 Positive Specimen is positive for HPV type 18 DNA.
NEG HPV18	HPV18 Negative* HPV type 18 DNA was undetectable or below the pre-set threshold.
Invalid HPV18	Invalid HPV18 The result for HPV18 is invalid. For PreservCyt® specimens, the original specimen should be re-tested no more than two times to obtain valid results. If the results are still invalid a new specimen should be obtained. For SurePath™ specimens the original specimen should be retested if sufficient volume remains. If the results are still invalid a new specimen should be obtained.
Failed	No Result for Specimen Consult the cobas® 4800 System - User Assistance for instructions to review run flags and recommended actions. Original specimen should be re-tested to obtain valid results.

*A negative result does not preclude the presence of HPV infection because results depend on adequate specimen collection, absence of inhibitors, and sufficient DNA to be detected.

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑. ลงชื่อ..... ประธานกรรมการ

๒. ลงชื่อ..... กรรมการ

๓. ลงชื่อ..... กรรมการ

Table 2
Result Interpretation of the cobas® 4800 HPV Test for Patients with Cytological Abnormalities

Results	Interpretation
NEG Other HR HPV*, NEG HPV16, NEG HPV18	Very low likelihood of underlying \geq CIN2.
POS Other HR HPV*, NEG HPV16, NEG HPV18	Increased likelihood that underlying \geq CIN2 will be detected at colposcopy.
POS HPV16 and/or POS HPV18	Highest likelihood that underlying \geq CIN2 will be detected at colposcopy ^{34,35} .

*Other HR HPV DNA includes the following types: 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68.

NOTE: HPV negative results are not intended to prevent women from proceeding to colposcopy.

NOTE: In addition to the results tabulated above, invalid results for one or more combinations is also possible. If such a result is obtained, for example:

Other HR HPV NEG, HPV16 POS, HPV18 Invalid

The positive and negative results should be interpreted as shown in Table 1. In this example, HPV 18 results are invalid. The original specimen should be re-tested no more than two times to obtain valid results. If the results are still invalid a new specimen should be obtained.

NOTE: Negative results indicate HPV DNA concentrations are undetectable or below the pre-set threshold.

NOTE: Positive test results indicates the presence of any one or more of the high risk types, but since patients are often co-infected with low-risk types it does not rule out the presence of low-risk types in patients with mixed infections.

NOTE: Results of this test should only be interpreted in conjunction with information available from clinical evaluation of the patient and patient history.

LIST OF RESULT FLAGS

The following table lists common flags for the cobas® 4800 HPV Test which are relevant for result interpretation. Refer to the cobas® 4800 System - User Assistance for a full list of flags.

Table 3
List of flags for cobas® 4800 HPV Test

Flag code	Description	Recommended action
R20	Positive control is invalid.	Positive control values were invalid. 1. Repeat entire run with fresh reagents. 2. If the problem persists, contact Roche Service.
R21	Negative control is invalid.	Negative control values were invalid. To avoid carryover, use Good Laboratory Practice. 1. Repeat entire run with fresh reagents. 2. If the problem persists, contact Roche Service.
X3	Error: Clot was detected. Sample was not processed.	Make sure that the samples were handled according to the workflow description. 1. Check the sample for clots. 2. Rerun the sample.
X4	Error: Pipetting error occurred. Sample was not processed.	Insufficient sample volume or mechanical error during pipetting is the most likely reason. 1. Make sure that there is enough sample volume. 2. Check whether the tip eject plate is placed correctly. 3. Rerun the sample.

QUALITY CONTROL

One set of cobas® 4800 HPV Test Positive and Negative Controls are included in each run. For any run, valid results must be obtained for both the Positive and Negative Control for the cobas® 4800 Software to display the reportable cobas® 4800 HPV Test results from that run.

Positive Control

The HPV (+) Control result must be 'Valid'. If the HPV (+) Control results are consistently invalid, contact your local Roche office for technical assistance.

Negative Control

The (-) Control result must be 'Valid'. If the (-) Control results are consistently invalid, contact your local Roche office for technical assistance.

คณะกรรมการพิจารณาผลการตรวจราคาอิเล็กทรอนิกส์
 ๑. ลงชื่อ..... ประธานกรรมการ
 ๒. ลงชื่อ..... กรรมการ
 ๓. ลงชื่อ..... กรรมการ

PROCEDURAL PRECAUTIONS

As with any test procedure, good laboratory technique is essential to the proper performance of this assay. Due to the high analytical sensitivity of this test, care should be taken to keep reagents and amplification mixtures free of contamination.

PROCEDURAL LIMITATIONS

1. The **cobas**[®] 4800 HPV Test detects DNA of the high-risk types 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68. This test does not detect DNA of HPV low-risk types (e.g. 6, 11, 42, 43, 44) since there is no clinical utility for testing of low-risk HPV types³⁶.
2. The **cobas**[®] 4800 HPV Test for detection of human papillomavirus types 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68 is not recommended for evaluation of suspected sexual abuse.
3. The performance of the **cobas**[®] 4800 HPV Test has not been adequately established for HPV vaccinated individuals³⁷.
4. Prevalence of HPV infection in a population may affect performance. Positive predictive values decrease when testing populations with low prevalence or individuals with no risk of infection.
5. Infection with HPV is not an indicator of cytologic HSIL or underlying high-grade CIN, nor does it imply that CIN2-3 or cancer will develop. Most women infected with one or more high-risk HPV types do not develop CIN2-3 or cancer.
6. A negative high-risk HPV result does not exclude the possibility of future cytologic HSIL or underlying CIN2-3 or cancer.
7. Test only the indicated specimen type. The **cobas**[®] 4800 HPV Test has only been validated for use with cervical specimens collected in Roche Cell Collection Medium, **cobas**[®] PCR Cell Collection Media, PreservCyt[®] Solution and SurePath[™] Preservative Fluid.
8. Detection of high-risk HPV is dependent on the number of copies present in the specimen and may be affected by specimen collection methods, patient factors, stage of infection and the presence of interfering substances.
9. Beta-globin amplification and detection is included in the **cobas**[®] 4800 HPV Test to differentiate HPV negative specimens from those that do not exhibit HPV signal due to insufficient cell mass in the specimen. All HPV negative specimens must have a valid Beta-globin signal within a pre-defined range to be identified as valid negatives by the **cobas**[®] 4800 System.
10. Reliable results are dependent on adequate specimen collection, transport, storage and processing. Follow the procedures in this Package Insert and the **cobas**[®] 4800 System - User Assistance.
11. The addition of AmpErase enzyme into the **cobas**[®] 4800 HPV Master Mix enables selective amplification of target DNA; however, good laboratory practices and careful adherence to the procedures specified in this Package Insert are necessary to avoid contamination of reagents.
12. Use of this product must be limited to personnel trained in the techniques of PCR and the use of the **cobas**[®] 4800 System.
13. Only the **cobas** x 480 instrument and **cobas** z 480 analyzer have been validated for use with this product. No other sample preparation instrument or PCR system can be used with this product.
14. Due to inherent differences between technologies, it is recommended that, prior to switching from one technology to the next, users perform method correlation studies in their laboratory to qualify technology differences.
15. The effects of other potential variables such as vaginal discharge, use of tampons, douching, etc. and specimen collection variables have not been evaluated.
16. Though rare, mutations within the highly conserved regions of the genomic DNA of Human papillomavirus covered by the **cobas**[®] 4800 HPV Test's primers and/or probes may result in failure to detect the presence of the viral DNA.
17. The presence of PCR inhibitors may cause false negative or invalid results.
18. Cervical specimens often show visibly detectable levels of whole blood as a pink or light brown coloration. These specimens are processed normally on the **cobas**[®] 4800 System. If concentrations of whole blood exceeds 2% (dark red or brown coloration) in Roche Cell Collection Medium, **cobas**[®] PCR Cell Collection Media, or PreservCyt[®] Solution, or above 4% in SurePath[™] Preservative Fluid treated with **cobas**[®] Sample Prep Buffer, there is a likelihood of obtaining a false-negative result. See Interference results for details.
19. Use of the vaginal moisturizer Replens[®] has been associated with false-negative results in SurePath[™] Preservative Fluid.
20. Use of the RepHresh[®] vaginal hygiene products has been associated with false-negative results in Roche Cell Collection Medium and PreservCyt[®] Solution.
21. Removal of red blood cells from Roche Cell Collection Medium, PreservCyt[®], or SurePath[™] specimens through treatment with glacial acetic acid (GAA) has not been validated with the **cobas**[®] 4800 HPV Test. Any use of GAA treatments with the **cobas**[®] 4800 HPV Test must be validated by the testing laboratory.

PERFORMANCE CHARACTERISTICS

Performance Comparison to a CE Mark Comparator HPV Test

Clinical sensitivity and specificity to disease status (\geq CIN2) was determined for the **cobas**[®] 4800 HPV Test and a CE Mark comparator HPV test³⁸ in a population of women at least 21 years old with ASC-US cytology results, determined through routine cervical cancer screening. All testing was carried out using PreservCyt[®] Solution cervical specimens. A total of 1578 subjects with an Initial ASC-US cytology result underwent colposcopy and had valid HPV tests and cervical biopsy results. Disease status of the subjects was ascertained by a central pathology review panel from the biopsy specimens obtained at colposcopy. The results for an ASC-US population are summarized in Table 4 and indicate that the **cobas**[®] 4800 HPV Test performance was comparable to the comparator test.

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑.ลงชื่อ.....ประธานกรรมการ

๒.ลงชื่อ.....กรรมการ

๓.ลงชื่อ.....กรรมการ

Table 4
Comparison of the Performance of the cobas[®] 4800 HPV Test and a CE Mark Comparator HPV test in Detecting \geq CIN2 and \geq CIN3 in the ASC-US Population

	cobas [®] 4800 HPV Test		CE Marked HPV Test	
	Point Estimate	95% CI	Point Estimate	95% CI
\geq CIN2				
Sensitivity (%)	90.0 (72/80)	(81.5, 94.8)	87.2 (68/78) ¹	(78.0, 92.9)
Specificity (%)	70.5 (1,056/1,498)	(68.1, 72.7)	71.1 (1,056/1,485) ²	(68.8, 73.4)
PPV (%)	14.0 (72/514)	(12.8, 15.3)	13.7 (68/497)	(12.4, 15.1)
NPV (%)	99.2 (1,056/1,064)	(98.6, 99.6)	99.1 (1,056/1,066)	(98.3, 99.5)
Prevalence (%)	5.1 (80/1578)	(4.1, 6.3)	5.0 (78/1563)	(4.0, 6.2)
\geq CIN3				
Sensitivity (%)	93.5 (43/46)	(82.5, 97.8)	91.3 (42/46)	(79.7, 96.6)
Specificity (%)	69.3 (1,053/1,517)	(66.9, 71.5)	70.0 (1,062/1,517)	(67.7, 72.3)
PPV (%)	8.4 (43/514)	(7.6, 9.2)	8.5 (42/497)	(7.6, 9.4)
NPV (%)	99.7 (1,061/1,064)	(99.2, 99.9)	99.6 (1,062/1,066)	(99.0, 99.9)
Prevalence (%)	2.9 (43/1578)	(2.2, 3.9)	3.0 (46/1563)	(2.2, 3.9)

¹ Results for two subjects with a \geq CIN2 diagnosis could not be determined by the CE Mark Comparator HPV Test due to insufficient volume resulting from repeated testing.

² Results for thirteen subjects with a $<$ CIN2 diagnosis could not be determined by the CE Mark Comparator HPV Test due to insufficient volume resulting from repeated testing.

In women \geq 30 years with normal cytology, the risk of cervical disease (\geq CIN2) is 7.29 fold higher with a High Risk positive cobas[®] 4800 HPV Test result than with a negative cobas[®] 4800 HPV Test result. Relative risk estimates and their 95% confidence intervals are presented in Table 5.

In women 30 years or older, the cobas[®] 4800 HPV Test can be used to assess the presence or absence of HPV genotypes 16 and 18. The risk of cervical disease (\geq CIN2) is 13.71 fold higher with an HPV16 and/or HPV18 positive cobas[®] 4800 HPV Test result than with a negative result and the risk is 2.51 fold higher with an HPV16 and/or HPV18 positive cobas[®] 4800 HPV Test result compared to a positive cobas[®] 4800 HPV Test result for the 12 other high risk types. In all cases, the lower bound of the 95% confidence interval exceeds 1, suggesting a statistically higher risk of developing cervical disease with a positive HPV test result.

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์
 ๑.ลงชื่อ.....ประธานกรรมการ
 ๒.ลงชื่อ.....กรรมการ
 ๓.ลงชื่อ.....กรรมการ

Table 5
Relative Risk for Cervical Disease (\geq CIN2 by Central Pathology Review) in Women \geq 30 Years with Normal Cytology*

HPV Result	Relative Risk Estimate	95% CI*
Pos vs. Neg	7.29	(3.99, 22.11)
16+/18+ vs. Neg	13.71	(7.31, 41.92)
16+/18+ vs. 12 other HR+	2.51	(1.73, 3.61)
Note: 0.5 was added to a zero cell of the estimated number of diseased subjects in any of the 1000 bootstrap samples		
* 95% CI is 2.5 and 97.5 percentile of bootstrap CI based on 1000 bootstrap samples		

NILM (\geq 30 Years) Population - Performance Evaluation

For the NILM (\geq 30 years) population, estimates of sensitivity and specificity along with 95% CIs for HR HPV positive vs. HR HPV negative are presented in Table 6 for unadjusted results.

The unadjusted sensitivity and the specificity of the test for \geq CIN2 histology were 83.2% (109/131) with 95% CI: 75.9% to 88.6% and 60.4% (2492/4127) with 95% CI: 58.9% to 61.9%, respectively. The unadjusted sensitivity and specificity of the cobas[®] HPV Test for detecting \geq CIN3 histology were 90.0% (72/80) with 95% CI: 81.5% to 94.8% and 60.0% (2506/4178) with 95% CI: 58.5% to 61.5%, respectively.

Table 6
Performance of cobas[®] 4800 HPV Test In the NILM (\geq 30 years) Population (Unadjusted Estimates)

CPR Diagnosis	Performance	Estimate	95% CI
\geq CIN2	Sensitivity (%)	83.2 (109/131)	(75.9, 88.6)
	Specificity (%)	60.4 (2492/4127)	(58.9, 61.9)
	PPV(%)	6.3 (109/1744)	(5.8, 6.8)
	NPV(%)	99.1 (2492/2514)	(98.7, 99.4)
	Prevalence (%)	3.1 (131/4258)	(2.8, 3.6)
\geq CIN3	Sensitivity (%)	90.0 (72/80)	(81.5, 94.8)
	Specificity (%)	60.0 (2506/4178)	(58.5, 61.5)
	PPV(%)	4.1 (72/1744)	(3.8, 4.5)
	NPV(%)	99.7 (2506/2514)	(99.4, 99.8)
	Prevalence (%)	1.9 (80/4258)	(1.5, 2.3)

Overall (\geq 25 Years) Population - Comparison of Performance of HPV testing vs. Cytology

Clinical performance of the cobas[®] HPV Test and liquid based cytology (PreservCyt[®]) was determined in a population of 40,901 women 25 years and older, independent of cytology status (Overall Population). For the Overall (\geq 25 Years) population, estimates of sensitivity and specificity for the cobas[®] HPV Test vs. cytology for the detection of \geq CIN2 and \geq CIN3 are presented³⁸ in Table 7. The unadjusted sensitivities of the cobas[®] HPV Test and cytology for detection of \geq CIN2 were 88.2% (380/431) with 95% CI 84.8-90.9% and 51.5% (222/431) with CI 46.8-56.2%, respectively. The unadjusted sensitivities of the cobas[®] HPV Test and cytology for detection of \geq CIN3 were 92.0% (252/274) with 95% CI 88.1-94.6 and 53.3% (146/274) with CI 47.4-59.1%, respectively. Verification bias adjusted specificities of the cobas[®] HPV Test and cytology for detection of \geq CIN2 were 90.5% (36343/40163) with CI 90.2-90.8% and 94.1% (37811/40163) with CI 93.9-94.4%, respectively.

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์
 ๑.ลงชื่อ.....ประธานกรรมการ
 ๒.ลงชื่อ.....กรรมการ
 ๓.ลงชื่อ.....กรรมการ

Table 7
Comparison of performance of the cobas 4800
HPV Test and Cytology for the detection of \geq CIN2 and \geq CIN3 in the Overall (\geq 25 years) Population

	Cytology		cobas [®] HPV Test	
	% (n)	95% CI	% (n)	95% CI
\geq CIN2				
Sensitivity	51.5 (222/431)	(46.8-56.2)	88.2 (380/431)	(84.8-90.9)
Specificity	73.4 (5428/7392)	(72.4-74.4)	57.8 (4270/7392)	(56.6-58.9)
PPV	10.2 (222/2186)	(9.3-11.1)	10.9 (380/3502)	(10.4-11.3)
NPV	96.3 (5428/5637)	(95.9-96.6)	98.8 (4270/4321)	(98.5-99.1)
\geq CIN3				
Sensitivity	53.3 (146/274)	(47.4-59.1)	92.0 (252/274)	(88.1-94.6)
Specificity	73.0 (5509/7549)	(72.0-74.0)	58.9 (4299/7549)	(55.8-58.1)
PPV	6.7 (146/2186)	(6.0-7.4)	7.2 (252/3502)	(6.9-7.5)
NPV	97.7 (5509/5637)	(97.4-98.0)	99.5 (4299/4321)	(99.2-99.7)

Limit of Detection: PreservCyt[®] Solution and SurePath[™] Preservative Fluid

The limit of detection (LOD) of high risk HPV genotypes HPV16, HPV18 and HPV31 was determined for the cobas[®] 4800 HPV Test. The LODs were assessed using 1) plasmids of HPV31, HPV16 and HPV18 in the background of pooled HPV negative patient specimens collected in PreservCyt[®] Solution and SurePath[™] Preservative Fluid, and 2) HPV positive cell lines SiHa (HPV16) and HeLa (HPV18) in PreservCyt[®] Solution and SurePath[™] Preservative Fluid containing an HPV negative cell line (HCT-15) background. Plasmid and cell lines were diluted to concentrations below, above and at the expected LOD levels. A minimum of 60 replicates were tested for each plasmid or cell line level in both PreservCyt[®] Solution and SurePath[™] Preservative Fluid for each of 3 reagent lots. All testing in SurePath[™] specimen background was done using treatment with cobas[®] Sample Prep Buffer. The LOD is the level of HPV DNA in the sample that has positive test results at least 95% of the time. Tables 8 and 9 contain results from the reagent lot producing the most conservative (highest) LOD in the analysis for PreservCyt[®] Solution and SurePath[™] Preservative Fluid, respectively.

Table 8

Limit of Detection Levels for HPV Types 31, 16, 18 and Cell Lines SiHa (HPV16) and HeLa (HPV18) in PreservCyt[®] Solution

HPV Type	Titer (copies or cells/mL)	Number of Positive/Tested	% Positives	95% Confidence Interval	
				Lower	Upper
31	600	60/60	100%	94%	100%
31	300	59/61	97%	89%	100%
31	150	49/60	82%	70%	90%
16	1500	60/60	100%	94%	100%
16	600	60/60	100%	94%	100%
16	300	55/61	90%	80%	96%
18	1,500	60/60	100%	94%	100%
18	600	60/60	100%	94%	100%
18	300	42/61	69%	56%	80%
SiHa (HPV 16)	200	66/66	100%	95%	100%
SiHa (HPV 16)	100	64/65	98%	92%	100%
SiHa (HPV 16)	50	57/60	95%	86%	99%
HeLa (HPV 18)	80	60/60	100%	94%	100%
HeLa (HPV 18)	40	60/60	100%	94%	100%
HeLa (HPV 18)	20	56/60	93%	84%	98%

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑. ลงชื่อ..... ประธานกรรมการ

๒. ลงชื่อ..... กรรมการ

๓. ลงชื่อ..... กรรมการ

Table 9
Limit of Detection Levels for HPV Types 31, 16, 18 and Cell Lines SiHa (HPV16) and HeLa (HPV18) in SurePath™ Preservative Fluid

HPV Type	Titer (copies or cells/mL)	Number of Positive/Tested	% Positives	95% Confidence Interval	
				Lower	Upper
31	600	60/60	100%	94%	100%
31	300	59/59	100%	94%	100%
31	150	54/60	90%	80%	96%
16	600	60/60	100%	94%	100%
16	300	59/60	98%	91%	100%
16	150	40/60	67%	53%	78%
18	1,500	60/60	100%	94%	100%
18	600	60/60	100%	94%	100%
18	300	55/59	93%	84%	98%
SiHa (HPV 16)	400	60/60	100%	94%	100%
SiHa (HPV 16)	200	60/60	100%	94%	100%
SiHa (HPV 16)	100	55/60	92%	82%	97%
HeLa (HPV 18)	80	60/60	100%	94%	100%
HeLa (HPV 18)	40	59/60	98%	91%	100%
HeLa (HPV 18)	20	43/60	72%	59%	83%

Limit of Detection: Roche Cell Collection Medium

Dilution panels of HPV31 plasmid, HPV16 and HPV18 cell lines in the background of pooled HPV negative patient specimens collected in Roche Cell Collection Medium and PreservCyt® Solution were tested side-by-side. The limit of detection for the cobas® 4800 HPV Test was comparable.

Inclusivity Verification

To verify that the cobas® 4800 HPV Test is capable of accurately detecting all HPV high risk genotypes, the limit of detection (LOD) was determined (Tables 10 and 11) for genotypes 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68. The sensitivity of the cobas® 4800 HPV Test for HPV genotypes 16, 18 and 31 was determined in the Limit of Detection Study described above in this Package Insert. Quantified plasmid stocks of each HPV genotype were diluted into either PreservCyt® Solution or SurePath™ Preservative Fluid containing HPV-negative HCT-15 cells to concentrations below, above and at the expected LOD levels. One lot of reagents was used to produce a minimum of 48 replicates for each positive level in each media. For testing in SurePath™ Preservative Fluid using treatment with cobas® Sample Prep Buffer (Table 11), background material was prepared from cervical specimens collected in SurePath™ Preservative Fluid and tested as 24 replicates each with two lots of reagents. For each HPV type, the reported LOD was defined as the lowest testing concentration having a ≥ 95% positive hit rate with all higher concentrations having at least a 95% hit rate.

Table 10

Summary of High Risk Genotype Limit Of Detection for cobas® 4800 HPV Genotype Inclusivity Study (PreservCyt® Solution)

HPV DNA Type	LOD (copies/mL)	Number of Positive/Tested	Hit Rate	95% Confidence Interval	
				Lower	Upper
33	190	46/48	96%	86%	99%
35	480	48/48	100%	93%	100%
39	80	48/48	100%	93%	100%
45	190	48/48	96%	86%	99%
51	100	48/48	96%	86%	99%
52	2400	48/48	100%	93%	100%
56	1400	48/48	100%	93%	100%
58	480	47/48	98%	89%	100%
59	190	46/48	96%	86%	99%
66	640	48/48	100%	93%	100%
68	450	48/48	100%	93%	100%

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑. ลงชื่อ.....¹⁹ ประธานกรรมการ

๒. ลงชื่อ..... กรรมการ

๓. ลงชื่อ..... กรรมการ

Table 11
Summary of High Risk Genotype Limit Of Detection for cobas® 4800 HPV Genotype Inclusivity Study
(SurePath™ Preservative Fluid)

HPV DNA Type	LOD (copies/mL)	Number of Positive/Tested	Hit Rate	95% Confidence Interval	
				Lower	Upper
33	300	48/48	100%	93%	100%
35	600	47/48	100%	89%	100%
39	150	48/48	100%	93%	100%
45	300	48/48	100%	93%	100%
51	600	46/48	96%	86%	99%
52	4800	48/48	100%	93%	100%
56	1200	46/48	96%	86%	99%
58	600	48/48	100%	93%	100%
59	600	48/48	100%	93%	100%
66	1200	48/48	100%	93%	100%
68	300	48/48	100%	93%	100%

Precision: PreservCyt® Solution and SurePath™ Preservative Fluid

In-house Precision was examined using panel members prepared for the Limit of Detection Study described in this Package Insert. Levels at and above the limit of detection were used for the precision analysis. Panels were prepared by spiking plasmids of HPV31, HPV16 and HPV18 into the background of pooled HPV negative patient specimens collected in PreservCyt® Solution and SurePath™ Preservative Fluid. All testing in SurePath™ specimen background was done using treatment with cobas® Sample Prep Buffer.

The positive hit rates for panel members (PreservCyt® Solution and SurePath™ Preservative Fluid) at and above the LOD are shown in Tables 12 and 13, respectively. Hit rates were above 95% for all plasmid panel levels. The variance in Ct value for the test was analyzed, and contribution from reagent lot, systems, run-to-run, and within-run random factors were calculated and summarized in Table 14 for PreservCyt® Solution and Table 15 for SurePath™ Preservative Fluid. Table 16 shows the Ct value SD and %CV of components of variation in PreservCyt® Solution. Table 17 shows the Ct value SD and %CV of components of variation in SurePath™ Preservative Fluid.

Table 12
Summary of Hit Rates for cobas® 4800 HPV Precision Study At or Above LOD (in PreservCyt® Solution)

Target	Panel Level	Concentration (copies or cells/mL)	N Tests	N Pos	Hit Rate	95% CI for Hit Rate	
						Lower	Upper
HPV31	> LOD	600	186	186	100%	98%	100%
	= LOD	300	187	184	98%	95%	100%
HPV16	> LOD	1,500	186	186	100%	98%	100%
	= LOD	600	186	186	100%	98%	100%
HPV18	> LOD	1,500	186	186	100%	98%	100%
	= LOD	600	186	186	100%	98%	100%

Table 13
Summary of Hit Rates for cobas® 4800 HPV Precision Study At or Above LOD (in SurePath™ Preservative Fluid)

Target	Panel Level	Concentration (copies or cells/mL)	N Tests	N Pos	Hit Rate	95% CI for Hit Rate	
						Lower	Upper
HPV31	> LOD	300	180	180	100%	98%	100%
	= LOD	150	180	175	97%	94%	99%
HPV16	> LOD	600	180	180	100%	98%	100%
	= LOD	300	180	180	100%	98%	100%
HPV18	> LOD	1,500	180	180	100%	98%	100%
	= LOD	600	180	180	100%	98%	100%

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑.ลงชื่อ.....ประธานกรรมการ

๒.ลงชื่อ.....กรรมการ

๓.ลงชื่อ.....กรรมการ

Table 14
Analysis of Ct Value Variance Components for cobas® 4800 HPV Precision Study Panel Levels Prepared in PreservCyt® Solution

Target	Panel Level	N	Mean Elbow	Variance Components/Percent Contribution				
				Rgt Lot	System	Run	Random	Total
HPV16	> LOD	186	36.3	0.038 17%	0 0%	0.111 49%	0.079 35%	0.228 100%
	= LOD	186	37.5	0.025 11%	0 0%	0.042 18%	0.161 71%	0.228 100%
HPV18	> LOD	186	36.6	0.043 16%	0 0%	0.149 58%	0.067 26%	0.259 100%
	= LOD	186	37.8	0.027 10%	0 0%	0.050 19%	0.184 71%	0.261 100%
HPV31	> LOD	186	36.5	0.003 1%	0.002 1%	0.105 35%	0.187 63%	0.297 100%
	= LOD	187	37.6	0.020 3%	0 0%	0.157 24%	0.489 73%	0.666 100%

Table 15
Analysis of Ct Value Variance Components for cobas® 4800 HPV Precision Study Panel Levels Prepared in SurePath™ Preservative Fluid

Target	Panel Level	N	Mean Elbow	Variance Components/Percent Contribution				
				Rgt Lot	System	Run	Random	Total
HPV16	> LOD	180	37.2	0.014 7%	0 0%	0.039 18%	0.157 75%	0.209 100%
	= LOD	180	38.2	0 0%	0 0%	0.090 22%	0.316 78%	0.405 100%
HPV18	> LOD	180	36.3	0.011 5%	0 0%	0.119 58%	0.073 36%	0.204 100%
	= LOD	180	37.7	0 0%	0 0%	0.148 40%	0.219 60%	0.366 100%
HPV31	> LOD	180	37.2	0 0%	0 0%	0.099 20%	0.393 80%	0.493 100%
	= LOD	180	38.1	0.026 3%	0.015 2%	0.038 5%	0.684 90%	0.764 100%

Table 16
Analysis of Ct Value SD and %CV for cobas® 4800 HPV Precision Study Panel Levels Prepared in PreservCyt® Solution

Target	Panel Level	N	Mean Elbow	SD Components/%CV				
				Rgt Lot	System	Run	Random	Total
HPV16	> LOD	186	36.3	0.19 0.50%	0 0.00%	0.33 0.90%	0.28 0.80%	0.48 1.30%
	= LOD	186	37.5	0.16 0.40%	0 0.00%	0.20 0.50%	0.40 1.10%	0.48 1.30%
HPV18	> LOD	186	36.6	0.21 0.60%	0 0.00%	0.39 1.10%	0.26 0.70%	0.51 1.40%
	= LOD	186	37.8	0.16 0.40%	0 0.00%	0.22 0.60%	0.43 1.10%	0.51 1.30%
HPV31	> LOD	186	36.5	0.05 0.10%	0.05 0.10%	0.32 0.90%	0.43 1.20%	0.54 1.50%
	= LOD	187	37.6	0.14 0.40%	0 0.00%	0.40 1.10%	0.70 1.90%	0.82 2.20%

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑. ลงชื่อ..... ประธานกรรมการ

๒. ลงชื่อ..... กรรมการ

๓. ลงชื่อ..... กรรมการ

Table 17

Analysis of Ct Value SD and %CV for cobas® 4800 HPV Precision Study Panel Levels Prepared in SurePath™ Preservative Fluid

Target	Panel Level	N	Mean Elbow	SD Components/%CV				
				Rgt Lot	System	Run	Random	Total
HPV16	> LOD	180	37.2	0.12	0	0.20	0.40	0.46
				0.30%	0.00%	0.50%	1.10%	1.20%
	= LOD	180	38.2	0	0	0.30	0.56	0.64
				0.00%	0.00%	0.80%	1.50%	1.70%
HPV18	> LOD	180	36.3	0.11	0	0.34	0.27	0.45
				0.30%	0.00%	1.00%	0.70%	1.20%
	= LOD	180	37.7	0	0	0.38	0.47	0.61
				0.00%	0.00%	1.00%	1.20%	1.60%
HPV31	> LOD	180	37.2	0	0.02	0.32	0.63	0.70
				0.00%	0.10%	0.80%	1.70%	1.90%
	= LOD	180	38.1	0.16	0.12	0.20	0.83	0.87
				0.40%	0.30%	0.50%	2.20%	2.30%

Precision: Roche Cell Collection Medium

Panels were prepared by spiking HPV16 cell line DNA and HPV18 cell line DNA into a background of pooled HPV negative patient specimens collected in Roche Cell Collection Medium at and above the LOD. Testing of the panels prepared in Roche Cell Collection Medium demonstrated precision comparable to the precision with panels prepared in PreservCyt® Solution.

Analytical Specificity

A panel of bacteria, fungi and viruses, including those commonly found in the female urogenital tract, as well as several Human papillomavirus types classified as low or undetermined risk were tested with the cobas® 4800 HPV Test to assess analytical specificity. The organisms listed in Table 18 were spiked at high concentrations ($\geq 1 \times 10^3$ units/reaction) into HPV negative PreservCyt® Solution specimen background and into HPV negative PreservCyt® Solution specimen spiked with HPV 31, HPV16 and HPV18 plasmid DNA at 3 times the LOD. Organisms with an asterisk were also tested in SurePath™ Preservative Fluid specimen background under the same conditions. Organisms with a double asterisk were tested only in SurePath™ specimen background. All testing in SurePath™ specimen background was done using treatment with cobas® Sample Prep Buffer. Results indicated that none of these organisms interfered with detection of HPV31, HPV16 and HPV18 plasmid DNA or produced a false positive result in the HPV negative specimen.

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑.ลงชื่อ.....ประธานกรรมการ

๒.ลงชื่อ.....กรรมการ

๓.ลงชื่อ.....กรรมการ

Table 18
Microorganisms Tested for Analytical Specificity

<i>Achromobacter xerosis</i>	<i>Haemophilus ducreyi</i>	<i>Streptococcus agalactiae</i> *
<i>Acinetobacter calcaceticus</i>	Hepatitis B virus (HBV)	<i>Streptococcus anginosus</i>
<i>Acinetobacter Iwoffii</i>	Herpes simplex virus 1*	<i>Streptococcus faecalis</i> **
<i>Acinetobacter sp. Genospecies 3</i>	Herpes simplex virus 2*	<i>Streptococcus pyogenes</i> *
<i>Actinomyces israelii</i>	Human immunodeficiency virus (HIV-1)	<i>Streptococcus sanguis</i>
Adenovirus*	<i>Kingella kingae</i>	SV40
<i>Aerococcus viridans</i>	<i>Klebsiella pneumoniae ss ozaenae</i> *	<i>Treponema pallidum</i>
<i>Alcaligenes faecalis</i>	<i>Lactobacillus acidophilus</i> *	<i>Trichomonas vaginalis</i> *
<i>Bacillus thuringiensis</i>	<i>Lactobacillus crispus</i>	<i>Ureaplasma urealyticum</i>
<i>Bacteroides caccae</i> **	<i>Lactobacillus delbrueckii s. lactis</i>	<i>Veillonella parvula</i>
<i>Bacteroides fragilis</i>	<i>Lactobacillus jensenii</i>	<i>Vibrio parahaemolyticus</i>
<i>Bacteroides ureolyticus</i>	<i>Lactobacillus vaginalis</i>	<i>Weissella paramesenteroides</i>
<i>Bifidobacterium longum</i>	<i>Lactococcus lactis cremoris</i>	<i>Yersinia enterocolitica</i>
<i>Bifidobacterium adolescentis</i> *	<i>Legionella pneumophila</i>	HPV 6*
<i>Bifidobacterium brevi</i>	<i>Micrococcus luteus</i>	HPV 11*
<i>Campylobacter jejuni</i>	<i>Mobiluncus curtisii s. curtisii</i>	HPV 26*
<i>Candida albicans</i> *	<i>Moraxella osloensis</i>	HPV 30**
<i>Chlamydia trachomatis</i> *	<i>Morganella morganii</i>	HPV 34**
<i>Chromobacter violaceum</i>	<i>Mycobacterium avium</i>	HPV 40
<i>Citrobacter braakii</i>	<i>Mycobacterium smegmatis</i>	HPV 42
<i>Clostridium adolescentis</i>	<i>Mycoplasma genitalium</i>	HPV 53**
<i>Clostridium beijerinckii</i> **	<i>Mycoplasma hominis</i>	HPV 54
<i>Clostridium perfringens</i>	<i>Neisseria gonorrhoeae</i> *	HPV 55B
<i>Corynebacterium genitalium</i> **	<i>Neisseria meningitidis Serogroup A</i>	HPV 61
<i>Corynebacterium glutamicum</i>	<i>Pasteurella multocida</i>	HPV 62
Cytomegalovirus*	<i>Pediococcus acidilactici</i>	HPV 64
<i>Eikenella corrodens</i>	<i>Peptostreptococcus anaerobius</i> *	HPV 67*
<i>Enterobacter aerogenes</i> **	<i>Propionibacterium acnes</i>	HPV 69*
<i>Enterobacter cloacae</i>	<i>Proteus mirabilis</i> *	HPV 70*
<i>Enterococcus faecalis</i>	<i>Proteus vulgaris</i>	HPV 71
<i>Enterococcus faecium</i> *	<i>Providencia stuartii</i>	HPV 72
Epstein Barr Virus*	<i>Pseudomonas aeruginosa</i>	HPV 73*
<i>Erysipelothrix rhusiopathiae</i>	<i>Pseudomonas fluorescens</i> **	HPV 81
<i>Escherichia coli</i> *	<i>Ruminococcus productus</i>	HPV 82*
<i>Ewingella americana</i>	<i>Salmonella minnesota</i>	HPV 83
<i>Fusobacterium nucleatum</i>	<i>Serratia marcescens</i>	HPV 84
<i>Fusobacterium varium</i> **	<i>Staphylococcus aureus</i> *	HPV 85**
<i>Gemella morbillorum</i>	<i>Staphylococcus epidermidis</i> *	HPV 89 (CP6108)
<i>Gardnerella vaginalis</i>	<i>Staphylococcus saprophyticus</i>	

*Tested in both PreservCyt® and SurePath™ specimen background

**Tested only in SurePath™ specimen background

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑. ลงชื่อ..... ประธานกรรมการ

๒. ลงชื่อ..... กรรมการ

๓. ลงชื่อ..... กรรมการ

Interfering Substances

HPV positive and HPV negative cervical specimens as well as contrived specimens were used to assess the effects of endogenous and exogenous interfering substances that could potentially be present in cervical specimens. Testing materials used in these studies are described in Table 19. The concentrations of endogenous and exogenous substances tested represent conditions that could occur during specimen collection.

Whole blood, Peripheral Blood Mononuclear Cells (PBMC) and cervical mucus were tested as potential endogenous interfering substances found in cervical specimens. All testing in SurePath™ specimen background was done using treatment with cobas® Sample Prep Buffer. Levels of each potential interfering substance tested and performance observations are described in Table 20. No interference was seen for PBMC or cervical mucus at all levels tested. Whole blood showed no interference when present in visually detectable amounts of up to 2% in Roche Cell Collection Medium and PreservCyt® Solution. Whole blood showed no interference when present in visually detectable amounts of up to 4% in SurePath™ Preservative Fluid.

Table 19
Interference Testing Sample Descriptions

Sample type	Description
HPV Positive Cervical Specimens	10 individual HPV positive PreservCyt® Solution specimens were aliquoted for testing with and without endogenous interfering substances.
HPV Negative Cervical Specimens	10 individual HPV negative PreservCyt® Solution specimens were aliquoted for testing with and without endogenous interfering substances.
Contrived HPV Positive Cervical Specimen	HPV positive (channel 1) PreservCyt® Solution specimens were diluted with HPV negative specimen to an approximate level of 3 x LOD. HPV types 16 (channel 2) and 18 (channel 3) plasmids were then added at ~ 3 x LOD. HPV31 plasmid, HPV16 cell line DNA, and HPV18 cell line DNA were added at ~ 3X LOD to HPV negative specimens collected in Roche Cell Collection Medium.
3x LOD PreservCyt® Solution and SurePath™ Preservative Fluid pools	HPV types 31, 16, 18 plasmids were each diluted to 3 x LOD in pools of PreservCyt® Solution and SurePath™ Preservative Fluid negative specimen.

Table 20
Interference Testing Results with Endogenous Interferents

Interferent Tested	Concentrations Tested	Interference Observed	
		PreservCyt®	SurePath™
Whole Blood	1%, 1.5%, 2%, 3%, 4%, 6%, 8% v/v	Above 2%	Above 4%
PBMC	10 ⁴ , 10 ⁵ , 10 ⁶ cells/mL	None	None
Cervical Mucus	Mucus obtained from standard cervical cleaning procedure	None	None

A total of 18 over-the-counter (OTC) feminine hygiene and contraceptive products were tested as potential interfering substances. Types of potential interferents tested and performance observations in Roche Cell Collection Medium, PreservCyt® Solution, and SurePath™ Preservative Fluid 3 x LOD pools are described in Table 21.

Table 21
Interference Testing Results with Exogenous Interferents

Interferent Description	Interference Observed
Contraceptive Gels, Foams	None
Vaginal Lubricants	*Yes
Vaginal Douche	None
Anti-fungal creams containing 1% clotrimazole, Phenazopyridine Hydrochloride, 1% Hydrocortisone, 2% Miconazole nitrate, 6.5% Tioconazole Ointment, 20% Benzocaine	None

* Replens® (topical anti-dryness gel) produced negative results in replicates of the SurePath™ Preservative Fluid 3 x LOD pool. This interference was also seen when SurePath™ material was treated with cobas® Sample Prep Buffer. RepHresh® vaginal hygiene products produced negative results in replicates of the 3 x LOD pools of Roche Cell Collection Medium and PreservCyt® Solution.

๑. ลงชื่อ..... ประธานกรรมการ
 ๒. ลงชื่อ..... กรรมการ
 ๓. ลงชื่อ..... กรรมการ

SurePath™ Specimen Stability for 6 Weeks at 2-30°C using Treatment with cobas® Sample Prep Buffer

Three SurePath™ HPV negative specimen pools were spiked with HPV Type 51 positive SurePath™ specimen material to produce positive, high positive and low positive specimen pools. The low positive pool was at ~ Limit of Detection (LOD) for the test at Day 0 prior to treatment with cobas® Sample Prep Buffer. These pools were stored at 32°C and tested at intervals up to 6 weeks. Pooled materials were treated to reverse matrix-induced crosslinks followed by analysis with the cobas® 4800 HPV Test at each timepoint. All three pools maintained Ct averages below the clinical cutoff for HPV type 51 (40.0 for channel 1) through the 6 week storage period (see Table 22).

Table 22
SurePath™ Specimen Stability Results for 6 Weeks with Treatment using cobas® Sample Prep Buffer

SurePath™ Pools	Average Ct Values*				
	Day 0	Week 1	Week 3	Week 4	Week 6
High Positive	28.7	30.1	30.3	30.6	31.1
Positive	32.9	33.5	34.1	33.9	34.6
Low Positive (~LOD)	36.9	37.9	38.0	38.8	38.7

*Low Positive timepoints tested as 40 replicates; Positive timepoints as 30 replicates; High Positive timepoints as 20 replicates

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์
 ๑.ลงชื่อ.....ประธานกรรมการ
 ๒.ลงชื่อ.....กรรมการ
 ๓.ลงชื่อ.....กรรมการ

REFERENCES

1. Burd, Eileen M. 2003. Human Papillomavirus and Cervical Cancer. *Clinical Microbiology Reviews*. **16**:1-17.
2. zur Hausen, H. 2002. Papillomaviruses and Cancer: From Basic Studies to Clinical Application. *Nat Rev Cancer*. **2(5)**:342-50.
3. Walboomers, Jan M.M., Jacobs, Marcel V., Manos, M.M., et al. 1999. Human Papillomavirus is a Necessary Cause of Invasive Cervical Cancer Worldwide. *Journal of Pathology*. **189**:12-19.
4. Bernard HU. Review: The clinical importance of the nomenclature, evolution and taxonomy of human papillomaviruses. *J Clin Virol*. 2005; **32S**, S1-6.
5. Molijn A, Kleter B, Quint W, van Doorn, L. Review: Molecular diagnosis of human papillomavirus (HPV) infections. *J Clin Virol*. 2005; **32S**:S43-51.
6. zur Hausen H. Roots and perspectives of contemporary papillomavirus research. *J Cancer Res Clin Oncol*. 1996; **122**:3-13.
7. de Villiers EM, Fauquet C, Broker TR, Bernard HU, zur Hausen H. Classification of papillomaviruses. *Virology*. 2004; **324**:17-27.
8. Franco EL, Rohan TE, Villa LL. Epidemiologic evidence and human papillomavirus infection as a necessary cause of cervical cancer. *J Natl Cancer Inst*. 1999; **91**:506-511.
9. Lorincz AT, Reid R, Jenson AB, et al. Human papillomavirus infection of the cervix: relative risk associations of 15 common anogenital types. *Obstet Gynecol*. 1992; **79**:328-37.
10. Bosch, F.X., Manos, M.M., Munoz, N., et al. 1995. International Biological Study on Cervical Cancer (IBSCC) Study Group. Prevalence of Human Papillomavirus in Cervical Cancer: a Worldwide Perspective. *Journal of the National Cancer Institute*, Vol. 87, No. 11:796-802.
11. Bosch, F.X., A. Lorincz, N. Muñoz, C.J.L.M. Meijer, K.V. Shah (2002) "The causal relation between human papillomavirus and cervical cancer" *J Clin Path* **55**:244-265.
12. Muñoz N, F.X. Bosch, S. de Sanjosé, R. Herrero, X. Castellsagué, K.V. Shah, P.J.F. Snijders, and Chris J.L.M. Meijer, for the International Agency for Research on Cancer Multicenter Cervical Cancer Study Group. (2003) "Epidemiologic Classification of Human Papillomavirus Types Associated with Cervical Cancer" *N Engl J Med* **348(6)**:518-527.
13. Clifford GM, Smith JS, Plummer M, Muñoz N, Franceschi S. Human papillomavirus types in invasive cervical cancer worldwide: a meta-analysis. *Br J Cancer*. 2003; **88**:63-73.
14. Koutsky, L. 1997. Epidemiology of genital human papillomavirus infection. *American Journal of Medicine*. **102(5A)**:3-8.
15. Winer RL, Kiviat NB, Hughes JP, et al. Development and duration of human papillomavirus lesions, after initial infection. *J Infect Dis*. 2005; **191**:731-738.
16. Moscicki, A, Schiffman M, Kjaer S, Villa L. Updating the natural history of HPV and anogenital cancer. *Vaccine* 2006; **24(S3)**; 42-51.
17. Moscicki AB, Ellenberg JH, Farhat S, Xu J. Persistence of human papillomavirus infection in HIV-infected and -uninfected adolescent girls: risk factors and differences, by phylogenetic type. *J Infect Dis*. 2004 Jul 1; **190(1)**:37-45.
18. Palmer Castle PE, Schiffman M, Herrero R, Hildesheim A, Rodriguez AC, Bratti MC, Sherman ME, Wacholder S, Tarone R, Burk RD. A prospective study of age trends in cervical human papillomavirus acquisition and persistence in Guanacaste, Costa Rica. *J Infect Dis*. 2005 Jun 1; **191(11)**:1808-16.
19. Zielinski GD, Snijders PJF, Rozendaal I, et al. High-risk HPV testing in women with borderline and mild dyskaryosis; long term follow-up data and clinical relevance. *J Pathol* 2001; **195**:300-306.
20. Holowaty P, Miller AB, Rohan T, To T. Natural history of dysplasia of the uterine cervix, *J Natl Cancer Inst* 1999; **91**:252-58.
21. Nobbenhuis MA, Helmerhorst TJ, van den Brule AJ, Rozendaal L, Voorhorst FJ, Bezemer PD, Verheijen RH, Meijer CJ. Cytological regression and clearance of high-risk human papillomavirus in women with an abnormal cervical smear. *Lancet*. 2001; **358(9295)**:1782-1783.
22. Cuzick J, Clavel C, Petry KU, et al. Overview of the European and North American studies on HPV testing in primary cervical cancer screening. *Int J Cancer* 2006; **119**:1095-1101.
23. Whitlock EP, Vesco KK, Eder M, Lin JS, Senger CA, Burda BU. Liquid-based cytology and human papillomavirus testing to screen for cervical cancer: a systematic review for the U.S. Preventive Services Task Force. *Ann Intern Med* 2011; **155**:687,697, W214-5.
24. Myers, T.W. and Gelfand, D.H. 1991. Reverse transcription and DNA amplification by a *Thermus thermophilus* DNA polymerase. *Biochemistry*. **30**:7661-7666.
25. Davies, P., Kornegay, J., Iftner, T. 2001. Current methods of testing for human papillomavirus. *Best Practice and Research Clinical Obstetrics and Gynaecology*. Vol. 15, No. 5:677-700.
26. Birch, D.E., et al. 1996. The use of a thermally activated DNA polymerase PCR gives improved specificity, sensitivity and product yield without additives or extra process steps. *Nature*. Vol 381, No 6581:445-446.
27. Meng, Q., et al. 2001. Automated Multiplex Assay System for Simultaneous Detection of Hepatitis B Virus DNA, Hepatitis C Virus RNA, and Human Immunodeficiency Virus Type 1 RNA. *J Clin Microbiol*. Vol 39, No 8:2937-2945.

05641225001-21EN

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

Doc Rev. 19.0

๑.ลงชื่อ.....²⁶.....ประธานกรรมการ

๒.ลงชื่อ.....^{๑๗}.....กรรมการ

๓.ลงชื่อ.....^{๒๕}.....กรรมการ

28. Longo, M.C., Berninger, M.S. and Hartley, J.L. 1990. Use of uracil DNA glycosylase to control carry-over contamination in polymerase chain reactions. *Gene*. **93**:125-128.
29. Higuchi, R., Dollinger, G., Walsh, P.S., and Griffith, R. 1992. Simultaneous amplification and detection of specific DNA sequences. *Bio/Technology* **10**:413-417.
30. Heid, C.A., Stevens, J., Livak, J.K., and Williams, P.M. 1996. Real time quantitative PCR. *Genome Research* **6**:986-994.
31. Richmond, J.Y. and McKinney, R.W. eds. 1999. *Biosafety in Microbiological and Biomedical Laboratories*. HHS Publication Number (CDC) 93-8395.
32. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections*. Approved Guideline-Third Edition. CLSI Document M29-A3 Villanova, PA:CLSI, 2005.
33. International Air Transport Association. *Dangerous Goods Regulations*, 48th Edition. 2007.
34. Wheeler, C.M., Hunt, W.C., Joste, N.E., Key, C.R., Quint, W.G.V. and Castle, P.E. 2009. Human Papillomavirus Genotype Distributions: Implications for Vaccination and Cancer Screening in the United States. *J Natl Cancer Inst*. **101**:475-487.
35. Khan, M.J., Castle, P.E., Lorincz, A.T., Wacholder, S., Sherman, M., Scott, D.R., Rush, B.B., Glass, A.G. and Schiffman, M. 2005. The Elevated 10-Year Risk of Cervical Precancer and Cancer in Women With Human Papillomavirus (HPV) Type 16 and 18 and the Possible Utility of Type-Specific HPV Testing in Clinical Practice. *J Natl Cancer Inst*. **97**:1072-1079.
36. Wright TC, Jr., Massad LS, Dunton CJ, Spitzer M, Wilkinson EJ, and Solomon D. 2007. 2006 Consensus Guidelines for the Management of Women with Abnormal Cervical Cancer Screening Tests, *Am J Obstet Gynecol* **197** (4); 346-355.
37. Saslow D, Runowicz CD, Solomon D, Moscicki A-B, Smith RA, Eyre HJ, Cohen C. 2002, American Cancer Society Guideline for the Early Detection of Cervical Neoplasia and Cancer. *CA Can Jour Clin* **53**: 342-362.
38. Stoler MH, Wright TC Jr, Sharma A, Apple R, Gutekunst K, Wright TL; ATHENA(Addressing THE Need for Advanced HPV Diagnostics) HPV Study Group. High-risk human papillomavirus testing in women with ASC-US cytology: results from the ATHENA HPV study. *Am J Clin Pathol*. 2011 Mar;**135**(3):468-475.

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑.ลงชื่อ.....ประธานกรรมการ

๒.ลงชื่อ.....กรรมการ

๓.ลงชื่อ.....กรรมการ

Document Revision Information	
Doc Rev. 18.0 02/2018	<p>Added cobas® 4800 System - User Assistance.</p> <p>Removed cobas® 4800 System Operator's Manual and cobas® HPV Test Operator's Manual.</p> <p>Added table of result flags.</p> <p>Added RepHresh vaginal hygiene products to Procedural Limitations section and Table 21 footnote.</p> <p>Updated Digital Heater Block with alternative part number.</p> <p>Added stability claim for open bottle of cobas® Sample Prep Buffer (CSPB).</p> <p>Corrected Part number for Roche Cell Collection Medium Replacement Caps.</p> <p>Added a cervical collection device for Roche Cell Collection Medium specimens.</p> <p>Clarified instruction for treatment of SurePath™ specimens.</p> <p>Changed "Tris-HCl buffer" to "Tris buffer" as a reagent component.</p> <p>Please contact your local Roche Representative if you have any questions.</p>
Doc Rev. 19.0 01/2019	<p>Updated hazard warnings.</p> <p>Added Rx Only symbol and description to and updated descriptions of the harmonized symbol page.</p> <p>Please contact your local Roche Representative if you have any questions.</p>

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑.ลงชื่อ.....ประธานกรรมการ

๒.ลงชื่อ.....กรรมการ

๓.ลงชื่อ.....กรรมการ



Roche Molecular Systems, Inc.
1080 US Highway 202 South
Branchburg, NJ 08876 USA
www.roche.com



Roche Diagnostics (Schweiz) AG
Industriestrasse 7
6343 Rotkreuz, Switzerland

Roche Diagnostics
2, Avenue du Vercors
38240 Meylan, France

Distributed by

Roche Diagnostics GmbH
Sandhofer Strasse 116
68305 Mannheim, Germany

Distributore in Italia:
Roche Diagnostics S.p.A.
Viale G. B. Stucchi 110
20052 Monza, Milano, Italy

Roche Diagnostics, SL
Avda. Generalitat, 171-173
E-08174 Sant Cugat del Vallès
Barcelona, Spain

Distribuidor em Portugal:
Roche Sistemas de Diagnósticos Lda.
Estrada Nacional, 249-1
2720-413 Amadora, Portugal

Roche Diagnostics
201, boulevard Armand-Frappier
H7V 4A2 Laval, Québec, Canada
(For Technical Assistance call:
Pour toute assistance technique,
appeler le: 1-877-273-3433)

Roche Diagnostica Brasil Ltda.
Av. Engenheiro Billings, 1729
Jaguará, Building 10
05321-010 São Paulo, SP Brazil

Trademarks and Patents

See <http://www.roche-diagnostics.us/patents>

©2019 Roche Molecular Systems, Inc.
01/2019
Doc Rev. 19.0

05641225001-21



Roche Diagnostics GmbH
Sandhofer Str. 116
68305 Mannheim
Germany




คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑. ลงชื่อ..... ประธานกรรมการ

๒. ลงชื่อ..... กรรมการ

๓. ลงชื่อ..... กรรมการ

The following symbols are now used in labeling for Roche PCR diagnostic products.



Ancillary Software



In Vitro diagnostic medical device



Authorized representative in the European community



Manufacturer



Batch code



Store in the dark



Biological risks



Contains sufficient for <n> tests



Catalogue number



Temperature limit



Consult instructions for use



Test Definition File



Contents of kit



Use-by date



Distributed by



Global Trade Item Number

Rx Only

US Only: Federal law restricts this device to sale by or on the order of a physician.



This product fulfills the requirements of the European Directive 98/79 EC for *in vitro* diagnostic medical devices.

US Customer Technical Support 1-800-526-1247

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์
 ๑. ลงชื่อ..... ประธานกรรมการ
 ๒. ลงชื่อ..... กรรมการ
 ๓. ลงชื่อ..... กรรมการ

PreservCyt® Solution Sample Collection and Transport Medium for use with

HOLOGIC®

The ThinPrep® Pap Test

R_x only

Intended Use for Gynecologic Applications

PreservCyt® Solution is designed for use with the ThinPrep® 2000 System, the ThinPrep 3000 Processor, and the ThinPrep 5000 Processor. PreservCyt Solution is an alcohol-based, preservation solution that serves as a transport, preservative, and antibacterial medium for gynecologic samples.

The PreservCyt Solution component of the ThinPrep 2000 System is an alternative collection and transport medium for gynecologic specimens tested with the Cervista® HPV HR Test, the Cervista® HPV 16/18 Test, the Roche cobas® HPV Test and the Digene Hybrid Capture® System HPV DNA. Refer to the respective manufacturer's package inserts for instructions for using PreservCyt Solution for collection, transport, storage, and preparation of specimens for use in those systems.

The PreservCyt Solution component of the ThinPrep 2000 System is an alternative collection and transport medium for gynecologic specimens tested with the BD ProboTaq™ CT Qx Amplified DNA Assay, Hologic APTIMA COMBO 2^o CT/NG Assays and the Hologic APTIMA® Trichomonas vaginalis Assay. Refer to the respective manufacturer's package inserts for instructions for using PreservCyt Solution for collection, transport, storage, and preparation of specimens for use in those systems.

The PreservCyt Solution component of the ThinPrep 2000 System is also an alternative collection and transport medium for gynecologic specimens tested with the Roche Diagnostics COBAS® AMPLICOR CT/NG assay. Refer to Hologic's labeling (Document # MAN-02063-001) for instructions for using PreservCyt Solution for collection, transport, storage, and preparation of specimens and to the Roche Diagnostics COBAS AMPLICOR CT/NG package insert for instructions for use of that system.

Summary and Explanation of the ThinPrep System

The ThinPrep process begins with the patient's gynecologic sample being collected by the clinician using a cervical sampling device which, rather than being smeared on a microscope slide, is immersed and rinsed in a vial filled with PreservCyt Solution. The ThinPrep sample vial is then capped, labeled, and sent to a laboratory equipped with a ThinPrep Processor.

At the laboratory, the PreservCyt sample vial is placed into a ThinPrep Processor and a gentle dispersion step breaks up blood, mucus, non-diagnostic debris, and thoroughly mixes the cell sample. The cells are then collected on a ThinPrep Pap Test Filter specifically designed to collect diagnostic cells. The ThinPrep Processor constantly monitors the rate of flow through the ThinPrep Pap Test Filter during the collection process in order to prevent the cellular presentation from being too scant or too dense. A thin layer of cells is then transferred to a glass slide in a 20 mm-diameter circle, and the slide is automatically deposited in or sprayed with a fixative solution.

Reagents

For In Vitro Diagnostic Use. Contains buffered methanol. CAS 67-56-1

Warnings



Danger. Flammable. Contains Methanol.

- H301 – Toxic if swallowed.
- H311 – Toxic in contact with skin.
- H331 – Toxic if inhaled.
- H370 – Causes damage to organs.
- H226 – Flammable liquid and vapor.

NOT FOR EXTERNAL OR INTERNAL USE IN HUMANS OR ANIMALS. CANNOT BE MADE NON-POISONOUS. Use with adequate ventilation.

- P210 – Keep away from heat/sparks/open flames/hot surfaces.
- P233 – Keep container tightly closed.
- P264 – Wash hands thoroughly after handling.
- P280 – Wear protective gloves/protective clothing/eye protection/face protection.

First Aid

IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. See www.hologiccsds.com for the entire Safety Data Sheet.

Storage

เงื่อนไขเพิ่มเติมข้อ 4.3

Store PreservCyt Solution between 15°C (59°F) and 30°C (86°F). Do not use beyond the expiration date printed on the container.

Store PreservCyt Solution with cytologic sample intended for ThinPrep Pap testing between 15°C (59°F) and 30°C (86°F) for up to 6 weeks.

Disposal

Dispose in accordance with all applicable regulations.

Specimen Collection for Gynecologic Samples

Collect samples in the routine manner (e.g., CLSI guideline GP-15A3) and rinse in PreservCyt Solution following recommended technique (see reverse for appropriate sample collection technique).

Aliquot Removal (Optional)

Testing for certain sexually transmitted diseases (STD) and for Human Papilloma Virus (HPV) in conjunction with cytology may be performed using the residual specimen remaining in the PreservCyt sample vial after preparation of the ThinPrep Pap Test slide. Such testing may also be enabled by the removal of an aliquot of up to 4 mL (Aliquot Removal) from the PreservCyt sample vial before preparing the ThinPrep Pap Test slide.

Because cytology/HPV testing and STD testing address different clinical questions, Aliquot Removal may not be suitable for all clinical situations. Physicians and other persons responsible for ordering clinical tests should be familiar with the following:

- There is no evidence of degradation of cytology results by Aliquot Removal, however, this cannot be ruled out for all specimens. As with any subsampling step in anatomic pathology, chance misallocation of diagnostic cells may occur if they are very rare. If negative results from the specimen do not fit with the clinical impression, a new specimen may be necessary.
- Aliquot Removal from low-cellularity specimens may leave insufficient material in the PreservCyt sample vial for preparation of a satisfactory ThinPrep Pap Test slide.
- Aliquot Removal may leave insufficient material in the PreservCyt sample vial for performance of ancillary testing (e.g., reflexive HPV testing) using the residual specimen following preparation of a ThinPrep Pap Test slide.
- Co-collection of separate samples for the ThinPrep Pap Test and STD testing may be considered in lieu of Aliquot Removal.
- When opting for concurrent cytologic and STD testing, providers should consider risk and clinical history (e.g., disease prevalence, patient age, sexual history or pregnancy) as well as specimen suitability (e.g., exudates or bleeding) that can impact diagnostic reliability.

Sexually Transmitted Diseases Treatment Guidelines 2002 (Centers for Disease Control and Prevention, MMWR 2002: 51(No. RR-6)) provides clinical guidance for the management and treatment of individual patients, including use of Pap testing.

If ancillary testing is to be performed, laboratory personnel must follow specific instructions (refer to the Operator's Manual for the ThinPrep processor) to appropriately remove the desired aliquot volume and prepare the PreservCyt sample vial for the ThinPrep Pap Test. Adherence to these instructions must be maintained to ensure there is no adverse effect on the ThinPrep Pap Test result.

Materials Required but Not Supplied

- ThinPrep Pap Test Filters (clear) for gynecologic slide preparation
- Microscope slides
- Cover slips
- Staining reagents

คณะกรรมการพิจารณาผลการตรวจตรวจคัดกรองมะเร็งปากช่องคอก่อน

MAN-01328-4001 Rev. 001

๑. ลงชื่อ..... ประธานกรรมการ
 ๒. ลงชื่อ..... กรรมการ
 ๓. ลงชื่อ..... กรรมการ

Specific Performance Characteristics

When used and stored as described, PreservCyt Solution will preserve a cytologic sample for up to 6 weeks.

PreservCyt Solution is bactericidal. PreservCyt Solution has been shown to cause greater than 99.999 percent inactivation within 15 minutes for the following bacteria: *Candida albicans*, *Escherichia coli*, *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Mycobacterium tuberculosis* (determined according to the United States Pharmacopeia preservative antimicrobial effectiveness test, U.S.P. XXII,51). As with all laboratory procedures, universal precautions should be followed.

Limitations of Procedure

Always use good sampling techniques when collecting samples. Poor sampling techniques will produce inadequate samples. Sampling is limited to broom-like and endocervical brush/plastic spatula combination collection devices. Consult sampling instructions below for specific techniques.

Warranty

This product is warranted to perform as described in the labeling and in Hologic's literature. Hologic disclaims any implied warranty of merchantability of fitness for any other purpose and in no event will Hologic be liable for any consequential damages arising out of the aforesaid express warranty. Distributed by Hologic, Inc., 250 Campus Drive, Marlborough, MA 01752 USA, Phone: 1 (800) 442-9892 or (508) 263-2900.




Figure 1

Line on cap passes line on vial

SAMPLE COLLECTION TECHNIQUE USING BRUSH/SPATULA COLLECTION DEVICES:

- STEP 1 Sample ectocervix with a plastic spatula.
- STEP 2 Rinse spatula in the PreservCyt Solution vial by swirling vigorously 10 times. Place cap on vial until step 4. Discard collection device.
- STEP 3 Sample endocervix with an endocervical brush.
- STEP 4 Rinse the brush in the PreservCyt Solution by rotating the device in the solution 10 times while pushing against the PreservCyt vial wall. Swirl the brush vigorously to further release material. Discard the collection device.
- STEP 5 Tighten the PreservCyt sample vial cap so that the torque line on the cap passes the torque line on the vial. See Figure 1.

SAMPLE COLLECTION TECHNIQUE USING BROOM-LIKE COLLECTION DEVICE:

- STEP 1 Obtain a sample from the cervix using a broom-like device.
- STEP 2 Rinse the collection device into a PreservCyt Solution vial by pushing the brush into the bottom of the vial 10 times, forcing the bristles to bend apart to release the cervical material. As a final step, twirl the brush between the thumb and forefinger vigorously to further release cellular material. Discard the collection device.
- STEP 3 Cap the PreservCyt sample vial tightly. Tighten the cap of the vial so that the black torque line on the cap passes the black torque line on the vial. See Figure 1.

Symbols

The following symbols may appear on your product:

Symbol	Title	Description	Standard information
	Use-by date	Indicates the date after which the medical device is not to be used.	ISO 15223-1 Medical devices—Symbols to be used with medical device labeling and information to be supplied, ISO 15223-1:2012(E) Section 5.1.4
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1 Medical devices—Symbols to be used with medical device labeling and information to be supplied, ISO 15223-1:2012(E) Section 5.1.6
	In vitro diagnostic medical device	Indicates a medical device that is intended to be used as an <i>in vitro</i> diagnostic medical device	ISO 15223-1 Medical devices—Symbols to be used with medical device labeling and information to be supplied, ISO 15223-1:2012(E) Section 5.5.1
	Caution, consult instructions for use	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	ISO 15223-1 Medical devices—Symbols to be used with medical device labeling and information to be supplied, ISO 15223-1:2012(E) Section 5.4.4
	Consult instructions for use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1 Medical devices—Symbols to be used with medical device labeling and information to be supplied, ISO 15223-1:2012(E) Section 5.4.3
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1 Medical devices—Symbols to be used with medical device labeling and information to be supplied, ISO 15223-1:2012(E) Section 5.1.5
	Manufacturer	Indicates the medical device manufacturer, as defined in the EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.	ISO 15223-1 Medical devices—Symbols to be used with medical device labeling and information to be supplied, ISO 15223-1:2012(E) Section 5.1.1
	Temperature limitation	Indicates the upper and lower limit of temperature to which the medical device can be safely exposed.	ISO 15223-1 Medical devices—Symbols to be used with medical device labeling and information to be supplied, ISO 15223-1:2012(E) Section 5.3.7
	Acute Toxicity	Toxic if swallowed. Toxic in contact with skin. Toxic if inhaled.	United States Department of Labor Occupational Safety and Health Administration's (OSHA) Hazard Communication Standard (HCS), Appendix C to §1910.1200, Sections C.4.1, C.4.2, C.4.3 (Classified in Accordance with Appendix A.1)
	Flammable	Flammable liquid and vapor	OSHA's HCS, Appendix C to §1910.1200, Section C.4.19 (Classified in Accordance with Appendix A.1)
	Respiratory Sensitizer, Target Organ Toxicity	Causes damage to organs	OSHA's HCS, Appendix C to §1910.1200, Section C.4.11 (Classified in Accordance with Appendix A.1)

๑. ลงชื่อ.....ประธานกรรมการ

๒. ลงชื่อ.....กรรมการ

๓. ลงชื่อ.....กรรมการ



คุณลักษณะเฉพาะข้อ 2.3

Title: Technical Data Sheet Rovers[®] Cervex-Brush[®] (Sterile)
 (Reference: 380100331)(Roche P/N for ordering 08779040190)

Date: September 13th, 2019
 Doc. ID: 2019256C
 Version: 1.0

Author / Approved:


H. Vissers,
 Quality Assurance / Regulatory Affairs Manager


Signature:

THE INFORMATION CONTAINED HEREIN IS CONFIDENTIAL AND THE PROPERTY OF ROVERS MEDICAL DEVICES B.V., LEKSTRAAT 10, 5347 KV, OSS, THE NETHERLANDS. THE RECIPIENT IS RESPONSIBLE FOR SAFE-KEEPING, AND PREVENTION OF UNAUTHORIZED APPROPRIATION, USE, DISCLOSURE, AND COPYING.

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑.ลงชื่อ..... ประธานกรรมการ

๒.ลงชื่อ..... กรรมการ 1/10

๓.ลงชื่อ..... กรรมการ



1 Introduction

Rovers® Cell Sampling Devices (Sterile) for professional and home use are based on "thin-hair technology". Cell material collected with Rovers® Cell Sampling Devices can be processed with conventional as well as liquid based cytology(LBC) and also be used for bacteriological testing, virus testing (e.g. for Sexually Transmitted Diseases or HPV) and DNA analysis.

Cell collection with Rovers® Cell Sampling Devices is quick, simple, painless and reliable.

In accordance with the European Medical Device Directive 93/42/EEC (as amended by Directive 2007/47/EC) the Rovers® Cell Sampling Devices (Sterile) are Class Is Medical Devices with a low risk profile.

คณะกรรมการพิจารณาผลการประกวดราคาอิเล็กทรอนิกส์

๑.ลงชื่อ.....ประธานกรรมการ

๒.ลงชื่อ.....กรรมการ

๓.ลงชื่อ.....กรรมการ ๒/10