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Models L110, L111

1.17

- Provides an ImageReady™ MR Conditional pacing system* at 3T and 1.5T, full body, with no time limitations, with automatic MRI timeout feature to optimize workflow in the MR environment
- RF telemetry for wireless transmission of information and efficiency in the operating room and follow-up setting
- LATITUDE™ NXT Remote Patient Management enabled, offering the opportunity for wireless (RF) remote patient monitoring and follow-up
- PaceSafe™ RV and RA, providing dynamic adjustment of pacing outputs to ensure capture, to maximize efficiency and ease of use
- RightRate™ – MV sensor technology and the only MV sensor clinically proven to restore chronotropic competence¹
- AVSH+, designed to minimize unnecessary RV pacing without clinically significant pauses, therefore reducing the risk of HF development
- Enhanced features and diagnostics designed to provide you with greater insight into your patient's disease progression
- POST function to facilitate patient follow up with a fully automatic device and lead check
- EASYVIEW™ header with port identifiers designed to make the implant experience more efficient

*Please refer to MRI Technical Guide.

Mechanical Specifications

Model	Type	Size (cm) (W x H x D)	Mass (g)	Volume (cc)	Connector Type (RA RV LV)
L110	SR	4.45 x 4.81 x 0.75	23.6	13.2	RA/RV: IS1
L111	DR	4.45 x 5.02 x 0.75	24.8	13.7	RA: IS1 – RV: IS1

1.18

Projected Longevity (Years)

Pacing		SR	DR
50%	RA/RV 2.5V	10.4	9.3
100%	RA/RV 2.5V	9.7	8.2

Additional Longevity Information

- Settings: pacing pulse width 0.4ms, Impedance 750Ω, LRL 60bpm, Sensor On, EGM Onset On. These calculations also assume that the pulse generator spends 6 months in Storage mode during shipping and storage, the ZIP™ telemetry use for 1 hour at implant time and for 40 minutes annually for in-clinic follow-up checks. For longevity calculations based on different settings please contact Boston Scientific technical services or your local representative.
- Power Supply SR and DR models: lithium-carbon monofluoride cell; Boston Scientific; 402290.

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Models L110, L111

Pacing Therapy

Brady Modes	Normal DDD(RI-DDI(RI-VDD(RI-VVI(RI-AAI(RI-DOO-VOO-AGO-Off) Temporary: DDD-DDI-VDD-VVI-AAI-DOO-VOO-AGO-Off	1.3
AT/AF Management	ATR Mode Switch, Rate Smoothing	
Automaticity	Automatic Gain Control (AGC) for sensitivity Right Atrial Automatic Threshold (RAAT) Right Ventricular Automatic Capture (RVAC)	1.8, 1.9, 1.15
Rate Adaptive Pacing	Accelerometer, RightRate™ (Minute Ventilation) or Blended sensors with sensor trending function	1.1, 1.16
RV Pacing Reduction	AV Search +, AV Delay to 400 ms, RateHysteresis	1.13
Rate Management	Sudden Brady Response (SBR), PMT Termination, PVARP after PVC, Dynamic PVARP	
Pace/Sense Configuration	Unipolar, Bipolar, Bipolar/Unipolar, Unipolar/Bipolar, Unipolar/Off, Bipolar/Off, Lead Safety Switch	1.11, 1.12

Patient Diagnostics

Arrhythmia Logbook	Event Summary, Stored Electrograms with Annotation Markers (Intervals and approximately 14 minutes all multi-channel EGM, always with 10 seconds Onset and event storage prioritization), Implant activation of all available EGMs. On screen measurements of all stored signal, amplitudes and timing. Snapshot Function (up to 12 seconds trace of ECG/EGM display stored)
Histograms & Counters	Ventricular Tachy Counter, Brady Counter, Histograms, Intrinsic Promotion (Rate Hysteresis % successful and AVSH+ % successful)
Diagnostics	AT/AF Burden, A & V Arrhythmias, Weight and Blood Pressure*
DAILY TREND for last 365 Days	Events, AT/AF Burden, Heart Rate, Lead Impedance and Amplitude, RAAT Trend, RVAC Trend

*Weight and Blood Pressure are only available via LATITUDE NXT.

ImageReady™

MRI Lead Selection	Pulse Generator MR-conditional with all FINELINE™ II Sterox, FINELINE™ II Sterox EZ and INGEVITY™ Pacing Lead Models
MRI Conditions	Full body scan at 1.5T (≤SAR 2W/Kg) for all FINELINE™ II models* Full body scan at 3T and 1.5T (≤SAR 4W/Kg) for all INGEVITY™ MRI models*
MRI Mode	Pacing Mode: AOO,VOO,DOO,Off Protection Time Out: Off, 12,24,48 hours

*Please refer to the Pacing System MRI Technical Guide as the system is designated as MR Conditional in accordance with specified conditions.

Implant/In Clinic Follow Up

Implant Communication Mode	Programmable values: Enable use of ZIP™ telemetry (MICS) (Requires initial use of wand for device ID) or use wand for all telemetry Nominal: Enable use of ZIP™ telemetry (Requires initial use of wand for device ID)
In Clinic Follow Up	Snapshot Function up to 12 seconds trace of ECG/EGM display stored POST (Post-Operative System Test), provides an automatic device/lead check at a pre-determined time post-implant to help document proper system functionality without requiring manual system testing

Remote Follow Up

Remote Monitoring	This device is designed to be LATITUDE™ NXT enabled, LATITUDE NXT availability varies by region*
Thresholds	Automatic storage of last successful daily PaceSafe threshold test for all active chambers
Wireless	Remote follow-up for all devices (MICS)
Patient Triggered Monitor (PTM)	Triggers the storage of two minutes onset and one minute post – EGMs, intervals, and annotated marker data during a symptomatic episode by placing a magnet over the device

*LATITUDE™ NXT is not available for L100 and L101 models.

Safety Functions*

Safety Core	Is intended to provide life-sustaining therapy if certain non-recoverable or repeat fault conditions occur. Safety Core operates independently and acts as a backup to these components
Electrocautery Protection Mode	Provides asynchronous pacing at the programmed outputs and LRL when commanded by the programmer

*The Safety Functions do not have programmable parameters.

1. Chronotropic competence is defined by the Model of the Cardiac Chronotropic Response to Exercise, Willif B, Corey J, Blackburn G. A mathematical model of the cardiac chronotropic response to exercise. Journal of Electrophysiology 1969;3:176-180. Refer to the Physician's System Guide for more information on adaptive-rate therapy. Additional clinical performance was assessed using INSIGNIA™ Ultra clinical data with the AutoLifestyle™ feature programmed On. Boston Scientific. Data on file. ALTRUA™ Pacemaker System Guide, 2008.1.20-29, monthly Full Interrogations (scheduled remote follow-ups, and quarterly patient-interrogations)

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Guidant Corporation and Guidant Europe (wise are Boston Scientific companies. Information contained herein is for distribution outside the U.S. only. Information not intended for distribution in France.

CRM-251101-AA JUL 2014 Printed in the Netherlands by Gosling

Boston Scientific
Advancing science for life™

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
CAUTION: Federal law (USA)
restricts this device to sale by
or on the order of a physician
trained or experienced in
device implant and follow-up
procedures.

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a. If the Communication Mode is selected via the Utilities button on the PRM Startup screen, the Nominal setting within the ZOOMVIEW Programmer software application will correspond to the value chosen on the Startup screen.

Table A-2. Device Mode

Parameter	Programmable Values	Nominal
Device Mode	Exit Storage; Enable Electrocautery Protection; Enable MRI Protection ^a	Storage

a. Available in models with the MRI Protection Mode feature.

Table A-3. Pacing therapy parameters (specified into a 750 Ω load)

Parameter	Programmable Values	Nominal
Mode ^{a c}	DDD(R); DDI(R); DOO; VDD(R); VVI(R); VOO; AAI(R); AOO; Off; Temporary: DDD; DDI; DOO; VDD; VVI; VOO; AAI; AOO; Off (1.3)	Dual Chamber: DDD; Single Chamber: VVI
Lower Rate Limit (LRL) ^{a b c} (ppm)	30; 35; ...; 185 } (1.5)	60 (Tolerance \pm 5 ms)
Maximum Tracking Rate (MTR) ^{a c} (ppm)	50; 55; ...; 185 }	130 (Tolerance \pm 5 ms)
Maximum Sensor Rate (MSR) ^e (ppm)	50; 55; ...; 185	130 (Tolerance \pm 5 ms)
Pulse Amplitude ^{a c d} (dual chamber, atrium) (V)	Auto; 0.1; 0.2; ...; 3.5; 4.0; ...; 5.0; Temporary: 0.1; 0.2; ...; 3.5; 4.0; ...; 5.0 (1.6)	3.5 (Tolerance \pm 15% or 100 mV, whichever is greater)
Pulse Amplitude ^{a c d} (dual chamber, right ventricle) (V)	Auto; 0.1; 0.2; ...; 3.5; 4.0; ...; 7.5; Temporary: 0.1; 0.2; ...; 3.5; 4.0; ...; 7.5 (1.6)	3.5 (Tolerance \pm 15% or 100 mV, whichever is greater)
Pulse Amplitude ^{a c d} (single chamber) (V)	Auto; 0.1; 0.2; ...; 3.5; 4.0; ...; 7.5; Temporary: 0.1; 0.2; ...; 3.5; 4.0; ...; 7.5	3.5 (Tolerance \pm 15% or 100 mV, whichever is greater)
Pulse Amplitude Daily Trend ^f (independently programmable in each chamber that has the Pacesafe feature)	Disabled; Enabled	Enabled (ACCOLADE, PROPONENT, ESSENTIO, and ALTRUA 2 devices) Disabled (FORMIO, VITALIO, INGENIO, and ADVANTIO devices)
Pulse Width ^{a c d g} (atrium, right ventricle) (ms)	0.1; 0.2; ...; 2.0 (1.7)	0.4 (Tolerance \pm 0.03 ms at $<$ 1.8 ms; \pm 0.08 ms at \geq 1.8 ms)
Accelerometer ^h	On; Passive	Passive
Accelerometer Activity Threshold	Very Low; Low; Medium Low; Medium; Medium High; High; Very High	Medium
Accelerometer Reaction Time (sec)	10; 20; ...; 50	30
Accelerometer Response Factor	1; 2; ...; 16	8
Accelerometer Recovery Time (min)	2; 3; ...; 16	2
Minute Ventilation ^h	On; Passive; Off	Passive

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Table A-3. Pacing therapy parameters (specified into a 750 Ω load) (continued)

Parameter	Programmable Values	Nominal
Minute Ventilation Response Factor	1; 2; ...; 16	8
Minute Ventilation Fitness Level	Sedentary; Active; Athletic; Endurance Sports	Active
Patient's Age ^h	≤ 5; 6-10; 11-15; ...; 91-95; ≥ 96	56-60
Patient's Gender ^h	Male; Female	Male
Ventilatory Threshold (ppm)	30; 35; ...; 185	120 (Tolerance ± 5 ms)
Ventilatory Threshold Response (%)	Off; 85; 70; 55	70
Rate Hysteresis Hysteresis Offset ^e (ppm)	-80; -75; ...; -5; Off	Off (Tolerance ± 5 ms)
Rate Hysteresis Search Hysteresis ^a (cycles)	Off; 256; 512; 1024; 2048; 4096	Off (Tolerance ± 1 cycle)
Rate Smoothing (Up, Down) ^e (%)	Off; 3; 6; 9; 12; 15; 18; 21; 25	Off (Tolerance ± 1%)
Rate Smoothing Maximum Pacing Rate (ppm)	50; 55; ...; 185	130 (Tolerance ± 5 ms)
Sudden Brady Response (SBR) ^e	Off; On	Off
SBR Atrial Paces Before Therapy	1; 2; ...; 8	3
SBR Atrial Pacing Rate Increase (ppm)	5; 10; ...; 40	20
SBR Therapy Duration (min)	1; 2; ...; 15	2
SBR Inhibit During Rest	Off; On	On
Atrial Pace/Sense Configuration ^{a c} (dual chamber)	Unipolar; Bipolar; Bipolar/Unipolar; Unipolar/Bipolar; Unipolar/Off; Bipolar/Off (1.11), (1.12)	Bipolar
Right Ventricle Pace/Sense Configuration ^{a c} (dual chamber)	Unipolar; Bipolar; Bipolar/Unipolar; Unipolar/Bipolar (1.11), (1.12)	Bipolar
Pace/Sense Configuration ^{a c} (single chamber)	Unipolar; Bipolar; Bipolar/Unipolar; Unipolar/Bipolar	Bipolar
Safety Switch (independently programmable in each chamber)	Off; On	On
Automatic Lead Recognition	Off; On	On
Maximum Paced AV Delay ^{a c} (ms)	30; 40; ...; 400	180 (Tolerance ± 5 ms)
Minimum Paced AV Delay ^{a c} (ms)	30; 40; ...; 400	80 (Tolerance ± 5 ms)
Maximum Sensed AV Delay ^{a c} (ms)	30; 40; ...; 400	150 (Tolerance ± 5 ms)
Minimum Sensed AV Delay ^{a c} (ms)	30; 40; ...; 400	65 (Tolerance ± 5 ms)
AV Search + ^e	Off; On	Off
AV Search + Search AV Delay (ms)	30; 40; ...; 400	300 (Tolerance ± 5 ms)
AV Search + Search Interval (cycles)	32; 64; 128; 256; 512; 1024	32 (Tolerance ± 1 cycle)
RYTHMIQ ^a	AAI(R) with VVI Backup; Off	Off
Maximum A-Refractory (PVARP) ^{a c} (dual chamber) (ms)	150; 160; ...; 500	280 (Tolerance ± 5 ms)
Minimum A-Refractory (PVARP) ^{a c} (dual chamber) (ms)	150; 160; ...; 500	240 (Tolerance ± 5 ms)
Maximum V-Refractory (VRP) ^{a c} (dual chamber) (ms)	150; 160; ...; 500	250 (Tolerance ± 5 ms)

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Owner	manucd	gmpno	catno	offname	pdtname	desc	pageno	umdn	gmdn	RefItemNo
51119	44072		6501	Wrench Kit	Wrench Kit	Wrench Kit IMPLANTABLE LEADS and ACCESSORIES	4	12911		6344072001335
51119	44072		7840	Ingevity+ (Ingevity Plus) Pacing Lead	Ingevity+ (Ingevity Plus) Pacing Lead	Ingevity+ (Ingevity Plus) Pacing Lead IMPLANTABLE LEADS and ACCESSORIES	4	11458		6344072001336
51119	44072		7841	Ingevity+ (Ingevity Plus) Pacing Lead	Ingevity+ (Ingevity Plus) Pacing Lead	Ingevity+ (Ingevity Plus) Pacing Lead IMPLANTABLE LEADS and ACCESSORIES	4	11458		6344072001337
51119	44072		7842	Ingevity+ (Ingevity Plus) Pacing Lead	Ingevity+ (Ingevity Plus) Pacing Lead	Ingevity+ (Ingevity Plus) Pacing Lead IMPLANTABLE LEADS and ACCESSORIES	4	11458		6344072001338
51119	44072		7731	Ingevity MRI	Ingevity MRI	Ingevity MRI IMPLANTABLE LEADS and ACCESSORIES	4	11458		6344072001339
51119	44072		7732	Ingevity MRI	Ingevity MRI	Ingevity MRI IMPLANTABLE LEADS and ACCESSORIES	4	11458		6344072001340
51119	44072		7735	Ingevity MRI	Ingevity MRI	Ingevity MRI IMPLANTABLE LEADS and ACCESSORIES	4	11458		6344072001341
51119	44072		7736	Ingevity MRI	Ingevity MRI	Ingevity MRI IMPLANTABLE LEADS and ACCESSORIES	4	11458		6344072001342
51119	44072		7740	Ingevity MRI	Ingevity MRI	Ingevity MRI IMPLANTABLE LEADS and ACCESSORIES	4	11458		6344072001343
51119	44072		7741	Ingevity MRI	Ingevity MRI	Ingevity MRI IMPLANTABLE LEADS and ACCESSORIES	4	11458		6344072001344
51119	44072		7742	Ingevity MRI	Ingevity MRI	Ingevity MRI IMPLANTABLE LEADS and ACCESSORIES	4	11458		6344072001345

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INGEVITY™ + Pacing Lead


Active Fixation Models: 7840, 7841, 7842

The INGEVITY+ pacing leads are 6F (2.0 mm) steroid-eluting, endocardial pace/sense leads designed for permanent implantation for either atrial or ventricular applications.

INGEVITY+ is built on the proven INGEVITY platform, with nearly 700,000 INGEVITY leads sold worldwide with a 99.2% reliability at 7 years.¹

INGEVITY+ is specifically designed with three layers of insulation between conductors and a polyurethane lead body. The tri-filar inner coil design provides consistent, low, and repeatable turn counts when extending and retracting the helix².

These leads utilize an IS-1 bipolar connector. The tip features a flexible neck design and incorporates an IROX™ (iridium oxide) coating on the tip electrode.

Lead Specifications and Reimbursement Information

Product	INGEVITY+ Pacing Lead
Model/Length	7840 / 45 cm 7841 / 52 cm 7842 / 59 cm } 2.5
Type	Bipolar Atrial / Ventricular Straight
Connector	IS-1 BI } 2.4
Compatibility	Pulse generators with an IS-1 port, which accepts an IS-1 terminal
MRI Conditions of Use*	ImageReady™ MR-Conditional System when used with an MR-Conditional pulse generator - Full body scan 1.5T and 3T } 2.7
Introducer without guide wire	6F (2.0 mm) } 2.6
Introducer with guide wire	9F (3.0 mm)
Fixation	Extendable/retractable helix
Expected number of rotations to fully extend/retract the helix**	6 ± 2 turns with straight stylet 7 ± 3 turns with J stylet
Recommended maximum number of turns to extend/retract the helix**	30
Nominal fixation helix penetration depth	1.8 mm

Q3 2019 Boston Scientific Corporation Product Performance Report
Internal data on file.

* Refer to the MRI Technical Guide for a complete list of cardiology and radiology conditions of use.

** Use fluoroscopy markers for verification of full extension/retraction of the helix. The number of turns to extend or retract the helix may vary based on patient anatomy and implant conditions.

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INGEVITY™ + Pacing Lead

Active Fixation Models: 7840, 7841, 7842

Lead Specifications and Reimbursement Information (continued)

Product	INGEVITY+ Pacing Lead
Nominal Electrode: Fixation helix surface area Distance between electrodes Anode electrode surface area	4.5mm ² 10.7mm 20mm ²
Nominal Diameter: Insertion Anode electrode Lead body Fixation helix	2.0mm (6F) 2.0mm 1.9mm 1.2mm
Material: External insulation Internal insulation Terminal ring contact IS-1 terminal pin contact Tip electrode Anode electrode	Polyurethane (55D) Silicone rubber 316L stainless steel 316L stainless steel IROX™ (iridium oxide) coated Pt-Ir IROX (iridium oxide) coated Pt-Ir
Conductor Type	Tri-filar inner coil of MP35N™ and single-filar outer coil of MP35N with a silver core. ¹
Steroid	0.91 mg dexamethasone acetate (2.1)
Radiopaque Markers	Pt-Ir
Suture Sleeve	Radiopaque white silicone rubber
C-code	1898

¹MP35N is a trademark of SPS Technologies, Inc.

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
Active Fixation Models: 7840, 7841, 7842

Features

Lifetime Warranty: The INGEVITY+ pacing lead family is backed with a lifetime warranty.*

Lead Body Design: The isodiametric lead body consists of a coaxial design that includes a tri-filar inner coil and a single-filar outer coil. Both the inner and outer coils are designed for MR Conditional use in the MRI environment and provide robust flexural fatigue performance. In addition, the tri-filar inner coil provides consistent helix deployment performance. The conductors are separated by both a silicone rubber and Polytetrafluoroethylene (PTFE) lining. The outer coil is covered in Ethylene tetrafluoroethylene (ETFE) for extra insulation protection. The entire lead body is encompassed in a polyurethane outer insulation.

IROX™-coated Electrodes: The electrodes are coated with IROX to increase the microscopic surface area.

Steroid-eluting: Upon exposure to body fluids, the steroid elutes from the lead to help reduce tissue inflammation response at the distal electrode. The steroid suppresses the inflammatory response believed to cause threshold rises typically associated with implanted pacing electrodes.

2.1

Radiopaque Suture Sleeve: The radiopaque suture sleeve is visible under fluoroscopy and is used to secure, immobilize, and protect the lead at the venous entry site after lead placement. The window feature is designed to aid compression of the sleeve onto the lead during suturing.

*Limited lifetime warranty. For a full and complete description of the INGEVITY™+ warranty, please review the warranty card included with the product labeling.

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2.3

Active Fixation Features

Extendable / Retractable Fixation: The extendable/retractable helix design anchors the distal tip electrode to the endocardial surface without support of trabecular structures, offering various lead placement possibilities for the tip electrode in the right atrium and/or right ventricle. The helix serves as the cathode for endocardial pacing and sensing. The lead is designed with a tri-filar inner coil for consistent and repeatable turn counts when extending and retracting the helix. The helix is extended and retracted using the fixation tool.

Mapping: The lead helix is electrically conductive to allow mapping (measuring pacing and sensing thresholds) of potential electrode positions without extending the helix into the tissue. Mapping prior to lead fixation is recommended as it can reduce the potential need for multiple lead positionings.

Fluoroscopic Markers: radiopaque markers near the distal tip can be seen under fluoroscopy. These markers show when the helix is fully retracted or fully extended.

Packaged Accessories

- Vein Pick
- Fixation Tool
- Stylet Guide
- Stylets:

	Pre-loaded	Packaged
7840	45cm soft, long tapered	45cm soft, long tapered 45cm extra soft, tapered 45cm soft, atrial J 45cm soft, wide atrial J
7841	52cm soft, long tapered	52cm soft, long tapered 52cm extra soft, tapered 52cm soft, atrial J 52cm soft, wide atrial J
7842	59cm soft, long tapered	59cm soft, long tapered 59cm extra soft, tapered

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INGEVITY™ + and INGEVITY™ MRI Extendable/Retractable Fixation and Tined Fixation Pacing Leads

INDICATIONS

This Boston Scientific lead is indicated for use as follows:

(2.1), (2.2)

- intended for chronic pacing and sensing in the right atrium and/or right ventricle when used with a compatible pulse generator (INGEVITY+ and INGEVITY MRI extendable/retractable fixation)
- intended for chronic pacing and sensing in the right atrium (Bifurcated Atrial) or right ventricle (Straight) when used with a compatible pulse generator (INGEVITY MRI tined fixation)

CONTRAINDICATIONS

Use of these leads are contraindicated for the following patients:

- Patients with a hypersensitivity to a nominal single dose of 0.25mg dexamethasone acetate (for INGEVITY+ and INGEVITY MRI extendable/retractable fixation)
- Patients with a hypersensitivity to a nominal single dose of 0.5mg dexamethasone (for INGEVITY MRI tined fixation)
- Patients with mechanical tricuspid heart valves.

WARNINGS

Read the manual thoroughly before replacement to avoid damage to the pulse generator and/or lead. For single patient use only. Do not reuse, reprocess, or recondition. Always have external defibrillation equipment available during implant and electrophysiology testing. Ensure that an external defibrillator and medical personnel skilled in CPR are present during post-implant device testing should the patient require external rescue. Lead fracture, dislodgment, abrasion, or an incomplete connection constitutes a periodic or sustained loss of pacing or sensing or both. Although possible, the lead is not designed to tolerate excessive flexing, bending, or torsion. Do not kink, twist, or bend the lead with other leads. Implant of the system cannot be performed in an MRI site (Zone II) and higher. Take care to obtain appropriate electrode position. Failure to do so may result in suboptimal lead measurements. Unless all of the MRI Conditions of Use (as described in the MRI Technical Guide) are met, MRI scanning of the patient does not meet MRI Conditional requirements of the implanted system. Refer to the MRI Technical Guide for potential adverse events applicable when Conditions of Use are met or not met, as well as a complete list of MRI related Warnings and Precautions. Do not subject a patient with an implanted pulse generator and/or lead to diathermy.

PRECAUTIONS

For specific information on precautions, refer to the following sections of the product labeling: clinical considerations, sterilization and storage, handling, implantation, hospital and medical environments, and follow-up testing.

POTENTIAL ADVERSE EVENTS

Based on the literature and on pulse generator and/or lead implant experience, the following alphabetical list includes the possible adverse events associated with implantation of products described in this literature: An embolism, Allergic reaction, Arterial damage with subsequent stenosis, Bleeding, Bradycardia, Breakage/failure of the implant instrument, Cardiac perforation, Cardiac tamponade, Chiasm nerve damage, Component failure, Conductor coil fracture, Death, Electrolyte imbalance/dehydration, Elevated thresholds, Erosion, Excessive fibrotic tissue growth, Extracardiac stimulation (muscle/nerve stimulation), Fluid accumulation, Foreign body rejection phenomena, Formation of hematoma or seroma, Heart block, Hemorrhage, Hemothorax, Inability to pace, Inappropriate therapy (e.g., shocks and anti-tachycardia pacing (ATP) where applicable, pacing), Incisional pain, Incomplete lead connection with pulse generator, Infection including endocarditis, Lead dislodgment, Lead fracture, Lead insulation breakage or abrasion, Lead tip deformation and/or breakage, Malnutrition or skin burn due to fluoroscopic radiation, Myocardial trauma (e.g., tissue damage, valve damage), Myopotential sensing, Overensing/undersensing, Pericardial rub, effusion, Pneumothorax, Pulse generator and/or lead migration, Symploc, Tachyarrhythmias, which include acceleration of arrhythmias and early, recurrent atrial fibrillation, Thrombosis/thromboemboli, Valve damage, Vaso-vagal response, Venous occlusion, Venous trauma (e.g., perforation, dissection, erosion).

For a list of potential adverse events associated with MRI scanning, refer to the ImageReady™ MR Conditional Pacing System or Defibrillation System MRI Technical Guide.

4986202 (Rev B)

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Directions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

Boston Scientific

Advancing science for life™

Rhythm Management
 300 Boston Scientific Way
 Marlborough, MA 01752-1234
www.bostonscientific.com

Medical Professionals:
 1 800.CARDIAC (227.3422)
 Patients and Families:
 1 866.481.5268

© 2020 Boston Scientific Corporation or its affiliates. All rights reserved.

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
transmedic

Boston Scientific

The instructions for use (IFU) for this product are supplied in electronic form only. You have the responsibility of ensuring you have the correct version of the IFU. Please check the version (V) and the date to have the product label attached for reference.

Contact Information
Phone: +66-2-022-1899
Fax: +66-2-022-1898
Email: ASE.ANCOMPLANT@boston.com


Imported by Boston Scientific (Thailand) Limited
88 South Square Office Tower, 28th Floor, Unit 2807, 2811
North Sathorn Road, Sathon, Bangkok, 10500

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Ø 0.9 mm (insertion diameter)
Insertion site (mm)
Insertion site (mm)
Insertion site (mm)

CAUTION: Store at 20°-25°C (68°-77°F). Excursions are permitted between 0°C to 30°C (32°F to 86°F). Temperature spikes are permitted up to 30°C (86°F).
Please refer to the Boston Scientific website for more information.
Made in USA.

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
Manufactured at:
Goldstar Puerto Rico B.V., a subsidiary of
Boston Scientific Corporation
No. 12, Road 698
Dorado, Puerto Rico 00646

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TRANSMEDIC (THAILAND) CO.,LTD.
TRANSMEDIC (THAILAND) CO., LTD
Transmedic (S) Pte. Ltd. • 1111 Song • Philippines • Britain