

Clinically proven

Highly deliverable

Ultrathin 60 μ m
struts


Technical data /
ordering info

^{1,2}
Vascular Intervention // Coronary
Drug-Eluting Stent System

BIOTRONIK
excellence for life

Orsiro

Pongseth

Pongseth

Clinically proven

Extensive clinical program*


- >32,500 patients enrolled
- >50,500 patients planned in total
- >44 studies ongoing
- >55 studies planned in total

*status as of Feb 2017


Outstanding clinical results even in challenging subgroups

Orsiro has demonstrated consistently low target lesion failure (TLF) in all-comers trials compared to major modern drug-eluting stents (DES).


BIO-RESORT^{1,2} (n=3,514 patients)

BIO-SCIENCE^{3,4} (n=2,121 patients)

Sort-out VII^{5,6} (n=2,314 patients)

0 5 10
TLF at 12 months (%)

STEMI
5.4%
10.8% Xience
TLF* 24months
BIO-SCIENCE*

Diabetics
0.0%
ST* 60months
BIOFLOW-II*

Small vessels
0.0%
ST* 60months
BIOFLOW-II*

ST - Stent Thrombosis

>50,500
patients planned
in total

Pangatorn


Pangatorn

The new benchmark for DES

BIOFLOW-V 12-month clinical outcomes compared to Xience

In a post-hoc analysis of pooled patient-level data from three RCTs, Orsiro achieved a 96.9% probability of superiority* on TLF rate versus Xience.¹⁰

BIOFLOW-V / -IV / -II Bayesian Population (n=2,208)



*Posterior probability. Bayesian analytical methods were applied

Proven long term clinical outcomes


All stents implanted from 2007 until January 11, 2017
unadjusted (SCAAR)^{11,12}

Orsiro showed a lower restenosis rate than all DES out to five years.

3.8%
restenosis rate at five years



Highly deliverable

Designed for challenging cases, the Orsiro stent system provides better push and easier cross with a lower crossing profile.


Better push

Transmitting up to 57%¹² more force from hub to tip.¹³


Easier cross

Up to 68% less force^{13,14} needed to successfully cross demanding anatomies.

Lower crossing profile

Improved acute performance - up to 13% lower crossing profile.¹⁵

Ultrathin 60 µm struts

Thinner struts make the difference

Thinner struts create:

- Less disrupted flow¹⁸
- Less arterial injury¹⁸

Which leads to:

- Improved re-endothelialization¹⁸
- Reduced risk of restenosis and thrombosis¹⁸


The thinner the better, as long as the radial force can be maintained¹⁸

Up to 15% more radial strength^{19,20} for stronger scaffolding once implanted.


Orsiro
BIOTRONIK

Synergy
Boston Scientific

Xience Xpedition
Abbott

15%
more radial
strength

Strut thickness
in perspective¹⁸

Orsiro
BIOTRONIK
CoCr-SES

60 µm*

Synergy
Boston Scientific
PtCr-EES

74 µm

Ultimaster
Terumo
CoCr-SES

80 µm

Resolute Onyx
Medtronic
CoNi-ZES

81 µm

Xience Family
Abbott
CoCr/EES

81 µm

Promus
Boston Scientific
PtCr-EES

81 µm

BioMatrix
Biosensors
316L-BES

120 µm

* ø 2.25 - 3.0 mm

Orsiro

Vascular
Intervention
Coronary

Indicated for discrete de novo stenotic lesions and in-stent restenotic lesions.*

Technical Data		Stent
Stent material		Cobalt chromium, L-605 (2.1)
Passive coating		proBIO (Amorphous Silicon Carbide)
Active coating		BIOlute bioabsorbable Poly-L-Lactide (PLLA) eluting a limus drug (2.2)
Drug dose		1.4 µg/mm
Strut thickness		ø 2.25 - 3.0 mm: 60 µm (0.0024") (2.4) ø 3.50 - 4.0 mm: 80 µm (0.0031")
Delivery system		
Catheter type		Rapid exchange
Recommended guide catheter		5F (max. I.D. 0.035")
Lesion entry profile		0.017"
Guide wire diameter		0.016"
Usable catheter length		140 cm (2.5)
Balloon material		Semi-crystalline polymer material
Coating (distal shaft)		Hydrophilic coating
Marker bands		Two swaged platinum-iridium markers
Proximal shaft diameter		2.0F
Distal shaft diameter		2.6F ø 2.25 - 3.5 mm; 2.3F ø 4.0 mm
Nominal pressure (NP)		8 atm
Rated burst pressure (RBP)		16 atm

(1.2) Compliance Chart

		Balloon diameter x length (mm)					
		ø 2.25 x 9-40	ø 2.50 x 9-40	ø 2.75 x 9-40	ø 3.00 x 9-40	ø 3.50 x 9-40	ø 4.00 x 9-40
Nominal Pressure (NP)	atm**	8	8	8	8	8	8
	ø (mm)	2.25	2.50	2.75	3.00	3.50	4.00
Rated Burst Pressure (RBP)	atm**	16	16	16	16	16	16
	ø (mm)	2.50	2.77	3.05	3.33	3.88	4.44

**1 atm = 1.013 bar

Ordering Information

Stent ø (mm)	Catheter length 140 cm								
	9	13	15	18	22	26	30	35	40
(2.3) 2.25	364469	364475	364481	364487	364499	364505	364511	391234	391238
2.50	364470	364476	364482	364488	364500	364506	364512	391235	391239
2.75	364471	364477	364483	364489	364501	364507	364513	391236	391240
3.00	364472	364478	364484	364490	364502	364508	364514	391237	391241
3.50	364473	364479	364485	364491	364503	364509	364515	391018	391020
4.00	364474	364480	364486	364492	364504	364510	364516	391019	391021

1. van Boven et al. Very thin strut bioabsorbable polymer everolimus-eluting stents versus durable polymer zotarolimus-eluting stents in all-comers with coronary artery disease (BIO-RESORT), a three-arm, randomised, non-inferiority trial. The Lancet 2016; 10.1016/S0140-6736(16)01925-1 and presentation at ICT 2016; 2. FLF as a composite of cardiac death, target vessel-related myocardial infarction, or clinically indicated target lesion revascularization; 3. Pijnen et al. Ultrathin strut biodegradable polymer sirolimus-eluting stent versus durable polymer everolimus-eluting stent for percutaneous coronary revascularization (BIOSCIENCE), a randomised, single-blind, non-inferiority trial. The Lancet 2016; 10.1016/S0140-6736(16)0033-2; 4. FLF as a composite of cardiac death, target vessel myocardial infarction, and clinically indicated target lesion revascularization; 5. Jensen et al. Randomized comparison of a sirolimus-eluting Orsiro stent with a bioimus-eluting Nobori stent in patients treated with percutaneous coronary intervention: Rationale and study design of the Scandinavian Organization of Randomized Trials with Clinical Outcome VII trial. 10.1016/j.amj.2015.05.009; 6. Target Lesion Failure as a composite of cardiac death, myocardial infarction (not related to other than index lesion), or target lesion revascularization; 7. Piccolo R. Biodegradable polymer sirolimus-eluting stents vs. durable polymer everolimus-eluting stents in patients with STEMI: Two-year follow-up of the BIOSCIENCE oral presentation, EuroPCR 2016; 8. Definition of probable stent thrombosis per AHA definition; 9. Preliminary analysis based on non-included data - Ion Sirolimus, poster presentation, presented at ICT, November 2016; 10. Kardžan et al. Ultrathin Bioresorbable Polymer Sirolimus-Eluting Stents versus thin durable Polymer Everolimus-eluting stents in patients Undergoing Coronary Revascularization (BIOFLOW-V), a randomized trial. The Lancet 2017; 11. Adapted from SCAAR data (January 11, 2017) http://www.scaar.se/sites/default/files/scaar_forskning_scaar; 12. Compared to other DES included in SCAAR at five years; 13. Compared to Resolute Onyx; 14. The stent system is advanced through a model, to the point of blockage (simulating a total occlusion). The force at the proximal hub and the blockage is measured. Pushability is the force transmitted along the length of the catheter; 15. (P)3/2015 - (RBP)3/2014-2; 15. (P)3/2015 - (RBP)3/2014-2; 16. Compared to Xience Expedition; 16. The stent system is advanced through a stenosis model. Crossability is the mean resistance (mean force) registered by the stenosis during the complete passage of the stent delivery system; 16. (P)3/2015 - (RBP)3/2014-2; 17. Stefanini GG, Teiwazy M, Winducker S. Coronary stents: novel development. Heart doi:10.1136/heartjnl-2012-200522; 18. Foin et al. Impact of stent strut design in metallic stents and biodegradable scaffolds. Int J Cardiol 2014 Dec 20;179(3):309-2; 19. Compared to Xience Expedition; 21. Expanded 3.0mm diameter stents are radially compressed 15% at loading full length. The force required to compress the stent is radial strength. BIOTRONIK data on file.

Sinergy and Fraxard are registered trademarks of Boston Scientific/Biosolve. Integrity, Resolute Integrity and Resolute Onyx are registered trademarks of Medtronic/Aveia. Xience Prime and Xience Expedition are registered trademarks of Abbott Cardiovascular Systems. Nobori and Ultrastent are registered trademarks of Terumo. BioMark is a registered trademark of Biosolve.

BIOTRONIK AG
Aulendorfstrasse 1
8170 Bülach, Switzerland
Tel: +41 (0) 44 8645111
Fax: +41 (0) 44 8645105
info@biotronic.com
www.biotronic.com

© 2016 BIOTRONIK AG - All rights reserved.
No part of this document may be reproduced
without the prior written consent of Biotronic AG.

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