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
myMerlinPulse™ Mobile App Compatible

Product Highlights

- Bluetooth® Low Energy (LE) communication enabling smartphone connectivity through data encryption
- DeFT Response™ technology offers non-invasive programming options to optimize rescue therapy to each patient's unique physiology and changing conditions
- VF Therapy Assurance decreases time to treatment for arrhythmias in patients who are likely to be hemodynamically unstable
- Antitachycardia pacing (ATP) while charging and prior to charging in the VF zone further extends the programming options for terminating tachyarrhythmias without a high voltage shock
- ShockGuard™ technology with DecisionTx™ programming designed to reduce inappropriate therapy and minimize the need for programming adjustments at implant
 - SecureSense™ RV lead noise discrimination algorithm detects sustained lead noise and records short bursts of oversensing that would otherwise go unnoticed or potentially lead to one or more inappropriate shocks
 - Far Field MD™ morphology discrimination and chamber onset discrimination enhance SVT and VT discrimination for reduced inappropriate therapies
- SenseAbility™ sensing algorithm feature provides the flexibility to fine-tune programming around T-wave oversensing without decreasing sensitivity
- DynamicTx™ over-current detection algorithm automatically changes shock configurations to ensure delivery of high voltage therapy when high current is detected
- 2.2.3 • MRI Ready device tested in combination with MR Conditional leads for full-body scans using a 1.5T or 3T (Tesla) field strength MRI Scanner*
- 1.10 • New battery provides higher capacity than previous QHR® batteries to offer superior longevity/volume ratio
- DF4 connector designed to streamline defibrillation connections into a single terminal pin and reduce the number of set screws
- 1.8 • Cold can programmability provides an additional RV-SVC shock configuration to decouple the can from the shocking vector parameters
- Premature Atrial Contraction (PAC) Response to avoid pacing the atrium in a vulnerable zone
- Physiologic rate responsive AV Delay and PVARP
- Dual patient notification: audio notification through the device and visual notification via myMerlinPulse™ mobile app

Ordering Information

Contents: Cardiac Pulse Generator

MODEL NUMBER	DIMENSIONS (H - W - T, MM)	WEIGHT (G)	VOLUME (CC)	CONNECTOR DEFIBRILLATION	CONNECTOR SENSE/PACE
CDDRA300Q	69 × 51 × 12	71	31	DF4	IS-I; DF4

*See MRI Scan Parameters in MRI Ready user manual

Handwritten signature

Entrant™ Dual-Chamber ICD

IMPLANTABLE CARDIOVERTER DEFIBRILLATOR (ICD) DEVICE

PHYSICAL SPECIFICATIONS

Models	CDDRA300Q
1.1.1 Telemetry	Bluetooth® LE
1.9 Delivered/Stored Energy	36/39 J
Volume	31 cc
Weight	71 g
Size	69 × 51 × 12 mm
Defibrillation Lead Connections	DF4
Sense/Pace Lead Connections	IS-1; DF4
High Voltage Can	Electrically active titanium can
PARAMETER SETTINGS	
Sensing/Detection	
SenseAbility™ Sensing Algorithm Technology	Automatic Sensitivity Control of adjustment for atrial and ventricular events
Low Frequency Attenuation	On; Off
Threshold Start	Post-Sense, Atrial/Ventricular: 50; 62.5; 75; 100% Post-Pace, atrial 0.2 - 3.0 mV Post-paced, Ventricular: Auto; 0.2 - 3.0 mV Post-Sense: 0-220 ms Post-Pace: Auto; 0-220 ms
Decay Delay	125; 157 ms
Ventricular Sense Refractory	
1.7 Detection Zones	3 zone programming — 1 zone, 2 zones or 3 zones (VT-1, VT-2, VF)
1.2 SVT Discriminators	AV Rate Branch; Arrhythmia Onset (Chamber Onset or Sudden Onset); Interval Stability; AV Association; Morphology Discrimination (Far Field MD™ Morphology Discrimination) with Automatic Template Update
Monitor Mode	Detection, discrimination and diagnostics, no therapy delivery (VT or VT-1 zone)
Discrimination Modes	On; Passive; Off
SVT Threshold	150-240 min ⁻¹
SVT Timeout	0.25-8 min
Reconfirmation	Continuous sensing during charging
SecureSense™ RV Lead Noise Discrimination Algorithm	On; On with Timeout; Passive; Off
1.4, Antitachycardia Pacing Therapy	
1.7 ATP Configurations	Ramp; Burst; Scan; 1 or 2 schemes per VT zone
ATP in VF Zone	ATP While Charging; ATP Prior to Charging; Off
ATP Upper Rate Cutoff	150-300 min ⁻¹
Burst Cycle Length	Adaptive; Readaptive or Fixed
Min. Burst Cycle Length	150-400 in increments of 5 ms
Number of Bursts	1-15
Number of Stimuli	2-20
Add Stimuli per Burst	On; Off
ATP Pulse Amplitude	2.5V independent from Bradycardia and Post-Therapy Pacing
ATP Pulse Width	1.0 or 1.5 ms independently programmable from Bradycardia and Post-Therapy Pacing
High Voltage Therapy	
VF Therapy Assurance	On; Off
DynamicTV™ Over-Current Detection Algorithm	On; Off
DeFT Response™ Technology	Programmable pulse width for PL/P2 and tilt
High-Voltage Output Mode	Fixed Pulse Width; Fixed Tilt
Waveform	Biphasic; Monophasic
RV Polarity	Cathode (-); Anode (+)
1.8 Electrode Configuration	RV to Can; RV to SVC/Can; RV to SVC
Bradycardia Pacing	
Permanent Modes	Off; DDD(R); DDI(K); VVI(R); AAI(R)
Temporary Modes	Off; DDD; DDI; VVI; AAI; AAT; DDO; VOO; AOO
1.1, Rate-Adaptive Sensor	On; Off; Passive
1.1.3 Programmable Rate and Delay Parameters	Base Rate (min ⁻¹); Rest Rate (min ⁻¹); Maximum Tracking Rate (min ⁻¹); Off; Ataximum Sensor Rate (min ⁻¹); Paced AV Delay (ms); Sensed AV Delay (ms); Rate Responsive AV Delay (Atrial and RV) (ms); Hysteresis Rate (min ⁻¹); Rate Hysteresis with Search 0.25; 0.50; 4.0; 6.5; 7.5 V
Pulse Amplitude	0.05; 0.1; 0.2; 1.5 ms
Pulse Width	On; Off
Ventricular AutoCapture™ Pacing System	On; Monitor; Off
ACap™ Confirm Feature	Sensed/Paced AV delay
QuickOp™ Timing Cycle Optimization	Off; DDI(R); VVI(R)
Auto Mode Switch (AMS)	
Atrial Tachycardia Detection Rate	110-300 min ⁻¹
AMS Base Rate	40; 45; 155 min ⁻¹
Rate-Adaptive PVA/R/VREF	Off; Low, Medium, High (ms)
PAC Response	On; Off

PAC Response Interval	200-400 in increments of 10 ms
Auto PMT Detection/Termination	Atrial Pace on PMT; Off; Passive
Ventricular Intrinsic Preference (VIP™)	Off; On (50-200 ms)

1.1.4 Post-Therapy Pacing (Independently Programmable from Bradycardia and ATP)

Post-Shock Pacing Mode	Off; AAI; VVI; DDI; DDD
Post-Shock Base Rate	30-100 in increments of 5 min ⁻¹
Post-Shock Pacing Duration	Off; 0.5; 1; 2.5; 5; 7.5; or 10 min ⁻¹

Device Testing/Induction Methods

DC Fibrey™ Pulse Duration	0.5-5.0 sec
Burst Fiber Cycle Length	30-100 ms
Non-invasive Programmed Stimulation (NIPS)	2-25 stimuli with up to three extrastimuli

1.1.2 Patient Notifiers

Programmable Notifiers (On; Off)	Battery Assurance™ alert, possible HV circuit damage, HV charge timeout, long charge time for Capacitor Maintenance, device at ERL, ventricular pacing lead impedance out of range, high voltage lead impedance out of range, AT/AF Burden, SecureSense™ lead noise detection, non-sustained ventricular oversensing, ventricular pacing percentage greater than limit, ventricular capture threshold exceeded limit, and ventricular sense amplitude lower than limit
Device Parameter Reset	On
Entry into Backup VVI Mode	On
Auditory Duration	2; 4; 6; 8; 10; 12; 14; 16 sec
Number of Audio Alerts per Notification	2
Number of Notifications	1-16
Time Between Notifications	10, 22 hours

1.5 Electrograms and Diagnostics

Stored Electrograms	Up to 15 minutes (1 user programmable + discrimination channel), up to one minute programmable pre-trigger data per VT/VF electrograms; additional triggers include lead noise detection, non-sustained ventricular oversensing, morphology template updates, atrial episode, PMT termination, magnet reversion, noise reversion
Therapy Summary	Directory listing of up to 60 episodes with access to more details including stored electrograms
Episodes Summary	History of bradycardia events and device-initiated charging
Lifetime Diagnostics	Trend data and counts
AT/AF Burden Trend	Multi-Vector Trend Data
Ventricular HV Lead Impedance Trend Histograms	Event Histogram; AV Interval Histogram; Mode Switch Duration Histogram; Peak Filtered Rate Histogram; Atrial Heart Rate Histogram; Ventricular Heart Rate Histogram; AT/AF Burden; Exercise and Activity Trending; V Rates during AMS; DirectTrend™ reports up to 1 year
PMT Data	Information regarding PMT detections
Real-Time Measurements (RTM)	Pacing lead impedances; high voltage lead impedances; and signal amplitudes

MRI Settings	
Tachy Therapy	Disabled
MRI Mode	DOO; VOO; AOO; Pacing Off
MRI Base Rate	30-100 min ⁻¹
MRI Paced AV Delay	25-120 ms
MRI Pulse Amplitude	5.0 or 7.5 V
MRI Pulse Width	1.0 ms
MRI Pulse Configuration	Bipolar
MRI timeout	Off; 3; 6; 9; 12; 24 hours

2.2.3 MRI Scan Parameters*

LEAD MODEL	LEAD LENGTHS	RF TRANSMIT CONDITIONS	SCAN REGION
Durata™ Defibrillation Lead			
7130Q; 7122Q	58, 65 cm		
Optisure™ Lead			
LDA230Q; LDA210Q	58, 65 cm	Normal Operating Mode	Full-body
Tendril MRI™ Lead†			
LPAL260M	46, 52 cm		
Tendril™ STS Pacing Lead			
2088TC	46, 52 cm		

* For additional information about MRI conditions for ICDs and leads, including scan parameters, warnings, precautions, adverse conditions to MRI scanning, and potential adverse events, please refer to the Abbott MRI Ready Systems Medical manuals, www.abbott.com/mri, or our MRI Ready resources at mri.meritlin.net.


† Tendril MRI™ lead is AIP Conditional for LST MRI only.

Customer Support: 46-8-474-475e

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

Abbott
The Corporate Village
Da Vinci Lane II Box FJ
1935 Zaventem, Belgium
Tel: +32 2 774 68 11
Abbott.com

™ Indicates a trademark of the Abbott group of companies.
† Indicates a third-party trademark, which is property of its respective owner.
Bluetooth and the Bluetooth logo are registered trademarks of Bluetooth SIG, Inc.
© 2019 Abbott. All Rights Reserved.
MAT1900552 - Item approved for ID use only.


Ellipse™ DR

Dual-Chamber Implantable Cardioverter Defibrillator (ICD)

Implantable Cardioverter Defibrillator (ICD) Devices

Product Highlights

- MRI Ready device has been tested for safe performance of an MRI scan using a 1.5 T (Tesla) field-strength MRI scanner when used in combination with MRI Conditional leads^{1,2}
- Improved shape with reduced volume and thickness
- Parylene coating for improved abrasion resistance
- DynamicTx™ Over-Current Detection Algorithm automatically changes shock configurations to ensure delivery of high voltage therapy when high current is detected
- Cold Can programmability provides an additional RV-SVC Shock Configuration to decouple the can from the shocking vector parameters in cases of lead problems
- ShockGuard™ technology with DecisionTx™ programming designed to reduce inappropriate therapy and minimise the need for programming adjustments at implant
 - SecureSense™ RV lead noise discrimination detects sustained and short bursts of lead noise that would otherwise go unnoticed or potentially lead to one or more inappropriate shocks
 - Far Field MD™ morphology discrimination and Chamber Onset discrimination improve SVT and VT discrimination for reduced inappropriate therapies
- Low Frequency Attenuation filter designed to enhance sensing performance and may reduce the possibility of oversensing T-waves
- SenseAbility™ feature provides flexibility to fine-tune programming around T-wave oversensing without decreasing sensitivity
- DF4 connector designed to streamline defibrillation connections into a single terminal pin and reduce the number of set screws
- CorVue™ congestion monitoring feature monitors the intrathoracic impedance in multiple vectors for improved accuracy, and it provides the option for both patient and physician alerts
- Antitachycardia pacing (ATP) while charging and prior to charging in the VF zone further extends the programming options for terminating tachyarrhythmias without a high-voltage shock
- ST monitoring capability provides unprecedented, continuous insight into significant ST shift events and associated ventricular arrhythmias through enhanced monitoring of iEGM and ST-segment as a diagnostic tool to help guide appropriate clinical action
- 36 J delivered energy safety shock option can provide a greater DFT safety margin
- DeFT Response™ technology offers the most noninvasive options for managing high DFTs
- QHR™ chemistry battery provides greater capacity for enhanced longevity and improved charge time performance compared to previous SVO batteries

Merlin@home™
Transmitter
Compatible

Ordering Information

Contents: Dual-Chamber Implantable Cardioverter Defibrillator (ICD)

1, 6

Model Number	Dimensions (H x W x D mm)	Weight (g)	Volume (cc)	Connector Defibrillation	Connector Sense/Pace
CD2377-36C	69 x 51 x 12	65	31	DF1	IS-1
CD2377-36QC*	70 x 51 x 12	68	31	DF4	IS-1; DF4

*Indicates models that are MRI Conditional^{1,2}

Indications: The devices are intended to provide ventricular antitachycardia pacing and ventricular defibrillation for automated treatment of life-threatening ventricular arrhythmias.

Contraindications: Contraindications for use of the implantable cardioverter defibrillator include ventricular tachyarrhythmias resulting from transient or correctable factors such as drug toxicity, electrolyte imbalance or acute myocardial infarction.

Adverse Events: Implantation of the implantable cardioverter defibrillator, like that of any other device, involves risks, some possibly life-threatening. These include but are not limited to the following: acute hemorrhage, bleeding, or emboli, arrhythmia acceleration, cardiac or venous perforation, cardiogenic shock, cyst formation, erosion, exacerbation of heart failure, extrusion, fibrotic tissue growth, fluid accumulation, hematoma formation, histotoxic reactions, infection, keloid formation, myocardial irritability, nerve damage, pneumothorax, thrombemboli, vessel occlusion. Other possible adverse effects include mortality due to component failure,

device-programmer communication failure, lead abrasion, lead dislodgment or poor lead placement, lead fracture, inability to defibrillate, inhibited therapy for a ventricular tachyarrhythmia, interruption of therapy, electrical or magnetic interference, bleeding of energy from defibrillation paddles, system failure due to X-ray radiation. Other possible adverse effects include mortality due to inappropriate delivery of therapy caused by multiple counting of cardiac events including T-waves, P waves, or supraventricular premature signals. Adverse psychological effects of device implantation are imagined pulsing, dependency, fear of inappropriate pulsing and fear of losing pulse capability.

Refer to the User's Manual for detailed indications, contraindications, warnings, precautions and potential adverse events.

† QHR is a trademark of Greatbatch Medical

Ellipse™ DR

Lead-Driven or Any Implantable Cardioverter Defibrillator (ICD)

Implantable Cardioverter Defibrillator (ICD) Devices

Product Specifications

Physical Specifications

Models	CD2377-360	CD2377-360C
Topology	RF	RF
Delivered/Stored Energy (J)	35/79	35/79
Volume (cc)	31	31
Weight (g)	66	68
Size (mm)	59 x 51 x 12	70 x 51 x 12
Defibrillation Lead Connections	DF1	DF4
Sense/Pace Lead Connections	IS-1	IS-1, DF4
High-Voltage Can Coating	Electrically active titanium can	Electrically active titanium can
MRI Conditional	Parylene	Parylene Yes-MRI ready

Parameter Settings

Parameter	Settings
AF Management	
AF Suppression™ Pacing	On; Off
No. of Overdrive Pacing Cycles	15-40 in steps of 5
Maximum AF Suppression Rate	15-40 in steps of 5
Sensing/Detection	
SenseAbility™ Technology	Automatic Sensitivity Control adjustment for atrial and ventricular events
Low-Frequency Attenuation Threshold Start	On; Off (Post-Sensed, Atrial) 50; 62.5; 75; 100%; (Post-Paced, Atrial) 0.2-3.0 mV; (Post-Sensed, Ventricular) 50; 62.5; 75; 100%; (Post-Paced, Ventricular) Auto: 0.2-3.0 mV
Decay Delay	(Post-Sense/Post-Pace; Atrial/Ventricular) 0-225
Ventricular Sense Refractory (ms)	125; 157
Date, Item Zones	3 zone programming - 1 zone, 2 zones or 3 zones (V1-1, V1-2, V1)
SVT Discriminators	AV Rate Branch: Arrhythmia Onset (Chamber Onset or Sudden Onset); Interval Stability; AV Association; Morphology Discrimination (Far Field MD or Original MD) with Manual (Original MD only) or Automatic Template Update
Metric Mode	Detection, discrimination and diagnostics, no therapy delivery (VT or VT-1 zone)
Discrimination Modes	
SVT Threshold	On; Passive; Off
SVT Timeout	150-240 min
Recharge/Charge	0.25-5 min
Lead Noise Discrimination	Continuous sensing during charging; SecureSense™ RV lead noise discrimination (On, Off with Timeout; Passive; Off)
Anti-tachycardia Pacing Therapy	
ATP Configurations	Ramp, Burst, Scan, 1 or 2 schemes per VT zone
ATP in VF Zone	ATP While Charging; ATP Prior to Charging; Off
ATP Upper Rate Cutoff	150 - 300/min
Burst Cycle Length	Adaptive, Readaptive or Fixed
Min. Burst Cycle Length (ms)	150-400 in increments of 5
Number of Bursts	1-15
Number of Stimuli	2-20
Add Stimuli per Burst	On; Off
ATP Pulse Amplitude (V)	7.5 Independent from Bradycardia and Post-Therapy Pacing
ATP Pulse Width (ms)	1.0 to 1.5 Independent from Bradycardia and Post-Therapy Pacing
High-Voltage Therapy	
Dynamic Tx™ Algorithm	On; Off
DefT Responsiveness™ Technology	Programmable pulse width for P1/P2 and full
High-Voltage Output Mode	Fixed Pulse Width, Fixed TH
Waveform	Biphasic, Monophasic
RV Polarity	Cathode (-), Anode (+)
Electrode Configuration	RV to Can; RV to SVC/Can; RV to SVC
Bradycardia Pacing	
Permanent Modes	Off; DDD(R); DDI(R); VVIR(R); AA(R)
Temporary Modes	Off; DDD; DDI; VVI; AA; AAT; DDD; VDD; AOO
Rate-Adaptive Sensor	On; Off; Passive
Programmable Rate and Delay Parameters	Base Rate (min ⁻¹); Rest Rate (min ⁻¹); Maximum Tracking Rate (min ⁻¹); Off; Maximum Sensor Rate (min ⁻¹); Paced AV Delay (ms); Sensed AV Delay (ms); Rate Response AV Delay (Atrial and RV) (ms); Hysteresis Rate (min ⁻¹); Rate Hysteresis with Search
Ventricular AutoAdapt™ Pacing System	
Pacing System	On; Off
ACap™ Confirm	On; Monitor; Off
QuasiOpt™ Timing Cycle Optimisation	Sensed/Paced AV delay
Auto Mode Switch (AMS)	Off; On(R); VVIR(R)
Atrial Tachycardia Detection Rate (min ⁻¹)	110-300
AMS Base Rate (min ⁻¹)	20; 45; ...; 135
Auto PWT Detection/Termination	Atrial Pace on PWT; Off; Passive
Rate Response PVARP/VRET	Off; Low; Medium; High
Ventricular Intense Protection (VIP™)	Off; On (50-200)

Post-Therapy Pacing (Independently programmable from Bradycardia and ATP)

Post-Shock Pacing Mode	Off; AAI; VVI; DDI; DDD
Post-Shock Base Rate (min ⁻¹)	35-100 in increments of 5
Post-Shock Pacing Duration (min)	Off; 0.5, 1, 2, 5, 7.5 or 10

Device Testing/Induction Methods

DC Fiber™ Pulse Duration (sec)	0.5-5.0
Burst Fiber Cycle Length (ms)	70-100
Noninvasive Programmed Stimulation (NIPS)	2-25 stimuli with up to three extrastimuli

Patient Notifiers

Programmable Notifiers (On, Off)	Device or ERI; Charge Time Limit Reached; Possible HV Circuit Damage; Atrial Lead Impedance Out of Range; Ventricular Lead Impedance Out of Range; High-Voltage Lead Impedance Out of Range; ATAF Burden; V Rate During ATAF; % V pacing; CorVue Congestion Trigger; SecureSense™ lead noise detected; non-sustained lead noise detected; ST Episodes (Type 1) only
Device Parameter Reset	On
Entry into Backup VVI Mode	On
Vibration Duration (sec)	2; 4; 6; 8; 10; 12; 14; 16
Number of Vibrations per Notification	2
Number of Notifications	1-16
Time Between Notifications (hours)	10; 22

Electrograms and Diagnostics

Stored Electrograms	Up to 45 minutes including up to one minute programmable pre-trigger data per V/VF diagnosis/detection electrograms; triggers include diagnosis; detection; therapy; atrial episode; PMT termination; FC shock delivery; pause/resume; magnet event/sensing; morphology template verification; lead noise detected; non-sustained lead noise detected; NSVT/NSVF
Therapy Summary	Diagram of therapies delivered
Episodes Summary	Diagnostics listing of up to 50 episodes with access to more details including stored electrograms
Lifetime Diagnostics	History of bradycardia events and device-initiated charging
ATAF Burden Trend	Trend data and counts
Ventricular HV Lead Impedance Trend	Multi-Vector Trend Data
Histograms	Event Histogram; AV Interval Histogram; Mode Switch Duration Histogram; Peak Filtered Rate Histogram; Atrial Fibrillation Rate Histogram; Ventricular Fibrillation Rate Histogram; ATAF Burden; Exercise and Activity Trending; V Rates during AMS; Over-Trend™ reports up to 1 year
PMT Data	Information regarding PMT detections
Real-Time Measurements (RTM)	Pacing lead impedances; high-voltage lead impedances; and signal amplitudes
ST Monitoring	ST Histogram Data; Long-term ST Deviation Trend; ST Episode Log; ST Episode Details; 24-Hour ST & HR Trend; ST EGM Search & Signals prior to ST Episode; V/VF Inducement (Snapshots and 24-Hour trend at time of interrogation)
CorVue™ Congestion Monitoring	On; Off
CorVue Congestion Trigger	5-18 days

MRI Scan Restrictions

If the implanted system is comprised of a combination of leads that have differing whole body SAR, scan region and/or additional considerations, use the most restrictive of each to determine that overall set of scan conditions applicable for the total system.

Lead Model	Whole Body SAR	Scan Zone Restrictions
Durata™ Lead		
7120G (lead length: 58 cm; 76 cm)	≤ 2 W/kg	If MRI Mode = Pacing OFF: No scan zone restrictions
7122G (lead length: 1.5 m)	≤ 2 W/kg	
7122Q (lead length: 55 cm)	≤ 1.0 W/kg	
Optisure™ Lead		
LD42201 (lead length: 58 cm; 65 cm)	≤ 2 W/kg	If MRI Mode = VOO or DOO: Isocenter cannot be placed between C1 and L2
LD4210R (lead length: 58 cm)	≤ 2 W/kg	
LD4210G (lead length: 55 cm)	≤ 1.6 W/kg	
Tendril™ MRI Lead		
UP41206M (lead length: 16; 62; 58 cm)	≤ 2 W/kg	


- MRI Conditional Field Strength: 1.5 Tesla
- See MRI Procedure Information for approved MR Conditional Systems Device/Lead combinations and scan parameters.

Customer Support: +6-8-474-4766

Brief Summary: Prior to using these devices, please review the instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use. Devices depicted may not be available in all countries. Check with your St. Jude Medical representative for product availability in your country.

Unless otherwise noted, ™ indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST, JDE, MPTICAL and the diamond symbol are trademarks and service marks of St. Jude Medical, Inc. and its related companies. © 2015 St. Jude Medical, Inc. All Rights Reserved.

5AM-ELF-1114-0002(1) Item approved for international use only.

Tools ▾

1.3 Ventricle Output

MRI Device

V. Pulse Amplitude (V)

7.5

0.25

3.0
2.75
2.5 N
2.25
2.0
1.75

⌫ ⏪ ⏩ ⌛

Emergency Capture Settings

Backup Pulse Configuration, Sketch Protocol

Pulse Amplitude: 2.5 V

Pulse Width: 0.5 ms

Charge: ON

Charge Delay: OFF

Max Sensitivity: Auto

2020

10/10/2020

10:10:10

10/10/2020

10:10:10

10/10/2020

10:10:10

?

Tools ▾

1.3 Ventricle Pulse Width

MRI Device

Configure 1.3 Settings

- Monitor Off
- Pulse Amplitude 2.5 V 2.5 V
- Pulse Width 0.5 ms 0.5 ms
- Synchrony On Off
- Sensitivity Auto Auto

Secondary Capture Settings
 Refresh Pulse Configuration Search Interval

Searchability: Use Filter
 and Secondary Additional Settings

V. Pulse Width (ms)

1.5

0.7

0.6

0.5 N

0.4

0.3

0.2

0.05

X

◀

▶

[Handwritten signature]

Tools

MRI Device

<30 min⁻¹

Ready

1.7

ShockGuard™ Settings (Zone Configuration)

Zone Configuration

1 Zone 2 Zones **3 Zones** Off

Detection Criteria

VT-1	VT-2	VF
<p>150 <small>min⁻¹</small> / 400 ms</p> <p>18 intervals</p> <p>On</p>	<p>181 <small>min⁻¹</small> / 350 ms</p> <p>16 intervals</p> <p>On</p>	<p>214 <small>min⁻¹</small> / 280 ms</p> <p>12 intervals</p>

SVT Discrimination

Chamber Timeout **Dual** / **Off**

SVT Details...

ATP, VT Timeout Settings...

<p>ATP x 3</p> <p>ATP x 3</p> <p>36.0 J</p> <p>36.0 J x 2</p>	<p>ATP x 3</p> <p>ATP x 3</p> <p>36.0 J</p> <p>36.0 J x 2</p>	<p>ATP x 1</p> <p>36.0 J</p> <p>36.0 J x 4</p>
---	---	--

Secure Sense™ V Lead Noise Discrimination

SecureSense™ Configuration **On**

Additional Settings... **RVC-Cam**

Preview 154 Program

[Handwritten signature]

Table 20. Lead compatibility

Device	Lead compatibility
	LV low voltage: one IS4-LLLL lead connector (LV) lead

WARNING: Do not use another manufacturer's lead system without demonstrated compatibility as undersensing cardiac activity, and failure to deliver necessary therapy may result. **1.10**

Battery Information

Table 21. Battery information

Device	
Battery chemistry;	xxxxx300Q, xxxxx500Q, xxxxx600Q, xxxxx700Q, xxxxx800Q
Manufacturer; Model; Cells	Silver vanadium oxide/carbon monofluoride; Greatbatch Medical®; Model 3451; One cell QHR

DEFIBRILLATION LEADS

Durata™
Defibrillation Lead

Product Highlights

- The DF4 connector is designed to simplify implants by streamlining defibrillation connections into a single terminal pin and reducing the number of set screws
- Redundant conductors serve as a backup system in the unlikely event of a conductor failure
- Symmetrically aligned cables within the lead body and centrally located coil provide for additional protection to the inner coil
- Optim™ lead insulation – a chemical co-polymer that blends the best features of polyurethane and silicone for improved handling and increased durability
- Two innovative designs are intended to help prevent tissue ingrowth – flat-wire technology provides a low profile for the defibrillation coils, and silicone backfilling completely fills the shock coil space

Ordering Information

Contents: Defibrillation Lead

REORDER NUMBER	INSULATION	FIXATION	MINIMUM INTRODUCER (F)	SHOCK CONFIGURATION	SENSING	TIP-TO-PROXIMAL COIL (CM)	CONNECTOR ⁵	LENGTHS (CM)
7120	Optim	Ext/Ret Helix	8	Dual-coil	True bipolar	17	DF1; IS-1	60; 65; 75
7120Q	Optim	Ext/Ret Helix	7	Dual-coil	True bipolar	17	DF4	52; 58; 65
7121	Optim	Ext/Ret Helix	7	Dual-coil	True bipolar	21	DF1; IS-1	60; 65; 75
7121Q	Optim	Ext/Ret Helix	7	Dual-coil	True bipolar	21	DF4	52; 58; 65
7122	Optim	Ext/Ret Helix	7	Single-coil	True bipolar	N/A	DF1; IS-1	60; 65; 75
7122Q	Optim	Ext/Ret Helix	7	Single-coil	True bipolar	N/A	DF4	52; 58; 65
7170	Optim	Tines	7	Dual-coil	True bipolar	17	DF1; IS-1	60; 65; 75
7170Q	Optim	Tines	7	Dual-coil	True bipolar	17	DF4	52; 58; 65
7171	Optim	Tines	7	Dual-coil	True bipolar	21	DF1; IS-1	60; 65; 75
7171Q	Optim	Tines	7	Dual-coil	True bipolar	21	DF4	52; 58; 65
7172Q	Optim	Tines	7	Single-coil	True bipolar	N/A	DF4	52; 58; 65

Refer to the User's Manual for detailed indications, contraindications, warnings, precautions and potential adverse events.


PRODUCT SPECIFICATIONS

True Bipolar, Active-fixation Defibrillation Leads

Models	7120	7120Q	7121	7121Q	7122	7122Q
Fixation	Ext/Ret Helix	Ext/Ret Helix	Ext/Ret Helix	Ext/Ret Helix	Ext/Ret Helix	Ext/Ret Helix
Shock Configuration	Dual-coil	Dual-coil	Dual-coil	Dual-coil	Single-coil	Single-coil
Sensing Configuration	True Bipolar	True Bipolar	True Bipolar	True Bipolar	True Bipolar	True Bipolar
Min. Size Introducer	7 F	7 F	7 F	7 F	7 F	7 F
Lengths (cm)	60; 65; 75	52; 58; 65	60; 65; 75	52; 58; 65	60; 65; 75	52; 58; 65
Connector	DF1; IS-1	DF4	DF1; IS-1	DF4	DF1; IS-1	DF4
Body Diameter	6,8 F	6,8 F	6,8 F	6,8 F	6,8 F	6,8 F
Tip-to-Anode Spacing	11 mm	11 mm	11 mm	11 mm	11 mm	11 mm
Tip-to-Proximal Coil	17 cm	17 cm	21 cm	21 cm	N/A	N/A
Tip Electrode Area	6 mm ²	6 mm ²	6 mm ²	6 mm ²	6 mm ²	6 mm ²
Steroid Plug	Yes	Yes	Yes	Yes	Yes	Yes
Distal Shock Coil Area	367 mm ²	367 mm ²	367 mm ²	367 mm ²	367 mm ²	367 mm ²
Proximal Shock Coil Area	588 mm ²	588 mm ²	588 mm ²	588 mm ²	N/A	N/A

True Bipolar, Passive-fixation Defibrillation Leads

Models	7170	7170Q	7171	7170Q	7172Q
Fixation	Tines	Tines	Tines	Tines	Tines
Shock Configuration	Dual-coil	Dual-coil	Dual-coil	Dual-coil	Single-coil
Sensing Configuration	True Bipolar	True Bipolar	True Bipolar	True Bipolar	True Bipolar
Min. Size Introducer	7 F	7 F	7 F	7 F	7 F
Lengths (cm)	60; 65; 75	52; 58; 65	60; 65; 75	52; 58; 65	52; 58; 65
Connector	DF1; IS-1	DF4	DF1; IS-1	DF4	DF4
Body Diameter	6,8 F	6,8 F	6,8 F	6,8 F	6,8 F
Tip-to-Anode Spacing	11 mm	11 mm	11 mm	11 mm	11 mm
Tip-to-Proximal Coil	17 cm	17 cm	21 cm	21 cm	N/A
Tip Electrode Area	3.5 mm ²	3.5 mm ²	3.5 mm ²	3.5 mm ²	3.5 mm ²
Steroid Plug	Yes	Yes	Yes	Yes	Yes
Distal Shock Coil Area	367 mm ²	367 mm ²	367 mm ²	367 mm ²	367 mm ²
Proximal Shock Coil Area	588 mm ²	588 mm ²	588 mm ²	588 mm ²	N/A

Customer Support: 46-8-474-4756

™ Indicates a trademark of the Abbott group of companies.

† Indicates a third party trademark, which is property of its respective owner.

© 2020 Abbott. All Rights Reserved

MAT-2001642 | Item approved for Hong Kong, Taiwan, Thailand, Vietnam and Philippines use only.

Indications for Use ^{2.1, 2.2.1}

The Durata™ Models 7120Q, 7121Q, and 7122Q transvenous leads are indicated for use with compatible pulse generators (refer to the applicable defibrillator manual for system indications). They provide pacing and sensing and deliver cardioversion/defibrillation therapy to the heart.

A transvenous lead system may offer the patient the benefit of avoiding a thoracotomy for lead implantation. If the initial lead configuration is not effective, repositioning of the lead or other lead configurations should be attempted.

In some patients, a nonthoracotomy lead configuration may not provide reliable conversion of arrhythmias, and the use of subcutaneous or epicardial patch defibrillation leads should be considered.

Contraindications

Contraindications for use of the Durata™ leads with an implantable pulse generator include ventricular tachyarrhythmias resulting from transient or reversible factors such as drug toxicity, electrolyte imbalance, or acute myocardial infarction. Transvenous lead systems are contraindicated for patients with tricuspid valvular disease or a mechanical heart valve. Durata leads are contraindicated for patients for whom a single dose of 1.0 mg of dexamethasone sodium phosphate is contraindicated. Durata leads are contraindicated for use with extra firm (red color knob) stylets.

The lead is not designed, sold, or intended for use other than as indicated.

Lead Package Contents

All leads are packaged separately and supplied sterile. Package contents are as follows:

- 1 lead with stylet and suture sleeve in place
- 1 vein pick
- 2 clip-on tools
- Stylets of varying firmness
- Product documentation

Packages may also contain:

- 1 extra suture sleeve
- 1 lead cap

Potential Complications

Possible complications of the use of transvenous lead systems include, but are not limited to, supraventricular or ventricular arrhythmias, conduction disturbances, cardiac perforation, cardiac tamponade, loss of contractility, air embolism, heart wall rupture, myocarditis, post-operative heart failure, chronic mechanical stimulation of the heart, tricuspid valve dysfunction, lead fracture necessitating surgical removal, pneumothorax, hemothorax, infection, tissue necrosis, and erosion of the skin. Specific events and effects are summarized in Table 1.

A handwritten signature in black ink is written over a faint, circular stamp. The signature appears to be 'M. J. ...' and the stamp contains some illegible text.

Indications for Use

PACING LEADS

Tendril™ STS

Pacing Lead

Product Highlights - Pacing Leads

- The Tendril™ STS lead allows patients to undergo 1.5 T or 3 T MRI scans when used in conjunction with an MRI Ready device from St. Jude Medical*
- Soft silicone tip offers more compliance and less tip pressure at the lead tip-endocardium interface
- Small diameter lead offers improved ease of venous passage, reduced risk of venous thrombosis or rib-clavicle crush and ability to accommodate additional leads more easily
- Optim™ lead insulation — a chemical co-polymer that blends the best features of polyurethane and silicone for improved handling and increased durability
- Titanium nitride (TiN) fractal coating on the tip and ring electrodes is designed to promote precise sensing and to provide improved contact with the myocardium
- Lubricious Fast-Pass™ coating facilitates lead insertion through the introducer and veins to ease implantation
- Fits through a 6 F introducer

Ordering Information

MODEL NUMBER	DESCRIPTION	INSULATION	FIXATION	MIN. INTRODUCER (F)	CONNECTOR	LENGTH (CM)
2088TC	Tendril™ STS Pacing Leads	Optim™	Ext/Ret helix	6	IS-1 bipolar	46**; 52**; 58**; 65; 100

See MRI conditional parameters.

**Indicates lead lengths that are MRI conditional.

Indications: Tendril™ STS Lead is designed for permanent sensing and pacing in either the right atrium or the right ventricle, in combination with a compatible device. Active leads such as the Tendril STS lead may be indicated for patients where permanent fixation of passive leads is suspected to be unstable.

In atrial applications, the use of screw-in leads such as Tendril STS lead may be indicated in the presence of an abnormal, surgically altered or excised atrial appendage.

Contraindications: Tendril STS lead is contraindicated: in the presence of tricuspid atresia, for patients with mechanical tricuspid valves, in patients who are expected to be hypersensitive to a single dose of one milligram of dexamethasone sodium phosphate

Adverse Events: Potential complications associated with the use of Tendril STS lead are the same as with the use of other active fixation leads and include cardiac tamponade, diaphragmatic stimulation, embolism, excessive bleeding, induced ventricular ectopy, infection, loss of pacing and/or sensing due to dislodgment or mechanical malfunction of the pacing lead, phrenic nerve stimulation, thrombosis. Complications reported with direct subclavian venipuncture include pneumothorax, hemothorax, laceration of the subclavian artery, arteriovenous fistula, neural damage, thoracic duct injury, cannulation of other vessels, massive hemorrhage and, rarely, death.

Refer to the User's Manual for detailed indications, contraindications, warnings, precautions and potential adverse events.

Product Specifications

PHYSICAL SPECIFICATIONS

Models	208STC
3.5 Minimum Introducer Size	6 F
3.2 Type of Lead	Active-fixation, steroid-eluting, endocardial, straight pacing lead
3.3 Lead Connector	IS-1 bipolar
3.4 Lead Lengths	46*; 52*; 58*; 65; 100 cm
Fixation Mechanism	Extendable/Retractable helix
Typical Number of Rotations for Helix Extension	6-11 (straight stylet)
Lead Body Diameter	1.9 mm (max)
Tip-to-Ring Spacing	10 mm
Lead Tip Electrode (Cathode)	Active Titanium-nitride-coated Pt/Ir helix (2.0 mm extension)
Tip Electrode Surface Area	6.9 mm ²
Ring Electrode (Anode)	Titanium-nitride-coated Pt/Ir
Ring Electrode Surface Area	16 mm ²
Mapping	Capable with Titanium-nitride-coated Pt/Ir helix
3.1 Steroid	< 1 mg dexamethasone sodium phosphate
Inner Conductor/Outer Conductor	MP35N™ coil
Inner Insulation	Silicone rubber
Outer Insulation	Optini™ lead insulation
Lead Body Coating	Fast-Pass™ coating

In-Pack

Straight stylets	1 x-soft in lead, 1 x-soft, 1 soft
J-shaped stylets	2 soft
Helix extension/retraction clip-on tools	2 clip-on tools

Accessory Kits Available Separately	Model Number	Compatible Lengths	Description
Stylet Kit	DSO6002 with appropriate length designation	46; 52; 58; 65; 100 cm	1 fixation tool, 1 clip-on tool, 1 J-shaped soft, 1 x-soft, 1 soft, 1 firm, 1 x-firm
	DSO6003 with appropriate length designation	46; 52; 58; 65; 100 cm	1 clip-on tool, 1 J-shaped soft, 1 x-soft, 1 soft, 1 firm, 1 x-firm
Locator™ Plus Deflectable Stylet	1291 with appropriate length designation	46; 52; 58; 65 cm	Disposable implant tool to facilitate precise lead positioning
	1292 with appropriate length designation	52; 58; 65 cm	and manipulation with one hand

MRI Conditional Parameters

MRI scan parameters vary depending on the MRI Ready device used. Consult the MRI Ready System Manual for specific product combinations and associated MRI scan parameters.

MP35N is a trademark of SPS Technologies, LLC


Customer Support: 46-8-474-4756

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing indications, contraindications, warnings, precautions, potential adverse events and directions for use

™ Indicates a trademark of the Abbott group of companies.

© 2017 Abbott. All Rights Reserved

22780-SJM-TND-0914-0002(6) ; Item approved for international use only

Abbott

Limited Warranty

Entrant™ Tiered-Therapy Defibrillators

For device models:

CDDRA300Q

CDHFA300Q

CDVRA300Q

The purchaser must comply with the terms and conditions in this document for this Limited Warranty to apply.

This Limited Warranty provides credit for the replacement of the Abbott CRMD ("Abbott") Tiered-Therapy Defibrillator ("Device"), as indicated by the type of device below, to meet the written operating specifications due to defects in materials or workmanship, or due to battery depletion as specified below.

Cardiac Resynchronization Therapy-Defibrillator ("CRT-D") Devices

Subject to the terms and conditions of this Limited Warranty, Abbott will provide a full credit for replacement of the CRT-D Device:

- if the CRT-D Device fails to function within written operating specifications (when used as intended under normal surgical conditions and in conformance with the instructions for its use and maintenance) due to defects in materials or workmanship during the first 5 years (60 months) after date of implant, or
- if the CRT-D Device battery reaches its Elective Replacement Indicator (ERI) voltage (as defined in the Merlin™ PCS on-screen help) during the first 36 months after date of implant.

Subject to the terms and conditions of this Limited Warranty, Abbott will provide a partial credit for replacement of the CRT-D Device as follows: if the battery reaches its ERI voltage (as defined in the Merlin™ PCS on-screen help) during months 37-60 after date of implant, Abbott will issue a replacement credit for 50% of the net purchase price of the original CRT-D Device decreased on a sliding scale to the date of replacement over this 24-month period.

Dual Chamber (DR) Devices

Subject to the terms and conditions of this Limited Warranty, Abbott will provide a full credit for replacement of the DR Device:

- if the DR Device fails to function within written operating specifications (when used as intended under normal surgical conditions and in conformance with the instructions for its use and maintenance) due to defects in materials or workmanship during the first 6 years (72 months) after date of implant, or
- ⁴ if the DR Device battery reaches its Elective Replacement Indicator (ERI) voltage (as defined in the Merlin™ PCS on-screen help) during the first 60 months after date of implant.

Subject to the terms and conditions of this Limited Warranty, Abbott will provide a partial credit for replacement of the DR Device as follows: if the battery reaches its ERI voltage (as defined in the Merlin™ PCS on-screen help) during months 37-72 after date of implant, Abbott will issue a replacement credit for 50% of the net purchase price of the original DR Device decreased on a sliding scale to the date of replacement over this 36-month period.

Single Chamber (VR) Devices

Subject to the terms and conditions of this Limited Warranty, Abbott will provide a full credit for replacement of the VR Device:

- if the VR Device fails to function within written operating specifications (when used as intended under normal surgical conditions and in conformance with the instructions for its use and maintenance) due to defects in materials or workmanship during the first 7 years (84 months) after date of implant, or
- if the VR Device battery reaches its Elective Replacement Indicator (ERI) voltage (as defined in the Merlin™ PCS on-screen help) during the first 60 months after date of implant.

Subject to the terms and conditions of this Limited Warranty, Abbott will provide a partial credit for replacement of the VR Device as follows: if the battery reaches its ERI voltage (as defined in the Merlin™ PCS on-screen help) during months 49-84 after date of implant, Abbott will issue a replacement credit for 50% of the net purchase price of the original VR Device decreased on a sliding scale to the date of replacement over this 36-month period.

ABBOTT DOES NOT WARRANT THE SUITABILITY OF A DEVICE FOR ANY SPECIFIC PATIENT, SINCE FITNESS FOR USE IS A MEDICAL DECISION.

Other Important Terms, Conditions, and Limitations

The following additional terms, conditions, and limitations apply to any of the above warranty coverages under this limited warranty:

This warranty applies only for an Abbott device replacement in the original patient.

The credit may only be used toward the purchase of an equivalent or functionally equivalent Abbott device.

The Device must be implanted before the "USE BEFORE" date marked on the package.

The Device must be returned for analysis to Abbott Cardiac Rhythm Management, Attention: Product Analysis, 645 Almenor Avenue, Sunnyvale, CA 94085 USA or a local Abbott affiliate within forty-five (45) days after removal from the patient.

Abbott will analyze the Device and determine whether it qualifies for warranty coverage. No credit allowance will be made where there is evidence of improper handling, exposure to temperature above 60°C or below -20°C, improper storage, subjection to shock or electrical abuse, disassembly or alteration of the Device or any other usage not conforming with the instructions for its use and maintenance.

The purchaser must inspect the Device upon receipt. If the Device is received in a damaged condition, Abbott will replace it at no charge to the purchaser if the damage is reported to Abbott within fifteen (15) days of receipt of the Device and it is promptly returned to Abbott.

All Devices returned to Abbott become the property of Abbott.

In no event shall the amount of the replacement credit exceed the net invoiced price of the replacement Abbott device.

THIS WARRANTY IS IN LIEU OF, AND ABBOTT HEREBY DISCLAIMS AND EXCLUDES, ALL OTHER WARRANTIES, REPRESENTATIONS, OR CONDITIONS, WHETHER EXPRESS OR IMPLIED BY OPERATION OF LAW OR OTHERWISE, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY, NON-INFRINGEMENT, COURSE OF DEALING, QUIET ENJOYMENT, OR FITNESS FOR A PARTICULAR PURPOSE.

ABBOTT WILL NOT BE RESPONSIBLE FOR ANY INDIRECT, SPECIAL, COVER, PUNITIVE, CONSEQUENTIAL, INCIDENTAL OR OTHER DAMAGES DIRECTLY OR INDIRECTLY ARISING FROM THE USE OF OR RELATED TO THE DEVICE. ABBOTT'S SOLE LIABILITY AND THE PURCHASER'S SOLE REMEDY RELATED TO THE DEVICE IS THE CREDIT DESCRIBED IN THIS LIMITED WARRANTY.

IN NO EVENT SHALL ABBOTT BE LIABLE TO THE PURCHASER OR ANY THIRD PARTY FOR ANY CLAIM OR DAMAGES, HOWEVER ARISING, IN AN AMOUNT THAT EXCEEDS THE DEVICE PURCHASE PRICE.

SOME COUNTRIES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS OR THE EXCLUSION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THE PRECEDING LIMITATIONS OR EXCLUSIONS MAY NOT APPLY TO THE PURCHASER.

THE REMEDIES SET FORTH IN THIS WARRANTY WILL BE THE ONLY REMEDIES AVAILABLE TO ANY PERSON. NO PERSON HAS ANY AUTHORITY TO BIND ABBOTT TO ANY WARRANTY OR REPRESENTATION EXCEPT THOSE SPECIFICALLY CONTAINED HEREIN.

THE ABOVE TERMS, DISCLAIMERS AND LIMITATIONS SHALL BE CONSTRUED TO COMPLY WITH APPLICABLE LAWS AND THIS LIMITED WARRANTY SHALL BE REFORMED/REVISED ACCORDINGLY AS NECESSARY TO COMPLY WITH APPLICABLE LAWS OF THE COUNTRY WHERE THE PRODUCT IS IMPLANTED.

Handwritten signature and stamp

Handwritten signature

201404
AR1560010010 A

Entrant™ DR

Implantable Cardioverter Defibrillator

عمرس الرخمان وعموم نظم القلب العرس الررس (bg)Implantirpyer kardioverter-defibrilator
 (ca)Implantovateľný kardioverter defibrilátor (da)Implanterbar cardioverter-defibrillator (de)Implantierbare
 Kardioverter-Defibrillator (el)Επιχειρήσιμος καρδιοπαταριστής αυτοβύθισης (es)Defibrilador cardioversor
 implantable (et)Implanteeritav kardioverter-defibrilator (fi)Implantoitava sydäntaheräin (fr)Défibrillateur
 automatique implantable (gv)Imprantibla kardioverter-defibrilator (hu)Beültethető kardioverter defibrilátor
 (it)Defibrillatore cardioverter implantabile (ja)埋込型自動除細動器 (ko)심장자동제세동기 (nl)Implantatie, ommeas
 -cardioverters-defibrillatorus (pl)Implantyjams kardioverters defibrilators (pt)Bateria de cardioversores
 desfibriladores (ro)Implanteerbare cardioverterdefibrilator (ru)Имплаптербар кардиовертер-дефрллятор
 (sl)Vsežozepni kardioverter defibrilator (pt-br)Cardioversor desfibrilador implantável
 (tr)otCardioversor defibrilador implantev (ro)Defibrilator cardiac implantabil (ru)Имплаптерувальн
 кардиовертер-дефрллятор (sk)Implantovateľný kardioverter defibrilátor (sv)Utsatt kardioverter defibrilator
 (th)ฝังยาต้านจังหวะหัวใจด้วยเครื่องกระชกไฟฟ้า (cs)Implanterbar elkovy terenní defibrilator
 (plmpian) Edepori Kardioverter Defibrilator (zh)植入式自动体外除颤器

VVED DDR

UDI


(01)05415067032843(17)241130(21)810055083

	2024-11-30
REF	CDDRA300Q
SN	810055083

5

STERILE EO


0°C 50°F | 45°C 113°F | 22°C 72°F

MD CE 0123

1.15

2022-11-16

EC REP

St Jude Medical
 Cardiac Rhythm
 Management Division
 15800 Valley View Court
 Sylmar, CA 91342 USA
 +1 818 362 6822

EC REP

St Jude Medical
 Coordinator Center BVBA
 The Corporate Village
 Da Vinciiaan 11 Box F1
 1935 Zaventem, Belgium
 +32 2 774 68 11

MY

St Jude Medical
 Operations (M) Sdn Bhd
 Plot 102, Lebuhraya Kajang Jajar
 Bayan Lepas Industrial Zone
 11800 Praang, Malaysia

Product of Malaysia


Handwritten signatures and stamps

DURATA™

REF 7122Q

65 cm

5

100126996

2024-11-30

2021-12-08

SN BPA140836


STERILE EO

50°C 122°F
5°C 23°F

CE
0123

Made in USA

St. Jude Medical
Cardiac Rhythm
Management Division
15900 Valley View Court
Sydney, CA 91342 USA
+1 818 362 6822
Pat. <http://www.abbott.com/patents>

EC REP
St. Jude Medical
Coordination Center BVEA
The Corporate Village
De Vrieslaan 11 Box F1
1925 Zaventem, Belgium
+32 2 774 69 11

ST. JUDE MEDICAL

Handwritten signature

Handwritten signature

TENDRIL™ STS

Pacing Lead


2088TC

52 cm

Active-fixation
Bipolar
Steroid-eluting
Endocardial
Pacing lead

IS-1 BI

Extendable/retractable
Optim™ insulation with Fast-Pass™ coating
Fast-Pass™ coating improves the lubricity of leads

Minimum introducer size: 6 F / 2 mm

5 Contains less than 1.0 mg dexamethasone sodium phosphate

100140460

2024-11-30

2021-12-02

SN CNY978923

STERILE EO

50°C / 122°F
-5°C / 23°F

CE
0123

CAUTION

Not to be used when a single dose of 1.0 mg dexamethasone sodium phosphate is contraindicated.

CONTENTS

One lead, accessories and product literature.

Made in Malaysia

St. Jude Medical
Cardiac-Rhythm
Management Division
15900 Valley View Court
Syomar, CA 91347 USA
+1 818 362 6827
Pat. Reg. www.abtson.com/patents

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
"The Abbott NGQ ICD devices (Avant™ VR/DR/HF, Neutrino™ NxT VR/DR/HF, Gallant™ VR/DR/HF, and Entrant™ VR/DR/HF) and CRT-D (Cardiac Resynchronization Therapy Defibrillator) devices, along with compatible leads, constitutes the implantable portion of the NGQ IMD system. The IPG (Implantable Pulse Generator) works with the Merlin Patient Care System (PCS) programmer (Model 3650) with USB Bluetooth Dongle (IT Infrastructure component) and software model 3330 v25.0.2 or higher (reference concurrent change notification CN12-19) and the existing Merlin Inductive Telemetry Wand (Model 3630). The IPGs also work with the Merlin™ 2 PCS Programmer Model MER3700, and the Merlin™ 2 PCS Programmer Model MER3400 v1.1.2 software supports applications for the Unity, NGQ, and ICM device families, as well as the accessory PSA tool used to assess the integrity of the implanted lead system (refer to CN06-20 for Merlin 2 approval). The myMerlinPulse™ mobile application for patients and the Merlin.net web portal with MN5000 Device Data Translator (DDT) software (reference concurrent change notification CN13-19) constitute the Remote Care system for NGQ(Next Generation Quad) devices. The primary function of ICDs and CRT-Ds are to treat life-threatening ventricular tachyarrhythmias via anti-tachycardia pacing (ATP), or cardioversion/defibrillation. Implantable Cardioverter Defibrillators (ICDs) can either be a dual chamber (also known as DR if rate responsive) device that is connected to a right atrial (RA) pacing lead and right ventricular defibrillation lead or single chamber ICD (also known as VR if rate responsive) with only a defibrillation lead in the right ventricle (RV). The ICD, more generally known as an implantable pulse generator (IPG), is typically implanted in a subcutaneous device pocket located in the pectoral region and is attached to the pacing/defibrillator leads that are implanted through the transvenous technique with the distal tips positioned in the RA (dual chamber) and RV. Cardiac Resynchronization Therapy Defibrillator (CRT-D) systems consist of a right atrial (RA) pacing lead, a right ventricular (RV) defibrillation lead, as well as a third lead placed in the coronary venous system overlying the left ventricle (LV) to improve ventricular resynchronization through biventricular pacing for treatment of heart failure. CRT-D devices, more generally known as an implantable pulse generators (IPG), are typically implanted in a subcutaneous device pocket located in the pectoral region and are attached to the pacing/defibrillator leads that are implanted through transvenous technique"

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