
Figure 8. Mitral sizer

Surgical Guidelines

The actual choice of surgical technique, modified in accordance with the instructions described herein, is left to the discretion of the individual surgeon.

When implanting supra-annular valves, non-everting mattress sutures are recommended.

Avoid any contact between the implantation sutures or knot tying technologies and the leaflets.

Precautions

- Do not allow the valve tissue to dry. Place the valve in isotonic sterile saline rinse solution immediately upon removal from the valve storage solution. Once removed from this solution, the valve should be periodically irrigated during implantation.
- Do not use cutting edge needles, unprotected forceps, or sharp instruments as they may cause structural damage to the valve.
- Never handle the leaflet tissue.
- Do not lacerate the valve tissue. If a valve is damaged, the valve must be explanted and replaced.
- Do not attempt to repair a valve. Damaged valves must not be used.

Aortic Valve Implantation

The Epic™ Plus aortic valve is designed to allow intra-annular placement of the inflow edge of the valve with supra-annular placement of the sewing cuff. The Epic Plus Supra valve is designed for supra-annular implantation of both the valve and the sewing cuff.

1. After sizing, choose a valve of the appropriate size.


CAUTION: Position the valve so that the stent posts do not obstruct the coronary ostia.

2. To facilitate implantation, the valve holder handle may be removed from the valve holder by depressing the release button on the valve holder (Figure 9).

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Archie


Figure 9. Release the valve holder handle

1. Release button

3. To remove the valve holder from the valve, cut the three retaining sutures as shown in Figure 10, and pull the valve holder away from the valve. Examine the valve to ensure that there are no holder suture remnants.

Figure 10. Cut three sutures to remove the aortic holder

Mitral Valve Implantation


1. After sizing, choose a valve of the appropriate size.
CAUTION: Position the mitral valve in a manner to avoid commissure obstruction of the left ventricular outflow tract, and minimize any potential of commissure contact with the ventricular wall.
2. To facilitate insertion of the mitral valve into the annulus, the mitral valve stent posts may be temporarily deflected inward during implantation. To deflect the valve stent posts inward, rotate the valve holder handle in a clockwise direction (Figure 11) while

Notes
1/2/16
Notes

securing the valve holder. The valve stent posts will remain deflected until the valve retaining sutures are cut. To facilitate implantation, the valve holder handle may be removed from the valve holder by depressing the release button on the valve holder.

NOTE: Take care to avoid looping or entangling sutures around the commissural posts, as this may result in a compromise of leaflet function.


Figure 11. Rotate the valve holder handle to deflect the mitral stent posts inward

3. To remove the valve holder from the valve, make one cut where the three retaining sutures meet as shown in Figure 12, and pull the valve holder away from the valve. Examine the valve to ensure that there are no holder suture remnants.

NOTE: To facilitate removal of the holder, the sewing cuff may be stabilized in position with a gloved hand or a surgical instrument, such as forceps, while pulling the valve holder away from the valve.

Figure 12. Cut the sutures to remove the mitral holder

Intra-Operative Assessment

The suggested method for assessing competence of the Epic Plus and the Epic Plus Supra valve is with intra-operative Doppler echocardiography.

Disposal

This instructions for use is recyclable. Dispose of all packaging materials as appropriate. Dispose of valves and accessories per standard solid biohazard waste procedures.

Notes
disposal

Notes

Patient Registration

A medical device registration form is included with each device. After implantation, please complete all requested information, and return the original form to the address indicated on the medical device registration form. Tracking by manufacturers is mandatory within the United States.

Future Valve-in-Valve Considerations

A transcatheter heart valve implanted in a surgical aortic tissue valve (ViV) is emerging as a reasonable treatment option for severely symptomatic patients with a failed surgical valve who are judged by the heart team to be at high or prohibitive risk of surgical reoperation, and in whom improvement in hemodynamics is anticipated. The information in this section is intended to be helpful for future ViV procedures in Epic™ Plus valves.

CAUTION: When performing a ViV procedure, the instructions for use of the transcatheter heart valve must be followed. The safety and effectiveness of valve-in-valve procedures in an Epic™ Plus Valve or an Epic™ Plus Supra valve have not been established.

Epic Plus™ and Epic Plus™ Supra valves contain a radiopaque wire marker at the annulus to aid the interventional cardiologist when placing the transcatheter prosthesis within the Epic Plus valve. Physicians may use the annulus marker as a reference point, in conjunction with known transcatheter stent foreshortening properties, to determine the proper location to place the transcatheter valve prior to delivery system actuation. The Epic Plus marker can continue to be utilized while adjusting the transcatheter valve during the deployment process.

Surgical Valve Implantation in the Aortic Position

When implanting the Epic™ Plus or Epic™ Plus Supra valves in the aortic position, it is important to:

- Ensure that there is adequate clearance of the coronary artery openings
- Assess the suitability of the implanted surgical valve and anatomy for a future ViV intervention

The replica end of the sizer accurately represents the valve and should be used to assess the position of the valve relative to the encountered anatomy and thereby ensure proper coronary opening clearance and assess suitability of the valve for a future ViV intervention.

Surgical Valve Implantation in the Mitral position

When implanting the Epic™ Plus valve in the mitral position, it is important to:

- Ensure that there is no obstruction of the left ventricular outflow tract and minimize any potential of commissure contact with the ventricular wall
- Assess the suitability of the implanted surgical valve and anatomy for a future ViV intervention

Three-dimensional CT imaging is the recommended measurement method to assess a failed surgical valve dimensions when selecting the appropriate transcatheter heart valve size during a ViV procedure. The reference dimensions in Table 1 do not account for valve tissue thickening, calcium burden/distribution, or other factors that vary by patient and are important when selecting a transcatheter heart valve size.

Limited Warranty


Abbott Medical warrants that reasonable care has been used in the manufacturing of this device. THIS WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES NOT EXPRESSLY SET FORTH HEREIN, WHETHER EXPRESS OR IMPLIED BY OPERATION OF LAW OR OTHERWISE INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, since handling, storage, cleaning, and sterilization of this device as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Abbott's control directly affect this device and the results obtained from its use. ABBOTT SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE directly or indirectly arising from the use of this device other than the replacement of all or part of it. SJM neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.

Some states in the United States do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you. This limited warranty gives you specific legal rights, and you may have other rights which vary from jurisdiction to jurisdiction.

Descriptions of specifications, appearing in Abbott literature, are meant solely to generally describe the device at the time of manufacture and do not constitute any express warranties.




Symbols

The following symbols are used on the valve labels.

Symbol	Definition
	Stented Tissue Valve
	Rinse - 2 x 500 mL x 10 seconds
	Contents

Notes:
[Handwritten signatures and notes]

The following symbols are used on the valve labels.

AORTIC	Aortic
	Contains biological material of animal origin
MITRAL	Mitral
	Manufacturing facility
	Temperature Indicator: If red, no not use
	Do not use if package is damaged
	Do not re sterilize
STERILE A	Sterilized using aseptic processing techniques
STERILE LC	Sterilized by liquid chemical sterilant
SN	Serial number
	Temperature limitations
	Use by
	Date of Manufacture
REF	Catalog Number
	Follow the instructions for use
	Do not reuse
	MR Conditional
UDI	Unique device identification number
	Non-pyrogenic
 <small>medical.abbott/manufact</small>	Follow instructions for use on this website
R ONLY	USA only: Federal law restricts this device to sale by or on the order of a Physician
	Quantity, package contents
	Manufacturer
	Contains the presence of Formaldehyde

Aortic
Manufacturing facility
ONLY

Abbott Medical
177 County Road B, East
St. Paul, MN 55117 USA
+1 855 478 5833
+1 651 756 5833

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
André

Abbott

Abbott Medical (Thailand) Co., Ltd.
Rasa Tower II, 17th Floor, Unit 1704
555 Phaholyothin Road, Chatuchak
Chatuchak, Bangkok 10900, Thailand

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The Epic™ Plus valve is bioprosthetic heart valves manufactured from select porcine aortic valve cusps. The cusps are carefully matched for optimum leaflet coaptation and hemodynamics.

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**U.S. FOOD & DRUG
ADMINISTRATION**

CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

U.S. Food & Drug Administration
10903 New Hampshire Avenue
Silver Spring, MD 20993
www.fda.gov

Certificate No. 10964-6-2023

CERTIFICATE TO FOREIGN GOVERNMENT

In order to allow the importation of United States products into foreign countries, the U.S. Food and Drug Administration (FDA) certifies the following information concerning the product(s) to be exported listed below:

Name of Product(s)

Name of Manufacturer/Distributor, Address

See Attached List

See Attached List

(One Page)

(One Page)

The product(s) described above (and the manufacturing/distribution site(s) which produces/distributes it) is subject to the jurisdiction of the FDA under the Federal Food, Drug, and Cosmetic Act.

It is certified that the above product(s) may be marketed in, and legally exported from, the United States of America at this time. The manufacturing plant(s) in which the product(s) is produced is subject to periodic inspections. The last such inspection showed that the plant(s), at that time, appeared to be in substantial compliance with current good manufacturing practice requirements for the product(s) listed above.

Sincerely,

CDR Cesar A. Perez, PhD, Director
DRP2: Division of Establishment Support
Office of Regulatory Programs
Office of Product Evaluation and Quality
Center for Devices and Radiological Health
U.S. Food and Drug Administration, DHHS

This certificate is valid from June 26, 2023 to June 25, 2025.

Andres

Andres

**U.S. FOOD & DRUG
ADMINISTRATION**

CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

U.S. Food & Drug Administration
10903 New Hampshire Avenue
Silver Spring, MD 20993
www.fda.gov

Certificate No. 10964-6-2023

Certificate to Foreign Government - Name of Manufacturer/Distributor Attachment Page 1 of 1

Name of Manufacturer

Legal Manufacturer

Abbott Medical
177 County Road B East
Saint Paul, MN
USA 55117

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---END OF MANUFACTURER/DISTRIBUTOR LIST---

Arlos
Arlos

Arlos

**U.S. FOOD & DRUG
ADMINISTRATION**

CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

U.S. Food & Drug Administration
10903 New Hampshire Avenue
Silver Spring, MD 20993
www.fda.gov

Certificate No. 10964-6-2023

Certificate to Foreign Government - Name of Product(s) Attachment Page 1 of 1

Name of Manufacturer

Legal Manufacturer

Abbott Medical
177 County Road B East
Saint Paul, MN
USA 55117

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QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

This is to certify that:

Abbott Medical
177 County Road B East
St. Paul
Minnesota
55117
USA

Facility ID Number: F000415

Holds Certificate No: **MDSAP 719341**

The company listed on this certificate has been audited to and found to conform with ISO 13485:2016 including the following country specific requirements:

Australia: Therapeutic Goods (Medical Devices) Regulations, 2002, Schedule 3 Part 1 (excluding Part 1.6) - Full Quality Assurance Procedure

Brazil: RDC ANVISA n. 16/2013, RDC ANVISA n. 23/2012, RDC ANVISA n. 67/2009

Canada: Medical Devices Regulations - Part 1 - SOR 98/282

Japan: MHLW Ministerial Ordinance 169, Article 4 to Article 68, PMD Act

USA: 21 CFR 820, 21 CFR 803, 21 CFR 806, 21 CFR 807 - Subparts A to D, 21 CFR 821

Please see scope page.

For and on behalf of BSI:

Graeme Tunbridge, Senior Vice President Medical Devices

Original Registration Date: 2020-01-27

Effective Date: 2023-01-27

Expiry Date: 2026-01-26

BSI Group America Inc. is an MDSAP recognised auditing organization

Page: 1 of 4

...making excellence a habit.™

This certificate remains the property of BSI and shall be returned immediately upon request.
An electronic certificate can be authenticated [online](https://www.bsigroup.com/ClientDirectory). Printed copies can be validated at www.bsigroup.com/ClientDirectory
To be read in conjunction with the scope above or the attached appendix.

Americas Headquarters: BSI Group America Inc., 12950 Worldgate Drive, Suite 800, Herndon, VA 20170-6007 USA
A Member of the BSI Group of Companies.

Certificate No: **MDSAP 719341**

Registered Scope:

Design, Development, and Manufacture of finished Mechanical Heart Valves, Valved Grafts, Tissue Heart Valves, Annuloplasty Rings, Valve and Annuloplasty Ring Sizer Sets, Mechanical Valve Leaflet Testers, Holder Rotators, Transcatheter Heart Valves with associated Delivery and Loading Systems, Mitral and Tricuspid Valve Repair Systems and related Accessories, Occluders and their Accessories, Plugs, Delivery and Exchanges Systems, Catheters, Sheaths, Introducers, Guidewires, Sizing Balloons for minimal invasive treatment of cardiovascular defects and peripheral vasculature diseases for the area of cardiology and peripheral vascular medicine. Manufacture of Sizing Plates and intermediate components used in other medical devices. Manufacture of Vascular Closure Device.

Notary
[Faint circular notary seal]

Original Registration Date: 2020-01-27

Effective Date: 2023-01-27

Expiry Date: 2026-01-26

Page: 2 of 4

This certificate remains the property of BSI and shall be returned immediately upon request.
An electronic certificate can be authenticated [online](http://www.bsigroup.com/ClientDirectory). Printed copies can be validated at www.bsigroup.com/ClientDirectory.
To be read in conjunction with the scope above or the attached appendix.

Americas Headquarters: BSI Group America Inc., 12950 Worldgate Drive, Suite 300, Herndon, VA 20170-6007 USA
A Member of the BSI Group of Companies.

Notary

Certificate No: **MDSAP 719341**

Location	Registered Activities
Abbott Medical 3885 Bohannon Drive Menlo Park California 94025 USA Facility ID Number: F001775	Design, Development and Manufacture of sterile Mitral and Tricuspid Valve Repair Systems and related Accessories.
Abbott Medical Lot 20, B Street Caguas West Industrial Park Caguas, PR 00725 Puerto Rico Facility ID Number: F000417	The Manufacture of Annuloplasty Ring, Mechanical Heart Valves and Contract Manufacturing of Vascular Closure Device.
Abbott Medical Costa Rica Limitada Edificio #44 Calle 0, Ave. 2 Zona Franca Coyol El Coyol, Alajuela Costa Rica Facility ID Number: F002736	Manufacture of Occluders and their Accessories and Tissue Heart Valves.
Abbott Medical 1820 Bastian Court Westfield Indiana 46074 USA Facility ID Number: F006224	Manufacture of sterile Mitral Valve Repair Systems and related Accessories.
Abbott Medical 177 County Road B East St. Paul Minnesota 55117 USA Facility ID Number: F000415	Design, Development, and Manufacture of finished Mechanical Heart Valves, Valved Grafts, Tissue Heart Valves, Annuloplasty Rings, Valve and Annuloplasty Ring Sizer Sets, Mechanical Valve Leaflet Testers, Holder Rotators, Transcatheter Heart Valves with associated Delivery and Loading Systems, Mitral and Tricuspid Valve Repair Systems and related Accessories, Occluders and their Accessories, Plugs, Delivery and Exchanges Systems, Catheters, Sheaths, Introducers, Guidewires, Sizing Balloons for minimal invasive treatment of cardiovascular defects and peripheral vasculature diseases for the area of cardiology and peripheral vascular medicine. Manufacture of Sizing Plates and intermediate components used in other medical devices. Contract Manufacture of Vascular Closure Device.

Not for Distribution

Original Registration Date: 2020-01-27 Effective Date: 2023-01-27 Expiry Date: 2026-01-26

This certificate remains the property of BSI and shall be returned immediately upon request.
An electronic certificate can be authenticated [online](http://www.bsigroup.com/ClientDirectory). Printed copies can be validated at www.bsigroup.com/ClientDirectory
To be read in conjunction with the scope above or the attached appendix.

AMH

Certificate No: **MDSAP 719341**

Location

Registered Activities

St. Jude Medical Brasil Ltda.
Rua Professor Jose Vieira de Mendonca, 1301
Bairro Engenho Nogueira
Belo Horizonte
Minas Gerais
31.310-260
Brasil
Facility ID Number: F003466

The Manufacture and Final Inspection of Tissue made Heart Valves and Tissue Vascular Prostheses.

Handwritten signature
01/26/2026

Original Registration Date: 2020-01-27

Effective Date: 2023-01-27

Expiry Date: 2026-01-26

Page: 4 of 4

This certificate remains the property of BSI and shall be returned immediately upon request.
An electronic certificate can be authenticated [online](http://www.bsigroup.com/ClientDirectory). Printed copies can be validated at www.bsigroup.com/ClientDirectory
To be read in conjunction with the scope above or the attached appendix.

Americas Headquarters: BSI Group America Inc., 12950 Worldgate Drive, Suite 800, Herndon, VA 20170-6007 USA
A Member of the BSI Group of Companies.

Handwritten signature

Date 1 August 2022

Letter of Authorization of Distributor
Thailand


To Whom It May Concern:

This is to certify that St. Jude Medical, LLC ("St. Jude Medical"), having a location at One St. Jude Medical Drive, St. Paul, Minnesota, 55117, USA, authorizes its affiliate, Abbott Medical (Thailand) Co., Ltd (previously known as "St. Jude Medical (Thailand) Co. Ltd."), located at Rasa Tower II, 17th Floor, Unit 1704, 555 Phaholyothin Road, Kwang Chatuchak, Chatuchak, Bangkok 10900, in Thailand, to distribute products of St. Jude Medical in Thailand until further notice.

As such, Abbott Medical (Thailand) Co., Ltd, is authorized to undertake the necessary actions to participate in all auctions and tenders regarding the supply of medical equipment and products from any of the divisions of St. Jude Medical identified above, and, in connection with such activity, to conclude contracts for the supply of such products, to carry out deliveries according to the conditions stipulated in the contracts and tender offers, to appoint sub-distributors to distribute such products in Thailand, and also to carry out any maintenance service, training of staff and to fulfill guarantee obligation, as needed and appropriate.

St. Jude Medical's authorization is valid for all the product lines of the St. Jude Medical Divisions until St. Jude Medical provides further written notification confirming that the services of Abbott Medical (Thailand) Co., Ltd. for the above matters are no longer needed.



Yours sincerely,

Name: John A. McCoy Jr.
Position: Vice President and Treasurer

STATE OF ILLINOIS
COUNTRY OF LAKE

Signed and attested before me on 1 August, 2022 by John McCoy Jr.

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Beneficiary Name : Abbott Medical (Thailand) Co., Ltd.
Incorporated in U.S.A with Limited Liability