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
Avalus™
Bioprosthesis


**PROVEN
PLUS.**

With more than 40 years of heart valve innovations, we took proven valve design concepts and adapted them for excellent implantability for you and performance for your patients.

Medtronic
Further Together

17180


Oh


INTRODUCING
THE AVALUSTM
AORTIC VALVE
BY MEDTRONIC.

YOU WANT THE VERY BEST FOR YOUR PATIENTS. SO DO WE.


We designed the next generation bovine pericardial valve for better overall performance, improved implant experience, and a contemporary design to facilitate future valve-in-valve (ViV).

Interior-mounted leaflets minimize damaging contact with the frame — a design platform for long-term durability.

AOA™ tissue treatment* to mitigate calcification — over 20 years of clinical use on the Medtronic surgical tissue valve portfolio.^{1,2}

3.2

- Supra-annular design to enhance hemodynamics³
- Three laser cut bovine pericardial leaflets matched for thickness and deflection to provide consistent performance
- Two-part polymer frame minimizes stress zones on leaflets
- Sewing markers facilitate suture placement and valve orientation


*No clinical data is available which evaluates the long-term impact of AOA treatment in patients.

Ch
J. 12.1



PLUS

Ease of Implant for You

- Soft and pliable sewing cuff facilitates needle penetration, suture placement, and valve seating for an improved implant experience
- Lower valve profile and narrow commissure posts expand ostia clearance and give you more space for knot tying
- Streamlined valve holder improves visibility in both standard and minimally invasive approaches
- Simple one-cut release


of the valve holder is designed to provide a clear view of the valve and the sewing cuff, and to facilitate the placement of the valve holder in the heart.


Q

15.80


Designed to achieve 100% coaptation and minimize central regurgitation.


Flexible support frame with firm base designed to maintain circularity and consistent hemodynamic performance.

Performance and Lifetime Management for Your Patients

- Valve dimensions and geometry enable future ViV replacements
- PEEK base frame impregnated with barium sulfate provides for radiopacity and visibility
- Polymer frame mitigates the risk of potential metal on metal corrosion with transcatheter stent materials
- MRI Safe in all MR environments without conditions


3.6

Ch

↓ 152

Ordering and Specifications

Avalus Valve Order Number	Valve Size	Stent Diameter (TAD)	Internal Orifice Diameter	External Sewing Ring Diameter	Valve Profile Height	Aortic Protrusion
		(1)	(2)	(3)	(4)	(5)
40019	19 mm	19 mm	17.5 mm	27.0 mm	13.0 mm	11.0 mm
40021	21 mm	21 mm	19.5 mm	29.0 mm	14.0 mm	12.0 mm
40023	23 mm	23 mm	21.5 mm	31.0 mm	15.0 mm	13.0 mm
40025	25 mm	25 mm	23.5 mm	33.0 mm	16.0 mm	14.0 mm
40027	27 mm	27 mm	25.5 mm	36.0 mm	17.0 mm	15.0 mm

TAD – Tissue Annulus Diameter

Accessories

Order Number	Description
7420	Valve Handle
7400S	Avalus Sizers
T7400	Tray, Accessory, Avalus

References

1. Medtronic Freestyle® Aortic Root Bioprosthesis was first implanted clinically in August 1992. Freestyle Aortic Root Bioprosthesis 15-Year Clinical Compendium. ©2016 Medtronic.
2. Jamieson WR, Riess FC, Raudkivi PJ, et al. Medtronic Mosaic porcine bioprosthesis: assessment of 12-year performance. *J Thorac Cardiovasc Surg.* August 2011;142(2):302-7.
3. Ruzicka DJ, Hettich I, Hutter A, et al. The complete supraannular concept. *Circulation* 2009;120[suppl 1]:S139-S145.

For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use.

Medtronic

Europe
 Medtronic International Trading Sàrl
 Route du Molliou 31
 Case Postale
 CH-1131 Tolochenaz
 Switzerland
 Tel: +41 (0)21 802 70 00
 Fax: +41 (0)21 802 79 00

UC201608716aEE © 2017 Medtronic.
 All rights reserved.

medtronic.com

de


ha

Bioprosthesis

1. Device description

2.3
 The Avalus™ bioprosthesis, Model 400, consists of a polyester-covered base frame and trileaflet support frame structure shaped to the geometry of a trileaflet valve. The base frame and trileaflet support frame are injection-molded using polyetheretherketone (PEEK) material. The base frame PEEK material is impregnated with barium sulphate to allow for radiographic visualization. The leaflets are laser cut from bovine tissue that has been cross-linked in buffered glutaraldehyde. The leaflets are inserted between the trileaflet support frame and the base frame, and then all components are securely sutured together. A sewing ring, fabricated from polyester cloth, is integrated into the inflow base frame cover to allow the user to suture and seat the Avalus bioprosthesis in the supra-annular position. Sewing markers are located on the sewing ring in the mid sinus area of each cusp to provide guidance for even spacing of the implant sutures (Figure 1). A disposable holder is attached to the outflow of the valve to facilitate implantation. The aortic bioprosthesis (valve) is treated with an alpha amino oleic acid (AOA™) antimineralization process that has been shown to mitigate leaflet calcification in animal studies. The disposable holder is designed to fit the reusable Medtronic valve handle, Model 7420. The holder features a single cut point to remove the holder from the valve. Double-ended sizers are used to select the appropriate valve size. The barrel end represents the valve orifice. The replica end imitates the prosthesis geometry.

Table 1. Available sizes and sewing ring diameters (nominal values in millimeters)

Size	19 mm	21 mm	23 mm	25 mm	27 mm
(1) Stent diameter (TAD)	19	21	23	25	27
(2) Internal orifice diameter	17.5	19.5	21.5	23.5	25.5
(3) External sewing ring diameter	27.0	29.0	31.0	33.0	36.0
(4) Valve profile height	13.0	14.0	15.0	16.0	17.0
(5) Aortic protrusion	11.0	12.0	13.0	14.0	15.0

Note: TAD—tissue annulus diameter

2. Indications for use

1, 2.1
 The Avalus bioprosthesis is indicated for the replacement of diseased, damaged, or malfunctioning native or prosthetic aortic valves.

3. Contraindications

None known.

4. Warnings and precautions

4.1. Warnings

- This device was designed for single use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Do not resterilize the valve by any method. Exposure of the valve and container to irradiation, steam, ethylene oxide, or other chemical sterilants will render the valve unfit for use.
- Check the shipping temperature indicator inside the carton. If the shipping temperature indicator window is black, the valve is not suitable for clinical use.
- Do not use the valve in any of the following circumstances:
 - The valve has been dropped, damaged, or mishandled in any way
 - The Use-by date has elapsed
 - All tamper strips are damaged
 - The serial number tag does not match the number on the container label
 - The glutaraldehyde storage solution does not completely cover the valve
 - Fluid is leaking from the packaging
 - The valve shows signs of deterioration
- Do not handle the tissue portion of the valve with instruments. Extreme care must be taken to prevent damage to the delicate valve tissue. Even a minor perforation could enlarge in time to significantly impair valve function.
- Do not attempt to repair a damaged valve. A damaged valve must not be used.
- Do not add antibiotics or any other substance to the storage solution or the rinse solution.

De

✓ TAD

It is recommended that prophylactic antibiotic therapy be given to patients undergoing dental or other procedures which are potentially bacteremic in order to minimize the risk of endocarditis.

Encourage patients to carry the Implanted Device Identification Card, provided by Medtronic, with them at all times.

11. How supplied

11.1. Packaging

3.5

The Avalus bioprosthesis is chemically sterilized and is supplied sterile in a buffered 0.2% glutaraldehyde solution. Sterility is compromised if the glass jar-and-lid container is opened or damaged. The outside of the container is not sterile. Do not place it in the sterile field.

11.2. Storage

Store the Avalus bioprosthesis between 5°C and 25°C (41°F and 77°F). Refrigeration is not required, and freezing could damage the valve. Room temperature storage up to 25°C (77°F) is satisfactory. Do not expose the valve to sunlight or other ultraviolet light sources, or placed where significant temperature fluctuations could occur.

Maintain appropriate inventory control so that bioprostheses with earlier Use-by dates are implanted first to avoid expiration dates.

11.3. Return of explanted bioprostheses

Medtronic would like to obtain recovered Avalus bioprostheses. When determined to be appropriate, explants will be studied by a consulting pathologist. A written report summarizing the findings will be returned to the physician. Product return kits, including an explant information form, are available by contacting Medtronic distribution centers or a Medtronic sales representative. It is important that the explant form is filled out completely. If a kit is not available, place the explanted bioprosthesis in a container of glutaraldehyde or 10% buffered formalin immediately after excision. For further instructions on the return of an explanted device, contact a Medtronic sales representative.

12. Patient information

12.1. Registration information

Note: Patient registration does not apply in countries where patient privacy laws conflict with providing patient information, including countries from the European Union.

A Patient Registration Form is included in each device package. After implantation, please complete all requested information. Return the original form to the Medtronic address indicated on the form. Provide the patient with the temporary identification card before discharge.

An Implanted Device Identification Card is provided to the patient. This card contains the name and telephone number of the patient's physician, as well as information that medical personnel would require in the event of an emergency.

13. Postoperative information

13.1. Magnetic resonance imaging (MRI) compatibility

The Avalus bioprosthesis Model 400, is magnetic resonance (MR) safe. The valve contains no metal and, therefore, poses no known hazards in all MR environments.

MR scanning at 3.0 T and 1.5 T can be performed immediately after implantation. The valve will not cause any harm to the patient when exposed to MR scanning.

13.2. Image artifact

MR image quality could be compromised if scanning on or near the position of the valve. If necessary, optimize the MR imaging parameters for the presence of the valve.

14. Accessories

Use only Avalus sizers, Model 7400S, and the Medtronic valve handle, Model 7420, to determine the appropriate Avalus bioprosthesis size.

Caution: Do not use the sizers or handles until they have been thoroughly cleaned and sterilized. Refer to the appropriate instructions for use for further instructions.






Caution: Do not use other manufacturers' valve sizers, or sizers from another Medtronic prosthesis, to size the Avalus bioprosthesis.

Al

3.5

Explanation of symbols on package labeling


Refer to the device labeling to see which symbols apply to this product.

	Do not reuse
	Use-by date
	Sterile LC: Device has been sterilized using liquid chemical sterilants according to EN/ISO 14160.
	Temperature limit
	Size
	Serial number
	Catalog number
	Date of manufacture
	Manufacturer
	MR Safe
	Quantity
	Nonpyrogenic
	Do not resterilize
	Do not use if indicator turns black
	Authorized representative in the European Community
	Manufactured in
	For US audiences only
	Model

Q

h.18

ZERTIFIKAT ◆ CERTIFICATE ◆ 認 證 證 書 ◆ CERTIFICADO ◆ CERTIFICAT

Product Service

EU Technical Documentation Assessment Certificate (MDR)

Pursuant to Regulation (EU) 2017/745 on Medical Devices, Annex IX Chapter II (Implantable Class IIb Devices and Class III Devices)

No. **G70 039709 1408 Rev. 00**

Manufacturer: **Medtronic, Inc.**
710 Medtronic Parkway
Minneapolis, MN 55432
USA

SRN Manufacturer: US-MF-000019977

Authorized Representative: Medtronic B.V.
Earl Bakkenstraat 10, 6422 PJ Heerlen, THE NETHERLANDS

The Certification Body of TÜV SÜD Product Service GmbH certifies that the manufacturer has drawn up and presented a Technical Documentation according to Annex II and III of the Regulation (EU) 2017/745 on medical devices. Details on devices covered by the Technical Documentation are described on the following page(s). The Report referenced below summarises the result of the assessment and includes reference to relevant CS, harmonized standards and test reports. The technical documentation assessment included an assessment of the clinical evaluation assessment. The conformity assessment has been carried out according to Annex IX chapter II of this regulation with a positive result.

Changes to the approved device, where such changes could affect the safety and performance of the device or the conditions prescribed for use of the device, shall require approval from the notified body TÜV SÜD Product Service GmbH. In order to place the devices on the market with CE-marking, an EU Quality Management System Certificate pursuant to Annex IX chapters I and III is necessary in addition to this EU Technical Documentation Assessment Certificate. All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with. For details and certificate validity see: www.tuvsud.com/ps-cert?q=cert:G70_039709_1408_Rev_00


Report No.: 713210950

Valid from: 2022-11-18

Valid until: 2027-11-17

Christoph Dicks
Head of Certification/Notified Body

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Product Service

EU Technical Documentation Assessment Certificate (MDR)

Pursuant to Regulation (EU) 2017/745 on Medical Devices, Annex IX Chapter II
(Implantable Class IIb Devices and Class III Devices)

No. G70 039709 1408 Rev. 00

Classification: III
 Device Group: P070301 - BIOLOGICAL CARDIAC VALVES
 Basic UDI-DI: 0763000B00007568C
 Intended Purpose: The intended purpose of the Avalus bioprosthesis is to restore unidirectional blood flow through the aortic annulus by replacing a diseased or damaged heart valve.
 Device(s): Avalus™ 400 Bioprosthesis
 40019
 40021
 40023
 40025
 40027
 40029

Classification: III
 Device Group: P070380 - HEART VALVES - ACCESSORIES
 Basic UDI-DI: 0763000B00010927L
 Intended Purpose: The intended purpose of the Avalus sizers 7400S is to facilitate the implantation of the Avalus bioprosthesis.
 Device(s): Avalus™ 7400S Sizers

The validity of this certificate depends on conditions and/or is limited to the following: ./.

TUV®

ZERTIFIKAT ◆ CERTIFICATE ◆ 證書 ◆ CERTIFICADO ◆ CERTIFICAT


Product Service

Certificate

No. Q6 039709 1376 Rev. 00

Holder of Certificate: Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA

Certification Mark:


Scope of Certificate: Production of :

- Disposable Medical Devices used in Cardiopulmonary Bypass Procedures, Disposable Products for Biopsy
- Vascular Devices, Subassemblies for Stent Graft Systems, Iliac Stents and Occluder Systems for Endovascular
- External Drainage Products
- Heart Valves
- Diagnostic and Ablation Catheters and Transseptal Introduction Products
- Annuloplasty Rings and Bands
- Medical Devices for Minimally Invasive Spinal and Orthopedic Procedures.
- Catheters-Based Technologies for Peripheral Artery Disease and for Neurovascular Revascularization
- Detachable Embolization Coils and Delivery System
- Detachable Embolization Devices and Delivery System
- Cerebral Protection Devices

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system (excluding subclause 7.3), which meets the requirements of the listed standard(s). All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with. For details and certificate validity see: www.tuvsud.com/ps-cert?q=cert:Q6_039709_1376_Rev_00

Report No.: 72166396

Valid from: 2022-02-11
Valid until: 2024-12-20

C.D.M.

Date, 2022-02-11

Christoph Dicks
Head of Certification/Notified Body

Ch
VSE